## Portavo® 908 MULTI

Quickstart guide .....3

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Instructions succinctes 113

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Other languages: [www.knick.de](http://www.knick.de)

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**EN** Copyright 2019 • Subject to change

Version: 2

This document was last updated on February 7, 2019

The latest documents are available for download on our website under the corresponding product description.

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<b>Measuring</b> .....	<b>31</b>
Toggling the Measured Value Display .....	31
Adjusting the Temperature.....	31
<b>Printing</b> .....	<b>32</b>
Printing Currently Measured Values .....	32
Printing a Calibration Record.....	33
<b>Data Logger</b> .....	<b>34</b>
Operating Modes of the Data Logger (Logger Type) .....	35
Configuring the Data Logger.....	37
Increasing the Battery Life .....	37
Starting/Stopping the Data Logger.....	39
Viewing the Logger Data.....	39
Deleting the Logger Data.....	39
<b>Error and Status Messages</b> .....	<b>41</b>
"Sensoface" Messages.....	42
Info and Help Texts .....	42
Error Messages .....	43
<b>Option 001 SOP</b> .....	<b>44</b>
<b>Option 002 Temp.Cal</b> .....	<b>44</b>
<b>How to Enter the Rescue PIN Code</b> .....	<b>48</b>
<b>Specifications</b> .....	<b>50</b>
<b>Index</b> .....	<b>56</b>

## Return of Products Under Warranty

Please contact our Service Team before returning a defective device.

Ship the cleaned device to the address you have been given.

If the device has been in contact with process fluids, it must be decontaminated/ disinfected before shipment. In that case, please attach a corresponding certificate, for the health and safety of our service personnel.

## Disposal

Please observe the applicable local or national regulations concerning the disposal of "waste electrical and electronic equipment".

## Registered Trademarks

The following names are registered trademarks. For practical reasons they are shown without trademark symbol in this manual.

- CaliMat®
- Calimatic®
- Memosens®
- Paraly®
- Portavo®
- Sensocheck®
- Sensoface®

Check the shipment for transport damage and completeness (see Package Contents).

### NOTICE!

Do not operate the device when one of the following conditions applies:

- the device shows visible damage
- the device fails to perform the intended function
- prolonged storage at temperatures above +70 °C / +158 °F
- severe transport stresses

In this case, a professional routine test must be performed.

This test should be carried out at our factory.

## Inserting the Batteries

With four AA batteries, the Portavo has an operating time of up to 500 h when operated in logger mode (see page 37).

Open the battery compartment on the rear of the device. Be sure to observe the correct polarity when inserting the batteries (see markings in the battery chamber). Close the battery compartment cover and screw it handtight.

A special lithium-ion battery (ZU 0925) suited to the battery compartment is available for the Portavo 908. Only this battery type can be charged directly from the USB port.

### A battery icon in the display indicates the battery power level:


	Icon fully filled	Batteries at full capacity
	Icon partially filled	Battery capacity is sufficient
	Icon empty	Battery capacity not sufficient; calibration is possible, no logging
	Icon blinks	Only a few operating hours remaining, measurement is still possible <b>NOTICE!</b> It is absolutely necessary to replace the batteries.

## Connecting a Sensor

You can connect a Memosens sensor for pH/ORP, conductivity or oxygen measurement to the Portavo 908. The Model SE 340 optical oxygen sensor can also be connected. The meter automatically recognizes the connected sensor and accordingly selects the corresponding process variable. Memosens is signaled in the display. Note that only **one** sensor may be connected to the meter at a time.

### Separate Temperature Probe

After power-on, a separate temperature probe is automatically recognized. If you want to replace the temperature probe, you must switch off the meter and then switch it on again.

### Connections

- a - Micro USB port
- b - M8, 4 pins, for Memosens lab cable
- c - Temperature probe GND
- d - Temperature probe
- e - M12, 8 pins, for flexible connecting cable for Memosens sensors or SE 340 sensor (optical oxygen)

Memosens sensors have a **cable coupling**, which allows convenient replacement of sensors while the cable remains connected to the meter. The connecting cable is connected to socket **b** or **e**.

## Switching On the Meter

You can use **meas** or **on/off** to switch the meter on:

### Analog sensors:

- When you press **meas** or **on/off**, the meter immediately switches to measuring mode.


### Memosens sensors:

- When you press **meas** or **on/off**, the meter displays selected sensor data before it switches to measuring mode.

## Icons


Important information about the state of the device:

## Display and Keypad

Display and keypad correspond directly via softkeys.

### Menu Icons

Information

Calibration

Data logger

Configuration

- Softkeys** Function is shown in the display above the key
- Arrow keys** Selecting / Adjusting entries
- enter** Confirming an adjustment
- on/off** Switching on / off
- meas** Switching on / Immediate access of measuring mode / Toggling the display / Displaying time and date

pH

ORP

Oxy

Cond

- 1) In measuring mode, press the **Menu** softkey.
- 2) Select "Information" and confirm by pressing **enter**.
- 3) Select the desired submenu and confirm by pressing **enter**.

The different submenus are described below.

## Calibration Record

Shows the data of the last calibration performed on the currently connected sensor. When the right softkey has been set to "Print" in the Configuration menu, it can be used to print out the calibration record (via micro USB port).

## Sensor Information

Shows the data of the currently connected sensor. When MemoLog has been activated (in the Configuration menu), you can save the sensor data in the device by pressing the **Save** softkey. The following table shows the sensor information depending on the respective sensor:

	pH/ pH/ORP**	Cond	Oxy	ISFET	ORP	Optical Oxy
Manufacturer	x	x	x	x	x	x
Ref. No.	x	x	x	x	x	x
Sensor serial no.	x	x	x	x	x	x
Cap serial no.						x
TAG	x	x	x	x	x	
SW version	x	x	x	x	x	x
HW version	x	x	x	x	x	
Calibration*	x	x	x	x	x	x
Zero point	x		x			x
Slope	x		x	x		x
ORP calibration* **	x					
Correction					x	
Nom. cell constant		x				
Temp. offset	x	x	x		x	
Sensor operating time	x	x	x	x	x	x
Cap operating time						x
Wear	x		x	x		
SIP	x	x	x	x	x	
CIP	x**	x				
Autoclaving	x**					
Cell constant		x				
Operating point				x		

\* latest calibration \*\* for pH/ORP combo sensor only

pH

Oxy

## Sensor Network Diagram (pH and Oxy only)


Provides single-glance information on the following parameters of the connected sensor:

- Slope
- Zero point (operating point for Memosens ISFET)
- Sensocheck (pH) or leakage current (ISFET and Oxy)
- Response time
- Calibration timer
- Wear (Memosens)

Parameters which cannot be checked are shown as inactive (gray) and are set to 100 %.

The parameter values should lie between the outer (100 %) and inner (50 %) hexagon. When a value enters the inner hexagon (< 50 %), the corresponding caption text flashes red (see example).

Example: Network diagram of a digital pH sensor (Memosens)

pH

ORP

Oxy

Cond

## Sensor Monitor

Shows the raw values available from the connected sensor:

<b>pH</b>	mV, temperature, glass impedance
<b>pH ISFET</b>	mV, leakage current, temperature
<b>ORP</b>	mV, temperature
<b>Cond</b>	Resistance, conductance, temperature
<b>Oxy</b>	Sensor current, leakage current, polarization voltage, partial pressure, air pressure, temperature
<b>Oxy, optical</b>	Partial pressure, temperature

## Messages

Shows all active error and status messages as well as supplementary help texts.

## MemoLog (Memosens only)

Displays the individual calibration records stored in the device. You have the possibility to delete individual entries or all entries. The following parameters are displayed:

- Sensor type
- Serial no.
- TAG
- Calibration date
- Zero point
- Slope
- Cell constant (Cond sensor)
- Operating point (ISFET sensor)

**Background:** The device provides a calibration data logger, which must be activated in the configuration menu. With "MemoLog" activated, up to 100 calibration records can be directly saved to the device. After every calibration, the complete Memosens data will be recorded. Convenient management of the calibration data is possible using the MemoSuite or Paraly SW 112 software.

MemoLog is not suitable for SE 340 (optical oxygen sensor).

pH

ORP

Oxy

Cond

## Device Info

Shows the following device information:

- Device name
- Serial number
- Software version
- Hardware version
- Air pressure
- Battery

## Device Test

A device self-test is automatically run in the background at regular intervals.

It checks the memory modules listed below. A green checkmark shows that the test was successful.

- FLASH program memory
- FLASH data memory
- FLASH parameter memory
- RAM (working memory)

## Display Test

- 1) Select "Display test" and press **enter**.
- 2) The display lights up red, green, blue and then white.
- 3) Press any key to stop the test.

## Keypad Test

- 1) Select "Keypad test" and press **enter**.
- 2) Press all nine keys one after the other.  
A green checkmark shows that a key functions properly.
- 3) Press any key to stop the test.

## pH

## pH Configuration

- 1) In measuring mode, press the **Menu** softkey.
- 2) Select "Configuration" and confirm by pressing **enter**.
- 3) Make the desired adjustments.

The following table gives you an overview.

Factory settings are shown in **bold print**.

"pH Configuration" menu selection – part 1

Language	<b>Deutsch</b>   English   Español   Italiano   Français   Português   Русский
Auto-off	<b>Off</b>   5 min   10 min   30 min   60 min
Temperature	°C   °F
Right softkey	Logger Start/Stop   <b>Hold value</b>   Print
+ pH sensor*	
Display format	<b>0.00 pH</b>   0.000 pH
Wear	<b>On</b>   Off
+ Calibration*	
Cal timer	<b>Off</b>   On
Interval	On: 00 ... 99 days
Cal mode	<b>Calimatic</b>   Manual   Data entry
Cal points	<b>Auto</b>   1-point   2-point   3-point
Buffer set	Mettler-Toledo    2.00/4.01/7.00/9.21 <b>Knick CaliMat</b> 2.00/4.00/7.00/9.00/12.00 Ciba                2.06/4.00/7.00/10.00 NIST technical    1.68/4.00/7.00/10.01/12.46 NIST standard    1.679/4.006/6.865/9.180 Hach                4.01/7.00/10.01/12.00 WTW                2.00/4.01/7.00/10.00 Hamilton         2.00/4.01/7.00/10.01/12.00 Reagecon         2.00/4.00/7.00/9.00/12.00 DIN 19267        1.09/4.65/6.79/9.23/12.75 Metrohm          4.00/7.00/9.00 User buffer 1**
MemoLog	<b>Off</b>   On
TAG	<b>Off</b>   On

\* "+" indicates that submenus can be opened by pressing **enter**.

\*\* Parameter can be configured using Paraly SW 112 software.

“pH Configuration” menu selection – part 2

	+ Time/Date*			
	Time format	<b>24 h</b>   12 h		
	Date format	<b>dd.mm.yyyy</b>   yyyy-mm-dd   dd/mm/yyyy   mm/dd/yyyy		
	Time	hh:mm:ss		
	Date	Date format as configured		
	+ Display*			
	Appearance	<b>Modern</b>   Retro		
	Lighting	<b>Permanent</b>   60 min   30 min   10 min   5 min   1 min   30 sec		
	Brightness	<b>Bright</b>   Standard   Dim		
	+ Data logger*			
	Meas.point	--		
	Note	--		
	Recording	<b>Non-circular</b>   Circular		
	Logger type	Shot		
		<b>Interval</b>	00.00.01...12:59:59   <b>00:02:00</b>	
		Difference	1st difference	<b>On</b>   Off
			Delta pH	pH 0.0...16.0   <b>pH 1.0</b>
			Delta mV	0 ... 2000 mV   1 mV
			2nd diff.	On   <b>Off</b>
			Delta °C	0...99.9 °C   <b>1.0 °C</b>
		Delta °F	0...450 °F   <b>1.0 °F</b>	
	Intv+Diff	Interval	See logger type: Interval	
		Difference	See logger type: Difference	
	Limit value	Interval	Basis/Event 00.00.01...12:59:59 <b>00:01:00/00:00:01</b>	
		Limit values	Min/Max, corresponding to permissible range (see Specifications)	
+ Printer				
Driver	<b>Auto</b>   HP PCL   Epson   Samsung   IBM			
Paper	<b>A4</b>   Small			
+ Options	001 SOP	Add-on function, enabled via TAN		
	002 Temp.cal			
Factory setting	Yes   <b>No</b> <b>Note:</b> Reset to factory settings will also erase all logger data!			

\* "+" indicates that submenus can be opened by pressing **enter**.

## ORP

## ORP Configuration

- 1) In measuring mode, press the **Menu** softkey.
- 2) Select "Configuration" and confirm by pressing **enter**.
- 3) Make the desired adjustments.

The following table gives you an overview.

Factory settings are shown in **bold print**.


"ORP Configuration" menu selection – part 1

Language	<b>Deutsch</b>   English   Español   Italiano   Français   Português   Русский
Auto-off	<b>Off</b>   5 min   10 min   30 min   60 min
Temperature	°C   °F
Right softkey	Logger Start/Stop   <b>Hold value</b>   Print
+ Calibration*	
MemoLog	<b>Off</b>   On
TAG	<b>Off</b>   On
+ Time/Date*	
Time format	<b>24 h</b>   12 h
Date format	<b>dd.mm.yyyy</b>   yyyy-mm-dd   dd/mm/yyyy   mm/dd/yyyy
Time	hh:mm:ss
Date	Date format as configured
+ Display*	
Appearance	<b>Modern</b>   Retro
Lighting	<b>Permanent</b>   60 min   30 min   10 min   5 min   1 min   30 sec
Brightness	<b>Bright</b>   Standard   Dim

\* "+" indicates that submenus can be opened by pressing **enter**.

“ORP Configuration” menu selection – part 2

\* "+" indicates that submenus can be opened by pressing **enter**.

## Cond

## Conductivity Configuration

- 1) In measuring mode, press the **Menu** softkey.
- 2) Select "Configuration" and confirm by pressing **enter**.
- 3) Make the desired adjustments.

The following table gives you an overview.


Factory settings are shown in **bold print**.

"Conductivity Configuration" menu selection – part 1

	Language	<b>Deutsch</b>   English   Español   Italiano   Français   Português   Русский
	Auto-off	<b>Off</b>   5 min   10 min   30 min   60 min
	Temperature	°C   °F
	Right softkey	Logger Start/Stop   <b>Hold value</b>   Print
	+ Cond sensor*	
	Conductivity	<b>S/cm</b>   S/m
	Range selection	<b>Auto</b>   0.000 µS/cm   00.00 µS/cm   000.0 µS/cm   0000 µS/cm   00.00 mS/cm   000.0 mS/cm   0000 mS/cm
	Compensation	<b>Off</b>   MΩ cm   TC   SAL   TDS   % by wt
	TC comp.	TC: <b>Linear</b>   NLF   NaCl   HCl   NH3   NaOH
	TC of solution	TC: 0 ... 20.0 %/K   <b>2.1 %/K</b>
	Ref. temp	TC: 0 ... 100.0 °C   <b>25 °C</b>   32 ... 212 °F   <b>77 °F</b>
	TDS factor	TDS: 0 ... 9.99   <b>1.00</b>
	Solution	% by wt: <b>NaCl</b>   HCl   NaOH   H2SO4   HNO3

\* "+" indicates that submenus can be opened by pressing **enter**.

“Conductivity Configuration” menu selection – part 2

\* "+" indicates that submenus can be opened by pressing **enter**.

Cond

“Conductivity Configuration” menu selection – part 3

	+ Data logger*
	Meas.point
	Note
	Recording
	Logger type
+ Printer	
Driver	
Paper	
+ Options	
Factory setting	

--		
--		
<b>Non-circular</b>   Circular		
Shot		
<b>Interval</b> 1... 12:59:59   <b>00:02:00</b>		
Difference	1st difference	On   Off
	Delta cond	0 ... 2000 mS/cm   <b>1.0 µS/cm</b>
	Delta conc	0 ... 9.99 %   <b>1.0 %</b>
	Delta MΩcm	0 ... 9.999 MΩcm   <b>1 MΩcm</b>
	Delta salinity	0 ... 45.0 g/kg   <b>1.0 g/kg</b>
	Delta TDS	0 ... 2000.0 mg/l   <b>1 mg/l</b>
	2nd diff.	On   <b>Off</b>
	Delta °C	0...99.9 °C   <b>1.0 °C</b>
	Delta °F	0...450 °F   <b>1.0 °F</b>
Intv+Diff	Interval	See logger type: Interval
	Difference	See logger type: Difference
Limit value	Interval	Basis/Event 00.00.01... 12:59:59 <b>00:01:00/00:00:01</b>
	Limit values	Min/Max, corresponding to permissible range (see Specifications)

<b>Auto</b>   HP PCL   Epson   Samsung   IBM	
<b>A4</b>   Small	
001 SOP	Add-on function, enabled via TAN
002 Temp.cal	
Yes   <b>No</b>	
<b>Note:</b> Reset to factory settings will also erase all logger data!	

\* “+” indicates that submenus can be opened by pressing **enter**.

## Oxygen Configuration

- 1) In measuring mode, press the **Menu** softkey.
- 2) Select "Configuration" and confirm by pressing **enter**.
- 3) Make the desired adjustments.

The following table gives you an overview.

Factory settings are shown in **bold print**.

"Oxygen Configuration" menu selection – part 1

	Language	<b>Deutsch</b>   English   Español   Italiano   Français   Português   Русский
	Auto-off	<b>Off</b>   5 min   10 min   30 min   60 min
	Temperature	°C   °F
	Right softkey	Logger Start/Stop   <b>Hold value</b>   Print
	+ Oxy sensor*	
	Display	<b>Saturation</b>   Concentration   Partial pressure
	Salinity	<b>0</b> ... 45.0 g/kg
	Pressure correction**	<b>Air pressure</b>   Manual
	Pressure	Manual: 0 ... 9999 mbar   <b>1013 mbar</b>
	Wear	<b>On</b>   Off
	+ Calibration*	
	Cal mode	<b>In air</b>   Zero point   Data entry
	Calibration timer	<b>Off</b>   On
	Interval	On: 0 ... 99 days
	MemoLog	<b>Off</b>   On
	TAG	<b>Off</b>   On
	+ Time/Date*	
	Time format	<b>24 h</b>   12 h
	Date format	<b>dd.mm.yyyy</b>   yyyy-mm-dd   dd/mm/yyyy   mm/dd/yyyy
	Time	hh:mm:ss
Date	Date format as configured	

\* "+" indicates that submenus can be opened by pressing **enter**.

\*\* The device provides an internal barometer.

Oxy

"Oxygen Configuration" menu selection – part 2

+ Display*
Appearance
Lighting
Brightness
+ Data logger*
Meas.point
Note
Right softkey
Recording
Logger type
+ Printer
Driver
Paper
+ Options
Factory setting

<b>Modern</b>   Retro
<b>Permanent</b>   60 min   30 min   10 min   5 min   1 min   30 sec
<b>Bright</b>   Standard   Dim

--		
--		
Logger Start/Stop	<b>Hold value</b>	
	<b>Non-circular</b>   Circular	
Shot		
<b>Interval</b>	00.00.01...12:59:59   <b>00:02:00</b>	
Difference	1st difference	<b>On</b>   Off
	Delta saturation	0 ... 200 %Air   <b>1% Air</b>
	Delta conc	0 ... 20 mg/l   <b>1 mg/l</b>
	Delta mbar	0 ... 999.99 mbar   <b>1 mbar</b>
	2nd diff.	On   <b>Off</b>
Delta °C	0...99.9 °C   <b>1.0 °C</b>	
	Delta °F   <b>1.0 °F</b>	
Intv+Diff	Interval	See logger type: Interval
	Difference	See logger type: Difference
Limit value	Interval	Basis/Event 00.00.01...12:59:59 <b>00:01:00/00:00:01</b>
	Limit values	Min/Max, corresponding to permissible range (see Specifications)

<b>Auto</b>   HP PCL   Epson   Samsung   IBM	
<b>A4</b>   Small	
001 SOP	Add-on function, enabled via TAN
002 Temp.cal	
Yes   <b>No</b>	
<b>Note:</b> Reset to factory settings will also erase all logger data!	

\* "+" indicates that submenus can be opened by pressing **enter**.

## pH Calibration

- 1) In measuring mode, press the **Menu** softkey.
- 2) Select "Calibration" and confirm by pressing **enter**.
- 3) Select the desired "Calibration mode" and confirm by pressing **enter**.
- 4) The "TAG" menu allows the sensor TAG to be edited.  
To do so, set "TAG" to **On** in the configuration menu (default setting: **Off**).
- 5) Perform the selected calibration as described on the following pages.  
Follow the instructions on the display.

**Note:** Calibration is not possible when the device is connected via USB with the Paraly SW 112 software.

## Calimatic Calibration

**(Automatic calibration with specification of the buffer solution used)**

- 1) Select the number of calibration points and the buffer set as shown in the table below and press the **Start** softkey.

Cal points	<b>Auto</b>	1-point	2-point	3-point
Buffer set	Mettler-Toledo		2.00/4.01/7.00/9.21	
	<b>Knick CaliMat</b>		2.00/4.00/7.00/9.00/12.00	
	Ciba		2.06/4.00/7.00/10.00	
	NIST technical		1.68/4.00/7.00/10.01/12.46	
	NIST standard		1.679/4.006/6.865/9.180	
	Hach		4.01/7.00/10.01/12.00	
	WTW		2.00/4.01/7.00/10.00	
	Hamilton		2.00/4.01/7.00/10.01/12.00	
	Reagecon		2.00/4.00/7.00/9.00/12.00	
	DIN 19267		1.09/4.65/6.79/9.23/12.75	
	Metrohm		4.00/7.00/9.00	
	User buffer 1		configurable via Paraly SW 112 software	

- 2) Immerse the sensor in the **1st/2nd/3rd** buffer solution and press **Continue** (repeat this step for each calibration point).
- 3) Finally, the calibration data will be displayed.  
You can **Apply** or **Discard** these values.

**Note:** To abort calibration, you can press **meas** at any time.

pH

## Manual Calibration

(Calibration with manual specification of the number of calibration points and the buffer solution)

- 1) Select the number of calibration points and press the **Start** softkey.
- 2) Adjust the temperature-corrected value (see buffer table) for the **1st**/2nd/3rd buffer solution and press **Continue** (repeat this step for each calibration point).

**Note:** When using sensors without temperature detector, you should adjust the temperature manually before starting calibration (see page 31).

- 3) Finally, the calibration data will be displayed.

You can **Apply** or **Discard** these values.

## Data Entry Calibration

(Calibration by entering known sensor values)

- 1) Press the **Start** softkey.
- 2) Enter the known sensor values for zero and slope.
- 3) Finally, you can **Apply** these values or **Cancel** the calibration.

pH

ORP

## Calibrating a pH/ORP Combo Sensor

The pH/ORP combo sensor can be calibrated as a pH sensor and/or as an ORP sensor.

### pH Calibration

Follow the instructions given for pH calibration, page 23.

### ORP Calibration

Follow the instructions given for ORP calibration, page 25.

**Note:** To abort calibration, you can press **meas** at any time.

## ORP Calibration

- 1) In measuring mode, press the **Menu** softkey.
- 2) Select "Calibration" and confirm by pressing **enter**.
- 3) Select the desired "Calibration mode" and confirm by pressing **enter**.
- 4) The "TAG" menu allows the sensor TAG to be edited.  
To do so, set "TAG" to **On** in the configuration menu (default setting: **Off**).
- 5) Enter the temperature-corrected setpoint of the calibration solution.
- 6) Immerse the sensor in the calibration solution and wait until the measured value is stable.
- 7) **Apply** or **Discard** the ORP setpoint.

**Note:** Calibration is not possible when the device is connected via USB with the Paraly SW 112 software.

**Note:** To abort calibration, you can press **meas** at any time.

## ISFET

## ISFET Calibration

- 1) In measuring mode, press the **Menu** softkey.
- 2) Select "Calibration" and confirm by pressing **enter**.
- 3) Select the desired "Calibration mode" and confirm by pressing **enter**.
- 4) The "TAG" menu allows the sensor TAG to be edited.  
To do so, set "TAG" to **On** in the configuration menu (default setting: **Off**).
- 5) Perform the selected calibration as described on the following pages.  
Follow the instructions on the display.

**Note:** Calibration is not possible when the device is connected via USB with the Paraly SW 112 software.

## Calibrating the ISFET Zero (Operating Point)

- 1) Select the "ISFET zero" calibration mode for setting the operating point for the first sensor calibration.

Calibration mode	<b>Calimatic</b>
	Manual
	Data entry
	ISFET zero (operating point)

- 2) Press the **Start** softkey.
- 3) Adjust the buffer value if required: default pH 7.00
- 4) Press the **Start** softkey.
- 5) Finally, you can **Apply** or **Discard** the calibration value for the operating point.  
When you apply the calibration value, the operating point will be stored in the device, but not in the sensor!  
Keep the sensor connected to the Portavo while performing the next calibration step. The operating point will be taken into account for the following calibration.

## Calimatic/Manual/Data Entry Calibration

Follow the instructions given for pH calibration, page 23.

If you disconnect the sensor before performing the calibration (e.g., Calimatic), you must set the operating point again as described above.

**Note:** To abort calibration, you can press **meas** at any time.

## Conductivity Calibration

- 1) In measuring mode, press the **Menu** softkey.
- 2) Select "Calibration" and confirm by pressing **enter**.
- 3) Select the desired "Calibration mode" and confirm by pressing **enter**.
- 4) The "TAG" menu allows the sensor TAG to be edited.  
To do so, set "TAG" to **On** in the configuration menu (default setting: **Off**).
- 5) Perform the selected calibration as described on the following pages.  
Follow the instructions on the display.

**Note:** Calibration is not possible when the device is connected via USB with the Paraly SW 112 software.

## Auto Calibration

**(Automatic calibration with specification of the calibration solution used)**

---

### NOTICE!

- Make sure that the values of the calibration solutions used correspond exactly to those specified in this manual.  
If not, the resulting cell constant will be incorrect.
  - When calibrating in a liquid, make sure that the sensor, the separate temperature probe (if present) and the calibration solution have the same temperature.  
Only this ensures that the cell constant is determined correctly.
- 

- 1) Select the calibration solution:
  - **NaCl 0.01 mol/l**
  - NaCl 0.1 mol/l
  - NaCl sat.
  - KCl 0.01 mol/l
  - KCl 0.1 mol/l
  - KCl 1 mol/l
- 2) Press the **Start** softkey.
- 3) Immerse the sensor in the solution and press **Continue**.
- 4) Finally, the calibration data record will be displayed.  
You can **Apply** or **Discard** these values.

**Note:** To abort calibration, you can press **meas** at any time.

## Cond

### “Entry of Solution” Calibration

(Calibration by entering the conductivity with display of the cell constant)

- 1) Press the **Start** softkey.
- 2) Immerse the sensor in the solution.
- 3) Enter the temperature-corrected conductivity value and press **enter**.
- 4) Finally, you can **Apply** these values or **Cancel** the calibration.

### Cell Constant / Cell Factor Calibration

(Calibration by entering the cell constant (cell factor) with display of the conductivity)

- 1) Press the **Start** softkey.
- 2) Immerse the sensor in the solution.
- 3) Modify the value of the cell factor (cell constant) until the temperature-corrected conductivity value is reached. Then press **enter**.
- 4) Finally, you can **Apply** these values or **Cancel** the calibration.

<b>Contacting conductivity sensor (conductive)</b>	<b>Cell constant</b>
SE 215 MS	1.00/cm $\pm$ 2%
<b>Toroidal conductivity sensor (inductive)</b>	<b>Cell factor</b>
SE 680 MS	6.4/cm

### Installation Factor Calibration

- 1) Make sure that the sensor is in normal mounting position in the medium.
- 2) Press the **Start** softkey.
- 3) Modify the installation factor until the correct conductivity value is displayed (reference measurement). Then press **enter**.
- 4) Finally, you can **Apply** these values or **Cancel** the calibration.

### Zero Calibration

- 1) Make sure that the sensor is outside the medium (in air).
- 2) Press the **Start** softkey.
- 3) Finally, you can **Apply** these values or **Cancel** the calibration.

**Note:** To abort calibration, you can press **meas** at any time.

## Oxygen Calibration

- 1) In measuring mode, press the **Menu** softkey.
- 2) Select "Calibration" and confirm by pressing **enter**.
- 3) Select the desired "Calibration mode" and confirm by pressing **enter**.
- 4) Select "Membrane module replacement" if you wish to save a change of membrane or electrolyte in the connected sensor. The digital, optical oxygen sensor automatically recognizes when its sensor cap has been replaced.
- 5) The "TAG" menu allows the sensor TAG to be edited.  
To do so, set "TAG" to **On** in the configuration menu (default setting: **Off**).
- 6) Perform the selected calibration as described on the following pages.  
Follow the instructions on the display.

**Note:** Calibration is not possible when the device is connected via USB with the Paraly SW 112 software.

## Calibration in Air

### (Calibrating the slope in air)

- 1) Place sensor in air and wait for a stable measured value.
- 2) Press the **Start** softkey.
- 3) Adjust the correct value for "Relative humidity". Then press **Continue**.  
Calibration will be performed.
- 4) Finally, you can **Apply** or **Discard** these values.

**Note:** To abort calibration, you can press **meas** at any time.

## Oxy

## Zero Calibration

(Zero calibration with oxygen-free medium, e.g., nitrogen 5.0)

- 1) Place sensor in oxygen-free medium and wait for a stable measured value.
- 2) Press the **Start** softkey. Calibration will be performed.
- 3) Finally, you can **Apply** these values or **Cancel** the calibration.

## Data Entry Calibration

(Calibration by entering known sensor values)

- 1) Press the **Start** softkey.
- 2) Adjust the known sensor values for zero and slope.
- 3) Finally, you can **Apply** these values or **Cancel** the calibration.

**Note:** To abort calibration, you can press **meas** at any time.

pH

ORP

Oxy

Cond

Once you have completed all preparations, you can start with the actual measurement.

- 1) Connect the desired sensor to the meter. Some sensors require a special preparation. Please proceed according to the operating instructions for the sensor.
- 2) Switch the meter on using the **on/off** or **meas** key.
- 3) Depending on the measurement method and the sensor used, immerse the sensing part of the sensor in the medium to be measured.
- 4) Watch the display and wait for the reading to stabilize.

**Note:** Measurement can also be controlled via the Paraly SW 112 software.

## Toggling the Measured Value Display

During measurement, you can toggle between display of primary / secondary measured values and clock by pressing **meas**.

## Adjusting the Temperature

When you connect a sensor without temperature detector, you can manually adjust the temperature for measurement or calibration:

- 1) Press **meas** to access measuring mode. The adjusted temperature will be displayed.
- 2) Set the desired temperature value using the ▼ or ▲ arrow. Holding the key depressed changes the temperature value at high speed.

pH

ORP

Oxy

Cond

## Printing

You can print out the currently measured values and the calibration record. Via the micro USB port, you can connect the Portavo 908 to a suitable printer. To do so, use the included USB adapter (A female to B male). The following printer protocols are supported: HP-PCL, Epson, Samsung, and IBM.

### Required Settings in the “Configuration” Menu

(depending on the measured variable;  
pH: page 14, Cond: page 18, Oxy: page 21)

- Select “Right softkey: Print”.
- From the “Printer” submenu, select the driver for your printer and the paper size:

+ Printer	
Driver	<b>Auto</b>   HP PCL   Epson   Samsung   IBM
Paper	<b>A4</b>   Small

## Printing Currently Measured Values

- 1) Use the included USB adapter (A female to B male) to connect a printer to the micro USB port of the meter and switch it on.
- 2) Switch the meter on using the **on/off** or **meas** key.

Press the right softkey: **Print**.

The currently measured values will be printed.

#### Note:

A4 printers do not eject the paper immediately.

The printer collects several values until it prints the page. If you want to eject the page before the print job is finished, you can press and hold the **Print** softkey (for approx. 2 seconds).

- 3) You will receive a printout of the currently measured values with date and time stamp (example):

```

908 MULTI          0003792
SE 533X/1-NMSN    2180694
26.11.2015 12:26   pH 7.137   3.8 mV   23.0 °C
  
```

pH

ORP

Oxy

Cond

## Printing a Calibration Record

- 1) Use the included USB adapter (A female to B male) to connect a printer to the micro USB port of the meter and switch it on.
- 2) Switch the meter on using the **on/off** or **meas** key.
- 3) Open the "Info" menu.
- 4) Select "Calibration record" (must have been created before).
- 5) Press the right softkey: **Print** (see previous page for configuration).

You will receive a printout of the calibration record with date and time stamp (example):

### Calibration record

```

Date                04.01.2016 09:19

Device information
Manufacturer        Knick
Serial no.          0003791
Device type         908MULTI
SW version          1.5.0 Build10931

Sensor information  Memosens
Sensor model        pH
Manufacturer        Knick
Order no.           SE515/1-MS
Serial no.          10000
SW version          1.0.6
HW version          0.5.2
TAG
Temp offset         0.0 K
Operating time      0 h
Wear                0 %
SIP                 0

Calibration data
Date of calibration 03.01.12 12:00
Zero point          pH 7.000 | 0.0 mV
Slope               100.000 % | 59.2 mV
Buffer 1            pH 4.010
Buffer 2            pH 9.210
    
```

pH

ORP

Oxy

Cond

## The Data Logger


The meter provides a data logger. **Prior to use**, it must be configured and then activated. You can choose from the following logger types:

- Shot (manual logging by pressing the **Save value** softkey)
- Interval (time-controlled logging at a fixed interval)
- Difference (signal-controlled logging of measured variable and temperature)
- Intv+Diff (combined time- and signal-controlled logging)
- Limit value (combined time- and threshold-controlled logging)

The data logger records up to 10,000 entries, which can be assigned to different points of measurement (TAGs) and notes. The following data will be recorded: meas. point, note, sensor ID, serial number of sensor (Memosens), primary value, temperature, time stamp, device status.

**It is always the currently selected process variable which is recorded.**

### Display: Icons related to the data logger

pH

ORP

Oxy

Cond


## Operating Modes of the Data Logger (Logger Type)

### Shot

In this mode, a measured value is recorded each time the **Save value** softkey is pressed. In the measuring mode (**meas**), it is always possible to hold a value and then save it.


### Interval (Time-Controlled)

In the "Interval" mode, the data are cyclically recorded.

### Difference

When the delta range (process variable and/or temperature) related to the last entry is exceeded, a new entry is created and the delta range is displaced upwards or downwards by the delta value. The first entry is automatically created when the data logger is started.

pH


ORP

Oxy


Cond

**Interval and Difference (Combined)**

When the delta range related to the last DIFF entry is exceeded, a new entry is created (example: entry **A**) and the delta range is displaced upwards or downwards by the delta value. As long as the measured value remains within the delta range, logging is performed at the preset interval. The first DIFF entry is automatically created when the data logger is started.

**Limit Value (Combined)**

When one of the two limit values (Min/Max) is exceeded, the data are logged as defined by the "event interval". Additionally, the last ten measured values before an event are recorded (pre-trigger). As long as the measured value remains within the limits, logging is performed at the preset "basic interval".

pH

ORP

Oxy

Cond

## Configuring the Data Logger

Prerequisite: Data logger is stopped.

The “Data logger” menu shows the number of occupied entries as well as the number of free entries. Configuration can also be done in the “Configuration” menu under “Data logger”.

- 1) Press **Menu** softkey.
- 2) Select “Data logger” and confirm by pressing **enter**.
- 3) Select “Configure data logger” and confirm by pressing **enter**.
- 4) Configure data logger as required (see table).
- 5) When you have completed the configuration, you can start the data logger!

## Increasing the Battery Life

To increase the battery life for logger operation, the time for the display lighting selected in the configuration should be as short as possible.

**Note:** When the selected time has expired, display and backlighting switch off automatically. They can be switched on again by pressing any key.

pH

ORP

Oxy

Cond

**Configuring the data logger (default in bold print)**

Meas.point	<b>Without</b>			
Note	<b>Without</b>			
Recording	<b>Non-circular</b>			
	Circular			
Logger type	Shot			
	<b>Interval</b>	Interval	00:00:01...12:59:59   <b>00:02:00</b>	
	Difference	1st difference	<b>On Off</b>	
		Delta pH	pH 0.0...16.0   <b>pH 1.0</b>	
		Delta mV	0 ... 2000 mV   <b>1 mV</b>	
		Delta cond	0 ... 2000 mS/cm   <b>1.0 µS/cm</b>	
		Delta conc	0 ... 9.99 % by wt   <b>1%</b>	
		Delta MΩcm	0 ... 9.999 MΩcm   <b>1.0 MΩcm</b>	
		Delta salinity	0.00 ... 45.0 g/kg   <b>1.0 g/kg</b>	
		Delta TDS	0.00 ... 2000.0 mg/l   <b>1 mg/l</b>	
		Delta saturation	0 ... 200% Air   <b>1% Air</b>	
		Delta conc	0 ... 20.0 mg/l   <b>1 mg/l</b>	
		Delta mbar	0 ... 1000 mbar   <b>1 mbar</b>	
		2nd difference	<b>On Off</b>	
		Delta °C	0...99.9 °C   <b>1.0 °C</b>	
	Delta °F	0...450 °F   <b>1.0 °F</b>		
	Intv+Diff	Interval	See logger type: Interval	
		Difference	See logger type: Difference	
	Limit value	Interval	Basis	00.00.01...12:59:59   <b>12:01:00 AM</b>
			Event	<b>00.00.01</b> ...12:59:59
Limit values		Min/Max corresponding to permissible range (see Specifications)		

pH

ORP

Oxy

Cond

## Starting/Stopping the Data Logger

With the data logger activated, automatic switch-off is disabled. Every time the meter has been switched off, the data logger must be restarted.

Depending on the assignment of the right softkey (see "Configuration", page 14) you can start/stop the data logger as follows:

Right softkey	
Start/Stop logger	1) Press right softkey <b>Start logger / Stop logger</b> .
Hold value	1) Press <b>Menu</b> softkey. 2) Select "Data logger" using the arrows and confirm by pressing <b>enter</b> . 3) Press <b>Start</b> or <b>Stop</b> softkey, resp.

## Viewing the Logger Data

In the "Data logger" menu you can view the recorded entries either individually or as curve characteristic (see examples).

You can also use the Paraly SW 112 software for managing the data logger.

- 1) Press **Menu** softkey.
- 2) Select "Data logger" using the arrows and confirm by pressing **enter**.
- 3) Select "View logger data" using the arrows and confirm by pressing **enter**.
- 4) Select filter ("Meas.point" or "Time + Meas.point" or "All values").
- 5) Select the parameter corresponding to the sensor.
- 6) Press **Menu** softkey.
- 7) Select the desired entries using the arrow (see example 1).
- 8) For display as curve characteristic, press **Graphic** softkey.

You can use the arrows to navigate between entries (see example 2).

## Deleting the Logger Data

To delete the recorded entries, proceed as follows:

- 1) Press **Menu** softkey.
- 2) Select "Data logger" using the arrows and confirm by pressing **enter**.
- 3) Select "Delete logger data" using the arrows and confirm by pressing **enter**.
- 4) Select deletion mode: "Complete", "Data", "Meas.point" or "Filter"  
(you can filter for measuring point, parameter or time).
- 5) Press **Delete** softkey. The data are deleted according to the configuration.
- 6) Press **Back** softkey to return to menu selection.

Example 1: Viewing the logger data

**Data Logger**

Primary measured value → **pH 7.000**    **25.0 °C** ← Measured temperature

☺ **15.08.2012 10:45:36** ← Date and time

Meas.point    ABC

Note

Sensor        Analog

Scroll using arrow keys → **1/3**

**Return**    **Graphic**

Example 2: Curve characteristic

**Data Logger**

7.50    25.8

pH    °C

6.50    24.3

Entry 2 of 3 → **2/3**

Select entry using arrows. →

**Return**    **Measured values**    **Date**    **Time**

pH 7.21    25.0 °C    24/08/2012    17:12:50


pH

ORP


Oxy

Cond

Error and status messages appear as plain text on the display. By pressing **enter** and **Help**, more detailed help texts can be displayed. Information on the sensor condition is indicated by the “Senseface” icon (friendly, neutral, sad) possibly accompanied by an info text.

Example of an error message: Press **enter** and **Help** to access the help text.

Help text for error 21

Senseface (the “smiley” icon) provides information on the sensor condition (maintenance request). Measurement can still be performed. After a calibration, the corresponding Senseface icon (friendly, neutral, sad) is shown together with the calibration data. Otherwise, Senseface is only visible in measuring mode.

pH

ORP

Oxy

Cond

## “Sensoface” Messages

The “Sensoface” icon provides information on the sensor condition:

### Sensoface Meaning

Sensor is okay

Calibrate the sensor soon

Calibrate or replace the sensor

## Info and Help Texts

When an error or status message appears on the screen, proceed as follows to view the corresponding info or help text:

- 1) Press **enter**.
- 2) Press the **Help** softkey.
- 3) The help text will be displayed. In most cases, you can remedy the cause of the error by yourself. Please refer to the following table for possible remedies.

Info	Message
Info 01	Cal timer expired
Info 02	Sensor wear
Info 03	Bad glass impedence
Info 05	Zero/Slope
Info 06	Response time too long
Info 07	Operating point (ISFET)
Info 08	Leakage current (ISFET)
Info 09	ORP offset
Info 10	Polarization

pH

ORP

Oxy

Cond

## Error Messages

Error	Message	Remedy
 blinks	Replace the batteries	Replace the batteries.
ERR 1	Primary variable range	Check whether the measurement conditions correspond to the adjusted measuring range.
ERR 2	ORP range	
ERR 3	Temperature range	
ERR 4	Zero point	Thoroughly rinse the sensor and recalibrate.
ERR 5	Slope	If this does not help, replace the sensor.
ERR 6	Cell constant too high/ too low	Enter nominal cell constant or calibrate the sensor using a known solution.
ERR 7	Air pressure range	Check if the opening for the pressure sensor located on the back of the device is blocked.
ERR 8	Identical buffers!	Use a buffer solution with a different nominal value before starting the next calibration step.
ERR 10	Buffers interchanged!	Repeat calibration.
ERR 11	Unstable value (Drift too high)	Leave the sensor in the liquid until the measured value is stable. If this does not help, replace the sensor.
ERR 14	Time and date invalid	Set the date and time.
ERR 18	System error	Restart, reset to factory settings, configure and calibrate. If the error occurs again, contact the Service.
ERR 19	Factory settings error	Data error, measurement with analog sensors no longer possible. Contact the Service team.
ERR 21	No sensor connected	Connect operational Memosens sensor.
ERR 25	Buffer distance	Re-enter the buffer table (Paraly SW 112).
ERR 30	Data logger full	Clear the logger completely or partially.
ERR 31	MemoLog full	Clear the MemoLog completely or partially.

pH

ORP

Oxy

Cond

## Option 001 SOP

The option meets specific requirements of the pharmaceutical and biotechnological industries.

### SOP Cal

Here, you specify which buffers are to be used in which sequence.

You can combine buffer solutions from different buffer sets. Please note that the minimum distance allowed between two buffer solutions is  $\Delta 2$  pH.

### User Management (Access Control)

You can create up to 4 users with different access rights for configuration or calibration.

### Sensor Verification

To make sure that only selected sensors can be operated on the meter, you can evaluate the sensor type and/or the "TAG" and "Group" data stored in the sensor.

The sensor will only be accepted if the data stored in the sensor corresponds to the data stored in the meter.

### Temperature Adjustment

For Memosens sensors, you can perform a 1-point calibration of the internal temperature detector.

## Option 002 Temp.Cal

(included in Option 001 SOP)

### Temperature Adjustment

For Memosens sensors, you can perform a 1-point calibration of the internal temperature detector.

## Enabling the Option 001 SOP

- 1) In measuring mode, press the **Menu** softkey.
- 2) Select "Configuration" and confirm by pressing **enter**.
- 3) Select Option "001 SOP" and enter your activation code.

## Adapt SOP Cal

Here, you specify which buffers are to be used in which sequence.

You can combine buffer solutions from different buffer sets. Please note that the minimum distance allowed between two buffer solutions is  $\Delta 2$  pH.

The "Configuration - Calibration" menu is extended as follows:


+ Calibration	
Cal mode	Calimatic   Manual   Data entry   <b>SOP cal</b>
Adapt SOP cal	
Cal points	1-point   2-point   3-point
Buffer 1	
Buffer set	Mettler-Toledo 2.00/4.01/7.00/9.21 <b>Knick CaliMat</b> 2.00/4.00/7.00/9.00/12.00 Ciba 2.06/4.00/7.00/10.00 NIST technical 1.68/4.00/7.00/10.01/12.46 NIST standard 1.679/4.006/6.865/9.180 Hach 4.01/7.00/10.01/12.00 WTW 2.00/4.01/7.00/10.00 Hamilton 2.00/4.01/7.00/10.01/12.00 Reagecon 2.00/4.00/7.00/9.00/12.00 DIN 19267 1.09/4.65/6.79/9.23/12.75 Metrohm 4.00/7.00/9.00 User buffer 1
Buffer	Select a buffer from the selected set
Buffer 2	Select buffer set 2 and buffer (see buffer 1)
Buffer 3	Select buffer set 3 and buffer (see buffer 1)
Verification	<b>Off</b>   On
Delta pH	<b>pH 0.05</b> (Enter maximum permitted deviation from verification buffer; exceeding this value generates an error message)
Verification buffer	Select buffer set and buffer (see buffer 1)

## Selecting SOP Calibration

- 1) In measuring mode, press the **Menu** softkey.
- 2) Select "Calibration" and confirm by pressing **enter**.
- 3) Select the "SOP cal" calibration mode and confirm by pressing **enter**.

## Performing an SOP Calibration

The sequence of buffers to be used is displayed as specified in the configuration. After each calibration step, the identified buffer is marked off in the display. The next operation procedure is displayed. Perform the calibration following the instructions given in the display.

**Note:** Calibration is not possible when the device is connected via USB with the Paraly SW 112 software.

pH

ORP

Oxy

Cond

## User Management (Access Control)

You can create up to 4 users with different access rights for configuration or calibration.

### Enabling the User Management

- 1) In measuring mode, press the **Menu** softkey.
- 2) Select "Configuration" and confirm by pressing **enter**.
- 3) Select "User management"/"Management – Enable".
- 4) Select  
User = ADMIN  
PIN code = 1989 (factory setting)
- 5) Press **enter**.
- 6) To create more users / assign more PIN codes: Press **Continue** softkey.

### Creating a User / Changing the PIN Code

1. Select a user (e.g., "User 1", default: ADMIN, PIN code 1989):

The first screenshot shows the 'Configuration' menu with 'Management - Disable' at the top. Below it are four options: '+ User 1', '+ User 2', '+ User 3', and '+ User 4'. At the bottom are 'Return' and 'Continue' buttons. An arrow labeled 'enter' points to the second screenshot. The second screenshot shows the same menu but with '- User 1' selected. Below it are fields for 'ADMIN', 'PIN code 1989', 'cal level Access', and 'conf level Access'. At the bottom are 'Return' and 'Continue' buttons.

You can allow or block access to configuration or calibration for each user.

2. Select ADMIN to open the editor for entering the user name:

The first screenshot shows the 'Configuration' menu with 'Management - Disable' at the top. Below it are '- User 1' and 'ADMIN' (selected). Below that are fields for 'PIN code 1989', 'cal level Access', and 'conf level Access'. At the bottom are 'Return' and 'Continue' buttons. An arrow labeled 'enter' points to the second screenshot. The second screenshot shows a text editor with 'End ADMIN |' at the top. Below it is a keyboard layout with letters and symbols. At the bottom are 'abc' and 'Finish' buttons.

### NOTICE!

If you lose the PIN code for the ADMIN user, the system access will be blocked. The manufacturer can generate a rescue PIN code.

pH

ORP

Oxy

Cond

## How to Enter the Rescue PIN Code

- 1) In measuring mode, press the **Menu** softkey.
- 2) Select "Configuration" using the cursor keys.
- 3) Press the ▼ and ▲ keys simultaneously.
- 4) Set User to "ADMIN".
- 5) Select "PIN code" and enter the rescue PIN code.  
Confirm by pressing **enter**.
- 6) Press the **Continue** softkey.

## Sensor Verification

To make sure that only selected sensors can be operated on the meter, you can evaluate the following data stored in the sensor:

- Model (sensor model)
- TAG (e.g., point of measurement)
- Group (e.g., facility)

With Option 001 enabled, the "Configuration" menu is extended as follows:

- Sensor verification	
Check model	<b>Off</b>   Info   Reject
Check TAG	<b>Off</b>   Info   Reject
Check group	<b>Off</b>   Info   Reject

You can select the following options:

- Off**      No verification.
- Info**     When a wrong sensor is connected, an error message will be displayed.  
Nevertheless, you can continue working with the sensor.
- Reject**   Here you specify values with which the sensor will be rejected.

pH

ORP

Oxy

Cond

## Enabling the Option 002 Temp.Cal (included in Option 001 SOP)

- 1) In measuring mode, press the **Menu** softkey.
- 2) Select "Configuration" and confirm by pressing **enter**.
- 3) Select Option "002 Temp.cal" and enter your activation code.

## Selecting Temperature Calibration

- 1) In measuring mode, press the **Menu** softkey.
- 2) Select "Calibration" and confirm by pressing **enter**.
- 3) Select the "Temperature" calibration mode and confirm by pressing **enter**.

For Memosens sensors, you can perform a 1-point calibration of the internal temperature detector. To do so, enter the reference temperature and confirm the temperature adjustment by pressing the **Apply** softkey:

The screenshot shows a 'Calibration' screen with a red header. The current temperature is 24.4 °C. Below this, it prompts the user to 'Enter the reference temperature.' The 'Temp. offset' is set to -0.3 K. The 'Reference temperature' is set to 24.1 °C. At the bottom, there are 'Cancel' and 'Apply' buttons.

<b>Calibration</b>	
24.4 °C	
Enter the reference temperature.	
Temp. offset	-0.3 K
Reference temperature	24.1 °C
Cancel	Apply

pH

ORP

Oxy

Cond

<b>Connections</b>	2x socket, 4 mm dia., for separate temp. detector 1x M8 socket, 4 pins, for flexible Memosens lab cable 1x micro USB-B for data transmission to PC or printer connection 1x M12 socket, 8 pins, for flexible connecting cable for Memosens sensors or SE 340 sensor (optical oxygen)	
<b>Air pressure measurement</b>	700 ... 1100 hPa	
<b>User interface</b>	Straightforward menu navigation with graphic icons and detailed operating instructions in plain text	
Languages	German, English, French, Spanish, Italian, Portuguese, Russian	
Sensoface	Status indication (friendly, neutral, sad)	
Status indicators	For battery power level, logger	
Graphic display	QVGA TFT display with white backlighting	
Keypad	[on/off], [meas], [enter], [◀], [▶], [▲], [▼], 2 context-sensitive softkeys	
<b>Data logger</b>	10 000 memory locations	
Recording	Manual, interval- or event-controlled, with management of tag numbers and notes	
<b>Cal data logger MemoLog (Memosens only)</b>	Up to 100 Memosens calibration records can be saved	
	Recording	Directly retrievable via MemoSuite or Paraly SW 112 (USB)
	Viewable on the display	Manufacturer, sensor type, serial no., zero, slope, calibration date
<b>Temperature input</b>	2 x 4 mm dia. for integrated or separate temperature detector	
Measuring ranges	NTC30 temperature detector	-20 ... +120 °C / -4 ... +248 °F
	Pt1000 temperature detector	-40 ... +250 °C / -40 ... +482 °F
Measuring cycle	Approx. 1 s	
Measurement error <sup>1,2,3)</sup>	< 0.2 K (Tamb = +23 °C / +73.4 °F); TC < 25 ppm/K	

1) according to EN 60746-1, at nominal operating conditions

2) ± 1 count

3) plus sensor error

pH

ORP

Oxy

Cond

<b>Communication</b>	USB 2.0
Profile	HID, driverless installation
Usage	Data exchange and configuration via Paraly SW 112 software

## Diagnostics functions

Sensor data (Memosens only)	Manufacturer, sensor type, serial number, wear, operating time, remaining lifetime, max. temperature, adaptive calibration timer, calibration and adjustment data, SIP, CIP and autoclaving counter
Calibration data	Calibration date; pH/Oxy: zero, slope; Cond: Cell constant
Device self-test	Automatic memory test (FLASH, EEPROM, RAM)
Device data	Device type, software version, hardware version

<b>Data retention</b>	Parameters, calibration data > 10 years
-----------------------	---

<b>EMC</b>	EN 61326-1 (General Requirements)
Emitted interference	Class B (residential environment)
Immunity to interference	Industrial environment EN 61326-2-3 (Particular Requirements for Transmitters)

<b>RoHS conformity</b>	According to directive 2011/65/EU
------------------------	-----------------------------------

<b>Power supply</b>	4 x AA alkaline batteries or 1x Li-ion battery, USB chargeable
---------------------	---

## Nominal operating conditions

Ambient temperature	-10 ... +55 °C / 14 ... +131 °F
Transport/ Storage temperature	-25 ... +70 °C / -13 ... +158 °F
Relative humidity	0 ... 95 %, short-term condensing allowed

## Housing

Material	PA12 GF30 (silver gray RAL 7001) + TPE (black)
Protection	IP66/67 with pressure compensation
Dimensions	Approx. (132 x 156 x 30) mm
Weight	Approx. 500 g

<b>Printer</b>	Printer protocols: HP-PCL, Epson, Samsung, IBM (ASCII texts) Connection via standard USB cable and USB adapter (A female to B male)
----------------	--

pH

ORP

<b>Memosens pH input (also ISFET)</b>	M8 socket, 4 pins, for flexible Memosens lab cable or M12 socket, 8 pins, for flexible connecting cable for Memosens sensors	
Display ranges <sup>4)</sup>	pH	-2.00 ... +16.00
	mV	-1999 ... +1999 mV
	Temperature	-50 ... +250 °C / -58 ... +482 °F
<b>Sensor standardization *</b>	pH calibration	
Operating modes *	Calimatic	Calibration with automatic buffer recognition
	Manual	Manual calibration with entry of individual buffer values
	Data entry	Data entry of zero and slope
Calimatic buffer sets *	-01- Mettler-Toledo	2.00/4.01/7.00/9.21
	-02- Knick CaliMat	2.00/4.00/7.00/9.00/12.00
	-03- Ciba (94)	2.06/4.00/7.00/10.00
	-04- NIST technical	1.68/4.00/7.00/10.01/12.46
	-05- NIST standard	1.679/4.006/6.865/9.180
	-06- HACH	4.01/7.00/10.01/12.00
	-07- WTW techn. buffers	2.00/4.01/7.00/10.00
	-08- Hamilton	2.00/4.01/7.00/10.01/12.00
	-09- Reagecon	2.00/4.00/7.00/9.00/12.00
	-10- DIN 19267	1.09/4.65/6.79/9.23/12.75
	-11- Metrohm	4.00/7.00/9.00
	-U1- (User)	loadable via Paraly SW 112
Permissible calibration range	Zero point	6 ... 8 pH
	With ISFET:	-750 ... +750 mV
	Operating point (asymmetry)	
	Slope	Approx. 74 ... 104 %
	(possibly restricting notes from Sensoface)	
<b>Calibration timer *</b>	Interval 1 ... 99 days, can be switched off	
<b>Sensoface</b>	Provides information on the sensor condition	
Evaluation of	zero/slope, response, calibration interval	

\* user-defined

4) ranges depending on Memosens sensor

---

<b>Memosens input</b>	M8 socket, 4 pins, for flexible Memosens lab cable or	
<b>ORP</b>	M12 socket, 8 pins, for flexible connecting cable for Memosens sensors	
Display ranges <sup>4)</sup>	mV	-1999 ... +1999 mV
	Temperature	-50 ... +250 °C / -58 ... +482 °F
Sensor standardization *	ORP calibration (zero adjustment)	
Permissible calibration range	$\Delta$ mV (offset)	-700 ... +700 mV

\* user-defined

4) ranges depending on Memosens sensor

## Cond

<b>Memosens conductivity input</b>	M8 socket, 4 pins, for flexible Memosens lab cable or M12 socket, 8 pins, for flexible connecting cable for Memosens sensors
Measuring range	SE 615/1-MS 10 $\mu\text{S}/\text{cm}$ ... 20 $\text{mS}/\text{cm}$ sensor
Measuring cycle	Approx. 1 s
Temperature compensation	Linear 0 ... 20 %/K, reference temperature adjustable nLF: 0 ... +120 °C / +32 ... +248 °F NaCl (ultrapure water with traces) HCl (ultrapure water with traces) NH <sub>3</sub> (ultrapure water with traces) NaOH (ultrapure water with traces)
Display resolution (autoranging)	Conductivity 0.001 $\mu\text{S}/\text{cm}$ ( $c < 0.05 \text{ cm}^{-1}$ ) 0.01 $\mu\text{S}/\text{cm}$ ( $c = 0.05 \dots 0.2 \text{ cm}^{-1}$ ) 0.1 $\mu\text{S}/\text{cm}$ ( $c > 0.2 \text{ cm}^{-1}$ ) Resistivity 00.00 ... 99.99 $\text{M}\Omega \text{ cm}$ Salinity 0.0 ... 45.0 g/kg (0 ... +30 °C / +32 ... +86 °F) TDS 0 ... 1999 mg/l (+10 ... +40 °C / +50 ... +104 °F) Concentration 0.00 ... 100 % by wt
<b>Concentration determination</b>	NaCl 0 – 26 wt% (0 °C / +32 °F) ... 0 – 28 wt% (+100 °C / +212 °F) HCl 0 – 18 wt% (-20 °C / -4 °F) ... 0 – 18 wt% (+50 °C / +122 °F) NaOH 0 – 13 wt% (0 °C / +32 °F) ... 0 – 24 wt% (+100 °C / +212 °F) H <sub>2</sub> SO <sub>4</sub> 0 – 26 wt% (-17 °C / -1.4 °F) ... 0 – 37 wt% (+110 °C / +230 °F) HNO <sub>3</sub> 0 – 30 wt% (-20 °C / -4 °F) ... 0 – 30 wt% (+50 °C / +122 °F) H <sub>2</sub> SO <sub>4</sub> 94 – 99 wt% (-17 °C / -1.4 °F) ... 89 – 99 wt% (+115 °C / +239 °F) HCl 22 – 39 wt% (-20 °C / -4 °F) ... 22 – 39 wt% (+50 °C / +122 °F) HNO <sub>3</sub> 35 – 96 wt% (-20 °C / -4 °F) ... 35 – 96 wt% (+50 °C / +122 °F) H <sub>2</sub> SO <sub>4</sub> 28 – 88 wt% (-17 °C / -1.4 °F) ... 39 – 88 wt% (+115 °C / +239 °F) NaOH 15 – 50 wt% (0 °C / +32 °F) ... 35 – 50 wt% (+100 °C / +212 °F)
<b>Sensor standardization</b>	Cell constant Input of cell constant with simultaneous display of conductivity value and temperature Input of solution Input of conductivity of the calibration solution with simultaneous display of cell constant and temperature Auto Automatic determination of the cell constant with KCl solution or NaCl solution

<b>Memosens input for amperometric oxygen</b>	M8 socket, 4 pins, for flexible Memosens lab cable or M12 socket, 8 pins, for flexible connecting cable for Memosens sensors	
Display ranges <sup>4)</sup>	Saturation	0.000 ... 200.0 %
	Concentration	000 µg/l ... 20.00 mg/l
	Partial pressure	0.0 ... 1000 mbar
Temperature meas. range <sup>4)</sup>	-20 ... +150 °C / -4 ... 302 °F	
<b>Sensor standardization</b>	Automatic calibration in air (100 % RH) Zero calibration	
<b>Storage</b>	in quiver with moistening sponge	
<b>Input optical oxygen</b>	M12 socket, 8 pins	
OXY meas. ranges at +20 °C / +68 °F	Saturation	0.000 ... 200.0 %
	Concentration	000 µg/l ... 20.00 mg/l
	Partial pressure	0.0 ... 1000 mbar
Response time	t90 < 30 s t99 < 60 s	
Measurement error <sup>1,2,3)</sup>	Zero signal < 0.1 % full saturation	
Temperature meas. range <sup>4)</sup>	0 ... +50 °C / +32 ... +122 °F	
Measurement error <sup>1,2,3)</sup>	Temperature ± 0.2 K	
<b>Sensor standardization</b>	Automatic calibration in air Zero calibration	
<b>Max. overpressure</b>	2.5 bar	
<b>Immersion depth</b>	Min. 60 mm Max. 25 m	
<b>Storage</b>	in quiver with moistening sponge	

1) according to EN 60746-1, at nominal operating conditions

2) ± 1 count

3) plus sensor error

4) ranges depending on Memosens sensor

**A**

AA batteries 6  
 Access control (Option 001 SOP) 47  
 Adapt SOP cal 45  
 ADMIN PIN code lost 47  
 Arrow keys 9  
 Auto, conductivity calibration 27  
 Automatic calibration (conductivity) 27  
 Automatic calibration (pH) 23  
 Automatic self-test 13  
 Auto, pH calibration 23

**B**

Basic interval (limit value logger type) 36  
 Battery capacity 6  
 Battery compartment 6  
 Battery icon 6  
 Battery life, increasing 37  
 Battery replacement 6  
 Buffer set (pH configuration) 14  
 Buttons and functions 9

**C**

Calibration (Cond), auto 27  
 Calibration (Cond), cell constant 28  
 Calibration (Cond), entry of solution 28  
 Calibration, conductivity 27  
 Calibration in air (Oxy) 29  
 Calibration, ISFET 26  
 Calibration, ORP 25  
 Calibration (Oxy), data input 30  
 Calibration, oxygen 29  
 Calibration (Oxy), zero point 30  
 Calibration (pH), Calimatic 23  
 Calibration (pH), data input 24  
 Calibration (pH), manual 24  
 Calibration, pH/ORP combo sensor 24  
 Calibration points (pH) 23  
 Calibration record 10  
 Calibration record, printout 33  
 Calibration, SOP Cal (Option) 44  
 Calibration, Temp.cal (Option) 49  
 Calibration, zero point (operating point), ISFET 26  
 Calimatic calibration 23  
 Cal mode 26  
 Cell constant, conductivity calibration 28  
 Charge level of batteries 6  
 CIP (sensor information) 10  
 Clearing the logger 39  
 Combo sensor, pH/ORP, calibration 24  
 Conductivity calibration 27  
 Conductivity configuration 18  
 Configuration (conductivity) 18  
 Configuration (ORP) 16  
 Configuration (oxygen) 21  
 Configuration (pH) 14  
 Configuring the data logger 37  
 Connecting cable for Memosens 7  
 Connections 7

Connection, USB (battery) 6  
 Control buttons and functions 9  
 Creating a user (Option 001 SOP) 47  
 Curve characteristic (data logger) 40

**D**

Data input, Oxy calibration 30  
 Data input, pH calibration 24  
 Data logger configuration 37  
 Data logger icons 34  
 Data logger, operating modes 35  
 Data logger, starting 39  
 Data logger, stopping 39  
 Data of the meter 50  
 Deleting the logger data 39  
 Delta range (data logger) 35  
 Device info (Information menu) 13  
 Device test (Information menu) 13  
 Difference (logger type) 35  
 Display 9  
 Display icons 8  
 Displaying the logger data 39  
 Display test 13  
 Disposal 5

**E**

Electrolyte, Oxy calibration 29  
 Enabling an option 45  
 Entry of solution, conductivity calibration 28  
 ERROR (error codes) 43  
 Error messages, overview 43  
 Event interval (limit value logger type) 36

**H**

Help texts 42

**I**

Icons 8  
 Icons for data logger 34  
 In air, Oxy calibration 29  
 Increasing the battery life 37  
 Information (menu) 10  
 Info texts 42  
 Inserting the batteries 6  
 Installation factor, calibration 28  
 Interval and difference (logger type) 36  
 Interval (logger type) 35  
 ISFET calibration 26

**K**

Keypad 9  
 Keypad test 13

**L**

Limit value (logger type) 36  
 Lithium-ion battery (start-up) 6  
 Logger type "Difference" 35  
 Logger type "Interval" 35  
 Logger type "Interval and difference" 36  
 Logger type "Limit value" 36

Logger type "Shot" 35  
Loss of PIN code 47

## M

Manual calibration (pH) 24  
meas, switch-on 8  
Measured values, printout 32  
Measuring 31  
Membrane module replacement 29  
MemoLog (Memosens only) 12  
Memosens connecting cable 7  
Memosens sensors, connection 7  
Menu icons 9  
Messages (Information menu) 12  
Micro USB port 7

## N

Network diagram 11  
Note (data logger) 38

## O

on/off, switch-on 8  
Operating modes of data logger 35  
Operating time of cap (optical oxygen) 10  
Operating time of sensor (sensor information) 10  
Option 001 SOP 44  
Option 001 SOP, enabling 45  
Option 002 Temp.cal 44  
Option 002 Temp.cal, enabling 49  
ORP calibration 25  
ORP configuration 16  
Overview of error messages 43  
Overview of status messages 41  
Oxygen calibration (Oxy) 29  
Oxygen configuration 21

## P

Passcode lost 47  
pH calibration 23  
pH configuration 14  
pH/ORP sensor, calibration 24  
PIN code, changing (Option 001 SOP) 47  
PIN code lost 47  
Ports 7  
Power-on 8  
Pressure correction (Oxy configuration) 21  
Pre-trigger (limit value logger type) 36  
Printing 32  
Printing a calibration record 33  
Printing measured values 32

## R

Rechargeable battery, Li-ion 6  
Redox calibration 25  
Registered trademarks 5  
Replacing the electrolyte (Oxy) 29  
Replacing the membrane (Oxy) 29  
Rescue PIN code entry 48

## S

Saving a membrane module replacement (Oxy) 29  
Self-test 13  
Sensoface messages 42  
Sensor connection 7  
Sensor information 10  
Sensor monitor 12  
Sensor network diagram 11  
Sensor verification (Option 001 SOP) 48  
Sensor without temperature detector 31  
Serial number of cap (optical oxygen) 10  
Serial number of sensor (sensor information) 10  
Setup, conductivity 18  
Setup, ORP 16  
Setup, oxygen 21  
Setup, pH 14  
Shot (logger type) 35  
SIP (sensor information) 10  
Slope calibration, oxygen 29  
Softkeys 9  
SOP calibration (Option 001) 46  
Specifications 50  
Starting the data logger 39  
Status messages, overview 41  
Stopping the data logger 39  
Switching on the meter 8  
Symbols in display 8

## T

Table of error messages 43  
Table of info texts 42  
TAG entry (Cond) 27  
TAG entry (ISFET) 26  
TAG entry (ORP) 25  
TAG entry (Oxy) 29  
TAG entry (pH) 23  
Technical data 50  
Temp.cal (Option) 49  
Temperature, manual adjustment 31  
Temperature probe, connection 7  
Toggling the measured value display 31  
Trademarks 5

## U

USB port (battery) 6  
USB port, micro 7  
User management (Option 001 SOP) 47

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The screenshot shows the main menu of the 'Datenlogger'. At the top is a green header with the title 'Datenlogger'. Below it, the main display area shows 'pH 7.000' and '25.0 °C'. A smiley face icon is to the left of the date and time '15.08.2012 10:45:36', which is circled in red. Below this, there are two columns of text: 'Messstelle ABC' and 'Notiz Sensor Analog'. At the bottom left, it says '1/3' with the text 'Mit Pfeiltasten blättern' and an arrow pointing to it. At the bottom are two green buttons labeled 'Zurück' and 'Grafik'. Red arrows point from external text labels to various elements: 'Gemessene Hauptgröße' to 'pH 7.000', 'Gemessene Temperatur' to '25.0 °C', 'Datum und Uhrzeit' to the circled date and time, and 'Mit Pfeiltasten blättern' to '1/3'.

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The screenshot shows the 'Grafik' (graph) view. It features a line graph with two data series: a red line for pH and a blue line for temperature. The y-axis has labels '7.50', '6.50', 'pH', and '°C'. The x-axis shows 'pH 7.21', '25.0 °C', '24.08.2012', and '17:12:50'. A red oval highlights a specific data point on the graph, with an arrow pointing to it from the text 'Eintrag mit Pfeiltasten auswählen'. At the bottom left, it says '2/3' with the text 'Eintrag 2 von 3' and an arrow pointing to it. At the bottom are two green buttons labeled 'Zurück' and 'Grafik'. Red arrows point from external text labels to '2/3', 'Messwerte' (pointing to the graph area), 'Datum' (pointing to the date), and 'Uhrzeit' (pointing to the time).

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