

Knick >

Protos II 4400 (X) **Process Analysis System**

User Manual

Protos MSU 4400(X)-180 **Communication Module** for Memosens Sensors





Read before installation. Keep for future use.



Supplemental Directives

READ AND SAVE THIS DOCUMENT FOR FUTURE REFERENCE. BEFORE ATTEMPTING TO ASSEMBLE, INSTALL, OPERATE OR MAINTAIN THE PRODUCT, PLEASE ENSURE A COMPLETE UNDERSTANDING OF THE INSTRUCTIONS AND RISKS DESCRIBED HEREIN. ALWAYS OBSERVE ALL SAFETY INFORMATION. FAILURE TO COMPLY WITH INSTRUCTIONS IN THIS DOCUMENT COULD RESULT IN SERIOUS INJURY AND/OR PROPERTY DAMAGE. THIS DOCUMENT IS SUBJECT TO CHANGE WITHOUT NOTICE.

These supplemental directives explain how safety information is laid out in this document and what content it covers.

Safety Chapter

This document's safety chapter is designed to give the reader a basic understanding of safety. It illustrates general hazards and gives strategies on how to avoid them.

Safety Guide

The separate safety guide is designed to give the reader a basic understanding of safety. It illustrates general hazards and gives strategies on how to avoid them.

Warnings

This document uses the following warnings to indicate hazardous situations:

Symbol	Category	Meaning	Remark	
A	WARNING	Designates a situation that can lead to death or serious (irreversible) injury.	The warnings contain information on how to	
A	CAUTION	Designates a situation that can lead to slight or moderate (reversible) injury.	avoid the hazard.	
None	NOTICE	Designates a situation that can lead to property or environmental damage.		

Related Documents

- Protos II 4400(X) Safety Guide
- Protos II 4400(X) Basic Unit User Manual

Table of Contents

MSU 4400(X)-180 Module

Intended Use	6
Package Contents	7
Safety	8
Operation in Hazardous Locations – MSU4400X-180 Module	8
Firmware Version	9
Terminal Plate	10
Terminal Assignments	11
Unical 9000(X) Connection	
Inserting the Module	13
Parameter Setting	14
Module Configuration: Process Variable	15
pH Parameter Setting	17
pH Calibration / Adjustment	
Function Check (HOLD) During Calibration / Adjustment	
Calibration Mode: Calimatic	34
Calibration Mode: Manual	
Calibration Mode: Product	
Calibration Mode: Data Entry	
Calibration Mode: ISFET Zero	
Calibration Mode: Temperature	42
pH Maintenance Functions	
Sensor Monitor	
pH Diagnostic Functions	44
ORP Parameter Setting	46
ORP Calibration / Adjustment	49
Function Check (HOLD) During Calibration / Adjustment	
Calibration Mode: ORP Data Entry	
Calibration Mode: ORP Adjustment	
Calibration Mode: ORP Check	
Calibration Mode: Temperature	
ORP Maintenance Functions	
Sensor Monitor	58
ORP Diagnostic Functions	59
pH, ORP Messages	61

Table of Contents

MSU 4400(X)-180 Module

Oxy Parameter Setting	65
Oxy Calibration / Adjustment	70
Recommendations for Calibration	
Function Check (HOLD) During Calibration / Adjustment	74
Calibration Mode: In Air	
Calibration Mode: In Water	
Calibration Mode: Data Entry	
Calibration Mode: Product	
Calibration Mode: Zero Point	
Calibration Mode: Temperature	
Oxy Maintenance Functions	
Sensor Monitor	
Membrane Body Replacement	85
Oxy Diagnostic Functions	86
Oxy Messages	88
Cond Parameter Setting	93
USP Function (Cond)	
Concentration (Cond)	
pH Value Calculation (Cond)	102
Cond Calibration / Adjustment	104
Function Check (HOLD) During Calibration / Adjustment	
Calibration Mode: Automatic	
Calibration Mode: Manual	
Calibration Mode: Product	
Calibration Mode: Installation Factor	
Calibration Mode: Data Entry	
Calibration Mode: Temperature	118
Cond Maintenance Functions	119
Sensor Monitor	119
Cond Diagnostic Functions	120
Cond Messages	122

Table of Contents

MSU 4400(X)-180 Module

Condl Parameter Setting	12/
USP Function (Condl)	
Concentration (Condl)	135
Condl Calibration / Adjustment	138
Function Check (HOLD) During Calibration / Adjustment	142
Calibration Mode: Automatic	
Calibration Mode: Manual	146
Calibration Mode: Product	
Calibration Mode: Zero Point	
Calibration Mode: Installation Factor	
Calibration Mode: Data Entry	
Calibration Mode: Temperature	153
Condl Maintenance Functions	154
Sensor Monitor	154
Condl Diagnostic Functions	155
Condl Messages	
Decommissioning	159
Specifications	160
•	
Buffer Tables	162
•	162
Buffer Tables Mettler-Toledo Buffer Table	
Mettler-Toledo Buffer Table	
Mettler-Toledo Buffer Table	162 162 163 164 165
Mettler-Toledo Buffer Table	
Mettler-Toledo Buffer Table	
Buffer Tables Mettler-Toledo Buffer Table Buffer Table, Knick CaliMat Buffer Table, DIN 19267 NIST Standard Buffer Table (DIN 19266: 2000-01) Buffer Table, NIST Technical Buffers Buffer Table, Hamilton Buffer Table, Kraft Buffer Table, Hamilton A	
Buffer Tables Mettler-Toledo Buffer Table Buffer Table, Knick CaliMat Buffer Table, DIN 19267 NIST Standard Buffer Table (DIN 19266: 2000-01) Buffer Table, NIST Technical Buffers Buffer Table, Hamilton Buffer Table, Kraft Buffer Table, Hamilton A Buffer Table, Hamilton B	
Mettler-Toledo Buffer Table	
Mettler-Toledo Buffer Table	
Mettler-Toledo Buffer Table Buffer Table, Knick CaliMat Buffer Table, DIN 19267 NIST Standard Buffer Table (DIN 19266: 2000-01) Buffer Table, NIST Technical Buffers Buffer Table, Hamilton Buffer Table, Kraft Buffer Table, Hamilton A Buffer Table, Hamilton B Buffer Table, Hamilton B Buffer Table, Hamilton B Buffer Table, HACH Buffer Table, Ciba Buffer Table, Reagecon	
Mettler-Toledo Buffer Table	

Intended Use

The MSU4400(X)-180 multiparameter module is a multifunctional communication unit that provides up to three RS-485 interfaces for simultaneous use. It enables the connection and operation of up to three Memosens sensors for simultaneous measurement of pH, oxidation-reduction potential (ORP), oxygen (add-on function FW4400-015), and conductivity (contacting, inductive) as well as the Unical 9000 electro-pneumatic controller.

An analog current input that supplies a pressure transmitter signal is used for oxygen sensor pressure correction during measurement and calibration.

The second and third sensor inputs (channels B and C) can be enabled via TAN.

- Channel B: add-on function FW4400-014
- Channels B+C: add-on function FW4400-018

Channel C is used for Unical actuation.

This user manual describes the use of the module in conjunction with Memosens sensors.

For a description of Unical actuation, see the Unical 9000/Protos II 4400 user manual.

The MSU4400X-180 module is intended for operation in hazardous locations that require equipment of Group II, device category 2(1), gas/dust. The MSU4400-180 module may not be used in hazardous locations.

Package Contents

All devices:

- · Measuring module
- · Installation guide
- Test report 2.2 according to EN 10204
- · Adhesive label with terminal assignments

Additionally for Ex version MSU4400X-180:

- Attachment to certificates (KEMA 03ATEX2530, IECEx DEK 11.0054)
- EU Declaration of Conformity
- · Control drawings

Safety

Operation in Hazardous Locations - MSU4400X-180 Module

The module is approved for operation in hazardous locations.

When installing the product in a hazardous location, observe the information in the supplements to the certificates and, if applicable, the relevant control drawings.

Observe all applicable local and national codes and standards for the installation of electrical equipment in hazardous locations. For orientation, please refer to IEC 60079-14, EU directives 2014/34/EU and 1999/92/EC (ATEX), NFPA 70 (NEC), ANSI/ISA-RP12.06.01.

A WARNING! Risk of impairment of explosion protection.

- Modules that have already been used must be subjected to a professional routine test before they may be operated in another type of protection.
- Before the product is commissioned, the operator must provide proof that
 the product is approved for connection to other equipment (including cables
 and wires).
- Connecting components designed for explosive atmospheres and those not designed for explosive atmospheres (mixed equipping) is not permitted.
- In hazardous locations, the device may only be cleaned with a damp cloth to prevent electrostatic charging.

Maintenance

Protos modules cannot be repaired by the user. For inquiries regarding module repair, please contact Knick Elektronische Messgeräte GmbH & Co. KG at www.knick.de.

Personnel Requirements

The operating company shall ensure that any personnel using or otherwise interacting with the product is adequately trained and has been properly instructed.

The operating company shall comply and cause its personnel to comply with all applicable laws, regulations, codes, ordinances and relevant industry qualification standards related to product. Failure to comply with the foregoing shall constitute a violation of operating company's obligations concerning the product, including but not limited to an unintended use as described in this document.

Firmware Version

MSU4400(X)-180 Module Firmware: Firmware version 01.xx.xx


Module Compatibility	MSU4400-180	MSU4400X-180
Protos II 4400 with FRONT firmware version 01.03.xx or higher	x	
Protos II 4400X with FRONT firmware version 01.03.xx or higher		x

Information on the firmware version history can be found at www.knick.de.

Query Current Device Firmware/Module Firmware


When the device is in measuring mode:

Press the **menu** key and open the Diagnostics menu: Device Description



Terminal Plate

MSU4400-180 Module Terminal Plate:



1	Name of manufacturer	8	Product number/serial number/ production year and week
2	Product name	9	Installation instructions
3	Designation of origin and address of the manufacturer	10	Terminal assignments
4	UKCA mark	11	Special conditions and danger points
5	CE mark	12	Serial number
6	WEEE mark	13	Model description
7	Permissible ambient temperature		

Note: The "Uniclean" function is currently unavailable.

Terminal Plate Adhesive Labels

The terminal plate adhesive labels for the lower modules can be attached to the inside door. This simplifies maintenance and service.



Terminal Assignments


	Memosens C	Cable Wiring		Wiring	
Terminal	Wire Color			Unical	Terminal
1	Brown (BN)	Power supply +	- ·: ›		
2	Transparent	Shield	en:		
3	Green (GN)	RS485 (A)	Channel A: Memosens		
4	Yellow (YE)	RS485 (B)	_ha ∕ler		
5	White (WH)	Power supply – (GND)	- 0 2		
6	Brown (BN)	Power supply +	∨		
7	Transparent	Shield	Channel B: Memosens		
8	Green (GN)	RS485 (A)	ַ עַעַ טָעָע זיין		
9	Yellow (YE)	RS485 (B)	⁻ha ∕ler		
10	White (WH)	Power supply – (GND)	0 2		
11		N.C.		7 V (Power)	20
12		N.C.	;	N.C.	
13	Brown (BN)	Power supply +	Channel C: Memosens/Unica	N.C.	
14	Transparent	Shield	Channel o	N.C.	
15	Green (GN)	RS485 (A)	ha	RS 485 B	18
16	Yellow (YE)	RS485 (B)	en C	RS 485 A	17
17	White (WH)	Power supply – (GND)	Σ	GND	19
18		+ Current input	·		
19		– 0(4) 20 mA			

The second and third sensor inputs (channels B and C) can be enabled via TAN.

- Channel B: Add-on function FW4400-014
- Channels B+C: Add-on function FW4400-018

Channel C (terminals 11 and 15 ... 17) is designed for Unical actuation. For wiring, see also the next page.

Unical 9000(X) Connection



For a detailed description of the Unical 9000(X) electro-pneumatic control, see the Unical 9000/Protos II 4400 user manual.

Inserting the Module

A CAUTION! Electrostatic discharge (ESD).

The modules' signal inputs are sensitive to electrostatic discharge. Take measures to protect against ESD before inserting the module and wiring the inputs.

NOTICE! Strip the insulation from the wires using a suitable tool to prevent damage.



- 01. Switch off the power supply to the device.
- 02. Open the device (loosen the 4 screws on the front).
- 03. Plug the module into the slot (D-SUB connector), as in the image.
- 04. Tighten the module's fastening screws.
- 05. Connect the sensor cable.
- 06. Check whether all connections are correctly wired.
- 07. Close the device and tighten the screws on the front.
- 08. Switch on the power supply.

A CAUTION! Risk of losing the specified ingress protection.

Fasten the cable glands and screw together the housing correctly. Observe the permissible cable diameters and tightening torques (see basic unit specifications).

Insert blanking plugs or sealing inserts if necessary.


A CAUTION! Incorrect parameter settings or adjustments can result in incorrect outputs.

Protos must therefore be commissioned by a system specialist, all its parameters must be set, and it must be fully adjusted.

NOTICE! The "function check" (HOLD) NAMUR contact is active during parameter setting. The behavior of the current outputs depends on the parameter setting, i.e., they may be frozen at the last measurement or set to a fixed value. The red "Alarm" LED blinks.

Measurement operations must not be carried out while Protos is in the function check (HOLD) mode, as this may put the user at risk due to unexpected system behavior.

For a description of NAMUR contacts, see the basic unit user manual.




For a description of Unical actuation, see the Unical 9000/Protos II 4400 user manual.

Module Configuration: Process Variable

Note: Function check (HOLD) active

The process variable, operating mode, and functionality are separately selected for each channel. If the Auto process variable is selected, the functionality settings are not shown.

Default: Process variable off




Module Configuration: Process Variable

Menu selection: Parameter Setting ▶ MSU ... Module ▶ Channel ...

Operating Mode: Memosens					
Process	Auto	Automatic selection of process variables in the measurement display			
Variable	рН	Functionality: pH, ISFET, pH/ORP, ORP			
	Conductivity	Functionality: 2-electrode / 4-electrode sensor			
	Conductivity (ind.)	Functionality: Condl			
	Oxygen	Functionality: Amperometric			
Operating Mode: SE 670, SE680K					
Process Variable	Conductivity (Ind.)	For use with SE 670, SE680K			

Note: Function check (HOLD) active




Menu selection: Parameter Setting ▶ MSU ... Module ▶ Memosens pH(/ORP)

Parameter	Default	Description, Options
		· · · · ·
Input Filter	T	
Pulse Suppression	Off	Enable/disable suppression of interference pulses.
Sensor Data		
Sensoface	On	Enable/disable display of Sensoface messages and icons.
Sensor Monitoring Details (see page 20)		pH sensors: Slope, Zero Point, Sensocheck, Response Time, Sensor Operating Time, Sensor Wear, SIP Counter
		ISFET sensors: Slope, ISFET Operating Point, ISFET Leakage Current, Response Time, Sensor Operating Time, Sensor Wear, SIP Counter
		pH/ORP sensors: Slope, Zero Point, ORP Offset, Sensocheck Reference Electrode, Sensocheck Glass Electrode, Response Time, Sensor Operating Time, Sensor Wear, CIP Counter, SIP Counter, Autoclaving Counter
Cal Presettings (see page 22)		
Calibration Mode	Calimatic	Presetting for calibration mode
		pH sensors: Calimatic, Manual, Product, Data Entry, Temperature
		ISFET: Calimatic, Manual, Product, ISFET Zero Point, Data Entry, Temperature
		pH/ORP sensors: Calimatic, Manual, Product, Data entry, ORP Data Entry, ORP Adjustment, ORP Check, Temperature
Buffer Set	Knick CaliMat	For Calimatic: Select buffer set
Cal Points	Auto	For Calimatic and manual: Auto, 1-point, 2-point, 3-point
Drift Check	Standard	Fine: 1.2 mV/min (stop after 180 s) Standard: 2.4 mV/min (stop after 120 s) Coarse: 3.75 mV/min (stop after 90 s)
Cal Timer Monitoring	Off	Off, Auto: 0168h, Individual
Adaptive Cal Timer	Off	Off, On
TC Process Medium (see page 24	1)	
Temperature Comp.	Off	Off, Linear, Ultrapure Water, Table
•		Linear: Enter Temp Factor +nn.nn%/K.
		Table: TC values specifiable, step size 5 °C / 9 °F

Menu selection: Parameter Setting ▶ MSU ... Module ▶ Memosens pH(/ORP)

Parameter	Default	Description, Options			
ORP / rH Value (for pH/ORP sens	ORP / rH Value (for pH/ORP sensors)				
Reference Electrode	Ag/AgCl, KCl 3mol	Ag/AgCl, KCl 1mol, Ag/AgCl, KCl 3mol, Hg, Tl/TlCl, KCl 3.5mol, Hg/Hg ₂ SO ₄ , K ₂ SO ₄ sat.			
ORP Conversion to SHE	No	Enable or disable ORP conversion to standard SHE hydrogen electrode.			
Calculate rH with Factor	No	Calculate rH with or without factor.			
Delta Function (see page 25)					
Delta Function	Off	pH sensors: Off, pH pH/ORP sensors: Off, pH, mV ORP, rH ISFET sensors: Off, pH			
Messages (see page 26) pH sensors, ISFET sensors:					
pH Value Messages	Off	Off, Max. Device Limits, Variable Limits			
Temperature Messages	Off	Off, Max. Device Limits, Variable Limits			
pH Voltage Messages	Off	Off, Max. Device Limits, Variable Limits			
pH/ORP sensors:					
pH Value Messages	Off	Off, Max. Device Limits, Variable Limits			
ORP Messages	Off	Off, Max. Device Limits, Variable Limits			
rH Value Messages	Off	Off, Max. Device Limits, Variable Limits			
Temperature Messages	Off	Off, Max. Device Limits, Variable Limits			
pH Voltage Messages	Off	Off, Max. Device Limits, Variable Limits			



CIP Counter¹⁾/SIP Counter

CIP/SIP cycles are used to clean or sterilize process-wetted parts in the process. Depending on the application, either one chemical (alkaline solution, water) or several chemicals (alkaline solution, water, acidic solution, water) are used.

- CIP temperature > 55 °C/131 °F
- SIP temperature > 115 °C/239 °F

The cleaning (cleaning in place) and sterilization (sterilization in place) cycles are counted to measure the load on the sensor, e.g., in biotechnology applications.

Note: If measurements are generally taken at high temperatures $(> 55 \, ^{\circ}\text{C}/131 \, ^{\circ}\text{F})$, the counters should be switched off.

When a CIP/SIP counter is switched on, a maximum number of cycles can be entered. A message can be set to indicate when a counter has reached a specified value.

Note: A CIP or SIP cycle is only entered into the logbook 2 hours after the start to ensure that the cycle is complete.


Autoclaving Counter¹⁾

Autoclaving cycles are counted to help measure the load on the sensor.

When the autoclaving counter is switched on, a maximum number of cycles can be entered. A message can be set to indicate when a counter has reached a specified value.

Menu

Display Action



The calibration presettings can be defined in the parameter settings or adjusted directly in the Calibration menu prior to the calibration.

Presettings for Calibration

Parameter Settings Memosens pH... Cal Presettings

Calibration Mode: Presetting of calibration mode, e.g., Calimatic,

Manual, Product, Data Entry, Temperature

If Calimatic automatic calibration is selected, the buffer set to be used must also be selected.

Cal Points: Selection of how many calibration points are to be used in the calibration (for Calimatic and manual calibration).

Selection: 1-, 2-, or 3-point

Drift Check: Setting of drift check sensitivity

Fine: 1.2 mV/min (stop after 180 s) Standard: 2.4 mV/min (stop after 120 s) Coarse: 3.75 mV/min (stop after 90 s)

Default: Standard

Calibration Timer

When a preset calibration interval expires, the calibration timer generates a message text to indicate the need for calibration.

"Auto": Interval is set to 168 h

"Individual": Any interval can be set

"Off": No monitoring of the calibration timer

Note: If Sensoface is enabled, (a) (neutral) is displayed once 80 % of the time interval has expired. On expiry of the complete time interval, (sad) is displayed and a corresponding message is generated (see p. 20). The NAMUR icon (a) is displayed. If the current outputs are configured accordingly, a 22-mA error signal is generated (see basic unit user manual).



Adaptive Calibration Timer: The time until the next calibration is automatically shortened depending on the temperature and pH value.

Old sensor = timer expires faster.

The following measuring conditions shorten the adaptive calibration timer interval:

- Temperatures above 30 °C / 86 °F
- pH ranges below pH 2 or above pH 12

The message test is displayed in the Diagnostics menu:

Diagnostics ▶ Message List

The calibration timer is reset after each calibration.

Menu Display Action TC Process Medium 7 You can select the following options: IBTC Process Medium (Administrator) Off Temperature Comp. ▼Table Linear (enter TC coefficient) TC at 0 °C (32 °F) 10.00% TC at 5 °C (41 °F) 10.00% Ultrapure Water TC at 10 °C (50 °F) 10.00% Table TC at 15 °C (59 °F) 10.00% Back

Linear Temperature Compensation of Process Medium

If the medium's pH value changes in linear fashion with the temperature, the temperature coefficient TC can be determined for temperature compensation in %/K as follows:

$$TC = (pH_{25} - pH_{T}) \cdot 100 / (25 \degree C - T) [\%/K]$$

TC Temperature coefficient [%/K]

pH₂₅ pH value at 25 °C

pH_T pH value at measuring temperature T

T Measuring temperature [°C]

Table


When using process media with a known pH value temperature response, the pH output value can be corrected using a table. The percentage deviation from the measured value in % can be entered for temperatures between 0 and 95 °C in steps of 5 °C. The pH output value is then corrected by the corresponding percentage deviation from the measured value in %, depending on the measuring temperature. Table values are linearly interpolated. If the temperature falls below or exceeds the specified value (< 0 °C or > 95 °C), the last value in the table is used for calculation.


The table must be completed with the following values in steps of 5 °C:

((pH25 / pHT) − 1) • 100 [%] pH25 pH value at 25 °C

pHT pH value at measuring temperature T

Note: If TC correction for the process medium is enabled, "TC" is shown on the display in measuring mode.





Back

Note: During calibration, the function check (HOLD) operating state is active for the corresponding module channel. The assigned current outputs and relay contacts behave as configured (BASE module).


The calibration data is saved in the Memosens sensor. This means that Memosens sensors can be cleaned, reconditioned, calibrated, and adjusted away from the measurement location, e.g., in a laboratory. Sensors in the system are replaced on-site by adjusted sensors.

Calibration: Detecting deviations without readjustment **Adjustment:** Detecting deviations with readjustment

Adjustment

Adjustment means that the values determined by a calibration are applied to the sensor. The zero point and slope values determined during calibration are entered in the adjustment record:

Diagnostics MSU ... Module Memosens pH Cal/Adj Record
These values are not used to calculate the process variables until the calibration has been terminated with an adjustment.



pH Calibration / Adjustment Explanations

Each pH sensor has a separate zero point and a separate slope. Both values change as a result of aging and wear. The voltage supplied by the pH sensor is corrected by Protos for the zero point and the electrode slope of the pH sensor, and displayed as the pH value.

During calibration, the sensor's deviation is initially determined (zero point, slope).

The sensor is immersed in buffer solutions with a precisely known pH value for this purpose. Protos measures the voltages of the sensors and the temperature of the buffer solution, using this information to calculate the sensor's zero point and slope.

Note: Make sure that the temperature of the sensor and the temperature of the buffer solution are not too far apart. The ideal temperature is 25 °C/77 °F.

Calibration Values Determined During Calibration

Zero Point The pH value at which the pH sensor supplies the voltage 0 mV.

The zero point is different for each sensor and changes with age

and wear.

Slope The slope of a sensor is the voltage change per pH unit.

With an ideal sensor, it is -59.2 mV/pH.

Temperature The temperature of the process solution must be logged, since

the pH measurement is temperature-dependent. Many sensors

feature an integrated temperature detector.

There are limits that are calculated during calibration when monitoring glass and reference impedances. The following limits apply to standard glass electrodes:

• Temperature range: 0 ... 80 °C/32 ... 176 °F

• Impedance range: 50 ... 250 M Ω (at 25 °C/77 °F)

Note: With Memosens sensors, the calibration data is stored in the sensor. This allows the use of precalibrated sensors.

When Protos is used for precalibrating sensors in the lab, you can use the calibration routines described below.

Calibration Methods


One-Point Calibration

The sensor is calibrated with just one buffer solution. A one-point calibration is useful and permissible if the measured values are close to the sensor zero point, such that the change in sensor slope is of negligible significance. The sensor's zero point is subsequently adjusted, while the slope remains unchanged.

Two-Point Calibration

The sensor is calibrated with two buffer solutions. This makes it possible to calculate the sensor's zero point and slope. The sensor's zero point and slope are subsequently adjusted. A two-point calibration is required in the following cases, for example:

- · The sensor has been replaced.
- The pH measured value covers a large range.
- · The pH measured value is far from the sensor zero point.
- The pH value needs to be determined with high precision.
- · The sensor is subject to heavy wear.




- 1 First point of first buffer solution
- 3 Recommended measuring range
- **2** Second point of second buffer solution
- 4 Result of ideal calibration at 25 °C/77 °F

Three-Point Calibration

The sensor is calibrated with three buffer solutions.

The zero point and slope are calculated using a line of best fit in accordance with DIN 19268. The sensor's zero point and slope are subsequently adjusted.



1	Value of first buffer solution	3	Value of third buffer solution
2	Value of second buffer solution	4	Rise

Temperature Compensation during Calibration

The pH sensor's slope is temperature-dependent. The measured voltage must therefore be corrected by the temperature influence.


The buffer solution's pH value is temperature-dependent. During calibration, the temperature of the buffer solution must therefore be known so that the actual pH value can be obtained from the buffer table.


Automatic Temperature Compensation

Protos measures the temperature of the buffer solution using the temperature detector integrated in the Memosens sensor.

Function Check (HOLD) During Calibration / Adjustment

Behavior of the signal and switching outputs during calibration / adjustment






Calibration Mode: Calimatic

During automatic calibration with Knick Calimatic, the sensor is immersed in one, two, or three buffer solutions. On the basis of the sensor voltage and the measured temperature, Protos automatically detects the nominal buffer value. The buffer solutions can be used in any order, but they must be part of the buffer set defined during parameter setting.

Calimatic accounts for the temperature dependence of the buffer value. All calibration data is converted to a reference temperature of 25 °C/77 °F.



Calibration Procedure

NOTICE! Use only new, undiluted buffer solutions that are part of the configured buffer set.

When you access calibration, the calibration values selected in the calibration presettings are used. These settings can still be changed in the Calibration menu. If you do not want to perform calibration, use the left softkey to go back one level and abort the calibration.

- 01. Select "Calimatic" calibration mode and press **enter** to confirm.

 Number of calibration points and buffer set are calibrated as in Cal Presettings; See Calibration Presettings, p. 22
- 02. If required, change the number of calibration points and the buffer set.
- 03. Take the sensor out of the medium and rinse it in deionized water.

▲ CAUTION! Risk of electrostatic discharge.


Do not wipe the sensor or dab it dry.

- 04. Immerse the sensor in the first buffer solution.
- 05. Start calibration with the *right softkey: Next*.
 - ✓ Calibration with first buffer is carried out.

The following parameters are displayed: Sensor Voltage, Calibration Temperature, Nominal Buffer Value, and Response Time.

The time for the measuring voltage to stabilize can be shortened with the *left softkey: Exit* (without drift check: reduced accuracy of calibration values). The response time indicates how long the sensor needs until the measuring voltage is stable. If the sensor voltage or the measured temperature fluctuate considerably, the calibration procedure is aborted after around 2 minutes. In this case, calibration needs to be restarted.

- 06. For one-point calibration: Exit calibration by pressing the softkey.
- 07. For two-point calibration: Rinse the sensor well with deionized water.
- 08. Immerse the sensor in the second buffer solution.
- 09. Start calibration with the right softkey: Next.
 - √ Calibration with second buffer is carried out.
- 10. Proceed as for one-point calibration.
- 11. For three-point calibration, the process uses the third buffer accordingly.
- ✓ By pressing the *right softkey: Adjust*, the calculated calibration values are used to calculate the measured values in Protos and also saved in the Memosens sensor.



Calibration Mode: Manual

During calibration with manual entry of buffer values, the sensor is immersed in one, two, or three buffer solutions.

Protos displays the measured temperature. The temperature-corrected buffer values must then be manually entered. For this purpose, take the buffer value that goes with the displayed temperature from the buffer table (e.g., on the bottle). Intermediate values must be interpolated.

All calibration data is converted to a reference temperature of 25 °C/77 °F.

Calibration Procedure

NOTICE! Use only new, undiluted buffer solutions that are part of the configured buffer set.

When you access calibration, the calibration values selected in the calibration presettings are used. These settings can still be changed in the Calibration menu. If you do not want to perform calibration, use the left softkey to go back one level and abort the calibration.

Select "Manual" calibration mode and confirm by pressing enter.

Number of calibration points are calibrated as in Cal Presettings;

See Calibration Presettings, p. 22

- 01. If required, change the number of calibration points.
- 02. Enter the first buffer value.
- 03. Continue with the *right softkey: Next*.
- 04. Take the sensor out of the medium and rinse it well in deionized water.

▲ CAUTION! Risk of electrostatic discharge.

Do not wipe the sensor or dab it dry.

- 05. Immerse the sensor in the first buffer solution.
- 06. Start calibration with the *right softkey: Next*.
 - ✓ Calibration with first buffer is carried out.

The following parameters are displayed: Sensor Voltage, Calibration Temperature, Nominal Buffer Value, and Response Time.

The time for the measuring voltage to stabilize can be shortened with the *left softkey: Exit* (without drift check: reduced accuracy of calibration values). The response time indicates how long the sensor needs until the measuring voltage is stable. If the sensor voltage or the measured temperature fluctuate considerably, the calibration procedure is aborted after around 2 minutes. In this case, calibration needs to be restarted.

- 07. For one-point calibration: Exit calibration by pressing the softkey.
- 08. For two-point calibration: Rinse the sensor well with deionized water.
- 09. Immerse the sensor in the second buffer solution.
- 10. Start calibration with the *right softkey: Next*.

 ✓ Calibration with second buffer is carried out.
- 11. Proceed as for one-point calibration.
- 12. For three-point calibration, the process uses the third buffer accordingly.
- ✓ By pressing the *right softkey: Adjust*, the calculated calibration values are used to calculate the measured values in Protos and also saved in the Memosens sensor.

Calibration Mode: Product

If the sensor cannot be removed – e.g., for sterility reasons – its zero point can be calibrated by "sampling". The currently measured value of the process is saved in the device for this purpose. A sample is taken directly afterward at the measuring point. The sample's pH value is measured in the laboratory.

The reference value is entered in the device.

Protos calculates the sensor's zero point from the difference between the measured and reference values. The slope is not changed in the process.

NOTICE! The sample's pH value is temperature-dependent. The reference measurement should be carried out at the sample temperature shown on the display. The sample should be transported in a vacuum flask.

The sample's pH value may also be falsified if volatile substances escape.

Calibration Procedure


When you access calibration, the calibration values selected in the calibration presettings are used. These settings can still be changed in the Calibration menu. If you do not want to perform calibration, use the left softkey to go back one level and abort the calibration.

- 01. Select "Product" calibration mode and press enter to confirm.
- 02. Prepare for sampling.
- 03. Start with the right softkey: "Next".

Product calibration is performed in 2 steps.

Step 1:

- 04. Take sample.
 - ✓ The measured value and temperature at the time of sampling are displayed.
- 05. Save with the *right softkey: Save*.



Step 2: Lab value has been measured.

- 06. Open the Product Calibration menu again.
- 07. Enter reference value ("Lab value").
- 08. Confirm with the *right softkey: Next* or repeat calibration with the *left soft-key: Cancel*.
- ✓ By pressing the *right softkey: Adjust*, the calculated calibration values are used to calculate the measured values in Protos and also saved in the Memosens sensor.

Exception: Sample value can be determined and entered on site: **Left softkey: Entry.**

Calibration Mode: Data Entry

Entry of values for the zero point and the slope of the sensor. The values must be known, e.g., determined beforehand in the laboratory.

Calibration Procedure

When you access calibration, the calibration values selected in the calibration presettings are used. These settings can still be changed in the Calibration menu. If you do not want to perform calibration, use the left softkey to go back one level and abort the calibration.

- 01. Select "Data Entry" calibration mode and press enter to confirm.
- 02. Remove the sensor and install the premeasured sensor.
- 03. Continue with the right softkey: Next.
- 04. Enter the zero point and slope measured values.
- ✓ By pressing the *right softkey: Adjust*, the calculated calibration values are used to calculate the measured values in Protos and also saved in the Memosens sensor.

Calibration Mode: ISFET Zero

When using Memosens ISFET sensors for pH measurement, the individual operating point of the sensor first needs to be determined, and should be in the pH 6.5 ... pH 7.5 range. The sensor is immersed in a buffer solution with a pH value of 7.00 for this purpose.

Calibration Procedure

When you access calibration, the calibration values selected in the calibration presettings are used. These settings can still be changed in the Calibration menu. If you do not want to perform calibration, use the left softkey to go back one level and abort the calibration.

- 01. Select the "ISFET Zero Point" calibration mode for setting the operating point for the first sensor calibration.
- 02. Continue with the right softkey: Next.
- 03. Adjust the buffer value if necessary: Default pH 7.00
- 04. Take the sensor out of the medium and rinse it well in deionized water.

A CAUTION! Risk of electrostatic discharge.

Do not wipe the sensor or dab it dry.

- 05. Immerse the sensor in buffer solution.
- 06. Start calibration with the *right softkey: Next*.
 - √ The ISFET operating point is calculated.
- 07. Finally, confirm the ISFET operating point with the *right softkey: Adjust*.

A pH calibration, e.g., Calimatic 2-point calibration, can be performed afterward.

Note: The operating point only needs to be determined once for each ISFET sensor.

Calibration Mode: Temperature

This function is used to adjust the individual temperature detector tolerances or cable lengths for the purpose of increasing the accuracy of the temperature measurement. The adjustment requires an accurate measurement of the process temperature using a calibrated reference thermometer. The measurement error of the reference thermometer should be less than 0.1 K. Adjustment without an accurate measurement may result in falsification of the displayed measured value.

Calibration Procedure

- 01. Select "Temperature" calibration mode and press enter to confirm.
- 02. Enter the measured process temperature and press *enter* to confirm.
 √ The temperature offset is displayed.
- 03. Adjust the temperature detector with the *right softkey: Save*.

The current adjustment and temperature offset data can be displayed in the Diagnostics menu:

▶ MSU ... Module ▶ Memosens ... ▶ Temp. Offset Log

Note: The pH/ORP combo sensor can be calibrated as a pH sensor and/or as an ORP sensor.


See also "ORP Calibration / Adjustment"; p. 49

pH Maintenance Functions


Note: Function check (HOLD) active

The current outputs and relay contacts behave in accordance with the parameter settings. Since the device is in function check (HOLD) mode, certain media can be used to validate the sensor and check the measured values without affecting the signal outputs.

To end the function check, return to measuring mode.



pH Diagnostic Functions



For a detailed description of general diagnostic functions, see the basic unit user manual.

Overview of pH Diagnostic Functions

In diagnostics mode, you can access the following submenus without interrupting the measurement:

Diagnostics MSU ... Module:

Module Diagnostics Protos periodically performs a self-test in the background.

The results can be displayed here.

Diagnostics MSU Module Memosens pH:			
Submenus	Description		
Sensor Information	The Sensor Information submenu shows data from the currently connected Memosens sensor, e.g., manufacturer, order no., serial no., firmware and hardware version, last calibration, operating time.		
Sensor Monitor	The raw measured values, such as pH voltage, glass impedance, reference impedance, etc., are displayed for diagnostic purposes.		


pH Diagnostic Functions


Submenus

pH Sensor Diagram

Description

The sensor diagram clearly indicates the status of the parameters in the connected sensor, including the calibration timer. Inactive parameters are shown in gray and set to 100 % (e.g., disabled calibration timer). The parameter values should lie between the outer (100 %) and inner (50 %) polygon. A warning signal flashes if a value drops below the inner polygon (< 50 %).





Critical range – "inner circle": Value outside tolerance

The tolerance limits (radius of the "inner circle") can be individually adjusted. See Parameter Setting > Sensor Data > Sensor Monitoring Details.


The calibration/adjustment record shows the data from the last calibration/adjustment performed on the currently connected sensor.

The temp. offset log shows the data from the last temperature equalization performed on the currently connected sensor.

pH Cal/Adj Record


Temp. Offset Log

Sensor Wear Monitor



The sensor wear monitor shows the sensor operating time and the maximum temperature during the operating time, as well as wear and the forecast remaining time.


Note: Function check (HOLD) active



Menu Selection: Parameter Setting MSU ... Module Memosens ORP

Parameter	Default	Description, Options
Input Filter		
Pulse Suppression	Off	Enable/disable suppression of interference pulses.
Sensor Data	Oli	Litable/disable supplession of interference pulses.
20.00.24.0		I
Sensoface	On	Enable/disable display of Sensoface messages and icons.
Sensor Monitoring Details (See next page)		Option to enter individual limits for monitoring ORP offset. Option to enter individual limits before a message for sensor operating time and SIP counter is triggered.
Cal Presettings		
Calibration Mode	ORP Data Entry	Presetting for calibration mode: ORP Data Entry, ORP Adjustment, ORP Check, Temperature
ORP Check	Test Period 10 s Test Difference 10 mV	Settings for test period in seconds and test difference in millivolts.
ORP / rH Value		
Ref El	Ag/AgCl, KCl 3 mol	Ag/AgCl, KCl 1 mol Ag/AgCl, KCl 3 mol Hg, Tl/TlCl, KCl 3.5 mol Hg/HgSO ₄ , K2SO ₄ sat.
ORP Conversion to SHE	No	Enable or disable ORP conversion to standard SHE hydrogen electrode.
Calculate rH with Factor	No	When using a pH sensor connected via another module at the same time: Calculate rH with or without factor
Delta Function		
Delta Function	Off	Display deviations from a preset value (delta value); see page 25
Messages		
ORP Messages	Off	Off, Max. Device Limits, Variable Limits
Temperature Messages	Off	Off, Max. Device Limits, Variable Limits

Menu Selection: Parameter Setting MSU ... Module Memosens ORP



Note: During calibration, the function check (HOLD) operating state is active for the corresponding module channel. The assigned current outputs and relay contacts behave as configured (BASE module).

The calibration and adjustment data is saved in the sensor. This means that Memosens sensors can be centrally calibrated, adjusted, reconditioned, and cleaned in a laboratory away from the measurement location. Only the sensors in the system are replaced on-site by calibrated / adjusted sensors.


Calibration: Detecting deviations without readjustment **Adjustment:** Detecting deviations with readjustment

Note: Make sure that the temperature of the sensor and the temperature of the buffer solution are not too far apart. The ideal temperature is 25 $^{\circ}$ C/77 $^{\circ}$ F.

Adjustment


Adjustment means that the values determined by a calibration are applied to the sensor. The zero point and slope values determined during calibration are entered in the adjustment record:


Diagnostics MSU ... Module Memosens ORP Cal/Adj Record
These values are not used to calculate the process variables until the calibration
has been terminated with an adjustment.



Function Check (HOLD) During Calibration / Adjustment

Behavior of the signal and switching outputs during calibration / adjustment





Calibration Mode: ORP Data Entry

Calibration by entering the ORP offset of a premeasured sensor.

Calibration Procedure

When you access calibration, the calibration values selected in the calibration presettings are used. These settings can still be changed in the Calibration menu. If you do not want to perform calibration, use the left softkey to go back one level and abort the calibration.

- 01. Select "ORP Data Entry" calibration mode and press *enter* to confirm.
- 02. Remove the sensor and install the premeasured sensor.
- 03. Continue with the right softkey: Next.
- 04. Enter the ORP offset value.
- √ The calibration record is displayed. By pressing the *right softkey: Adjust*,
 the calculated calibration values are used to calculate the measured values in
 Protos and also saved in the Memosens sensor.

Calibration Mode: ORP Adjustment

The sensor is immersed in an ORP buffer solution for ORP adjustment. Protos displays the measured temperature and the ORP. The temperature-corrected buffer values must then be manually entered. For this purpose, take the buffer value that goes with the displayed temperature from the buffer table (e.g., on the bottle). Intermediate values must be interpolated.

All calibration data is converted to a reference temperature of 25 °C/77 °F.

Calibration Procedure

When you access calibration, the calibration values selected in the calibration presettings are used. These settings can still be changed in the Calibration menu. If you do not want to perform calibration, use the left softkey to go back one level and abort the calibration.


- 01. Select "ORP Adjustment" calibration mode and press enter to confirm.
- 02. Continue with the right softkey: Next.
- 03. Take the sensor out of the medium and rinse it well in deionized water.

CAUTION! Risk of electrostatic discharge. Do not wipe the sensor or dab it dry.

- 04. Immerse the sensor in the ORP buffer solution and wait for the ORP measured value to stabilize.
- 05. Start calibration with the *right softkey: Next*.
- ✓ When the drift check is complete, the measured temperature and ORP are displayed.

The time for the measuring voltage to stabilize can be shortened with the *left softkey: Exit* (without drift check: reduced accuracy of calibration values). The response time indicates how long the sensor needs until the measuring voltage is stable. If the sensor voltage or the measured temperature fluctuate considerably, the calibration procedure is aborted after around 2 minutes. In this case, calibration needs to be restarted.

06. Enter the ORP setpoint (printed on bottle) of the buffer solution in the Calibration Mode submenu > ORP Adjustment > and press *enter* to confirm.



- 07. Exit calibration with the right softkey: Next.
- √ The calibration record is displayed. By pressing the *right softkey: Adjust*,
 the calculated calibration values are used to calculate the measured values in
 Protos and also saved in the Memosens sensor.

Calibration Mode: ORP Check

The sensor is immersed in a solution with a known ORP value for the ORP check. The test period and the permissible test difference are specified during parameter setting:

Parameter Setting ▶ MSU ... Module ▶ Memosens ORP ▶ Cal Presettings

Calibration Procedure

- 01. Select "ORP Check" calibration mode and press enter to confirm.
- 02. Take the sensor out of the medium and rinse it well in deionized water.

CAUTION! Risk of electrostatic discharge. Do not wipe the sensor or dab it dry.

- 03. Immerse the sensor in the ORP solution and wait for the ORP measured value to stabilize.
- 04. Start the ORP check with the *right softkey: Next*.
- ✓ When the drift check is complete, the measured temperature and ORP are displayed.
- ✓ If the specified test difference was not exceeded, the message "ORP Check Successful" appears.

If the specified test difference was exceeded, the message "ORP Check Unsuccessful" appears.

05. ORP adjustment should be carried out in the event of an unsuccessful ORP check.

Calibration Mode: Temperature

This function is used to adjust the individual temperature detector tolerances or cable lengths for the purpose of increasing the accuracy of the temperature measurement. The adjustment requires an accurate measurement of the process temperature using a calibrated reference thermometer. The measurement error of the reference thermometer should be less than 0.1 K. Adjustment without an accurate measurement may result in falsification of the displayed measured value.

Calibration Procedure

- 01. Select "Temperature" calibration mode and press enter to confirm.
- 02. Enter the measured process temperature and press *enter* to confirm.
 √ The temperature offset is displayed.
- 03. Adjust the temperature detector with the *right softkey: Save*.

The current adjustment and temperature offset data can be displayed in the Diagnostics menu:


▶ MSU ... Module ▶ Memosens ORP ▶ Temp. Offset Log

ORP Maintenance Functions


Note: Function check (HOLD) active

The current outputs and relay contacts behave in accordance with the parameter settings. Since the device is in function check (HOLD) mode, certain media can be used to validate the sensor and check the measured values without affecting the signal outputs.

To end the function check, return to measuring mode.



ORP Diagnostic Functions



The diagnostic functions are matched to NAMUR recommendation NE 107. For a detailed description of general diagnostic functions, see the basic unit user manual.

Overview of ORP Diagnostic Functions

In diagnostics mode, you can access the following submenus without interrupting the measurement:

Diagnostics MSU ... Module:

Module Diagnostics Protos periodically performs a self-test in the background.

The results can be displayed here.

Diagnostics MSU ... Module Memosens ORP:

Diagnostics / Mod in Module / Memosens on .		
Submenus	Description	
Sensor Information	The Sensor Information submenu shows data from the currently connected Memosens sensor, e.g., manufacturer, order no., serial no., firmware and hardware version, last calibration, operating time.	
Sensor Monitor	The raw measured values are displayed in the sensor monitor for diagnostic purposes.	
ORP Cal/Adj Record	The calibration/adjustment record shows the data from the last calibration/adjustment performed on the currently connected sensor.	

ORP Diagnostic Functions

Submenus	Description
Temp. Offset Log	The temp. offset log shows the data from the last temperature equalization performed on the currently connected sensor.
Sensor Wear Monitor	The sensor wear monitor shows the sensor operating time and the maximum temperature during the operating time, as well as wear and the forecast remaining time.

pH, ORP Messages

igotimes Failure igotimes Out of Specification igotimes Maintenance Required

See also "Decommissioning"; p. 159

No.	Message Type	Message / Notes	
P008	\otimes	Factory Settings: Switch device off (approx. 10 s). If the message persists, send in the device.	
P009	\otimes	Firmware Error: Switch device off (approx. 10 s). Reload the firmware. If the message persists, send in the device.	
P010	\otimes	pH Range: Measuring range exceeded. Possible causes: Sensor not or incorrectly connected, sensor faulty, cable faulty, incorrect temperature detector selected, temperature detector faulty.	
P011	\otimes	pH LO_LO: Value below configured monitoring limit	
P012	Â	pH LO: Value below configured monitoring limit	
P013	Â	pH HI: Value above configured monitoring limit	
P014	\otimes	pH HI_HI: Value above configured monitoring limit	
P015	\otimes	Temperature Range: Measuring range exceeded. Possible causes: Sensor incorrectly or not connected, cable faulty.	
P016	\otimes	Temperature LO_LO: Value below configured monitoring limit	
P017	À	Temperature LO: Value below configured monitoring limit	
P018	♠	Temperature HI: Value above configured monitoring limit	
P019	\otimes	Temperature HI_HI: Value above configured monitoring limit	
P020	\otimes	ORP Range: Measuring range exceeded. Possible causes: no ORP sensor connected, sensor incorrectly connected, sensor faulty, cable faulty.	
P021	\otimes	ORP LO_LO: Value below configured monitoring limit	
P022	A	ORP LO: Value below configured monitoring limit	
P023	À	ORP HI: Value above configured monitoring limit	
P024	\otimes	ORP HI_HI: Value above configured monitoring limit	
P025	À	rH Range: Measuring range exceeded. Possible causes: no pH/ORP combo sensor connected, sensor incorrectly connected, cable faulty.	
P026	\otimes	rH LO_LO: Value below configured monitoring limit	


No.	Message Type	Message / Notes
P027	Æ	rH LO: Value below configured monitoring limit
P028	À	rH HI: Value above configured monitoring limit
P029	\otimes	rH HI_HI: Value above configured monitoring limit
P045	\otimes	pH Voltage Range: Measuring range exceeded. Possible causes: Sensor not or incorrectly connected, sensor faulty, cable faulty.
P046	\otimes	mV LO_LO: Value below configured monitoring limit
P047	À	mV LO: Value below configured monitoring limit
P048	À	mV HI: Value above configured monitoring limit
P049	\otimes	mV HI_HI: Value above configured monitoring limit
P060	⊗,�	Sad Sensoface: Slope
P061	⊗,�	Sad Sensoface: Zero point
P062	User-defined	Sad Sensoface: Reference impedance
P063	User-defined	Sad Sensoface: Glass impedance
P064	User-defined	Sad Sensoface: Response time
P065	⇔	Sad Sensoface: Calibration timer
P069	⇔	Sad Sensoface: Calimatic (zero/slope)
P070	User-defined	Sad Sensoface: Sensor wear
P071	⇔	Sad Sensoface: ISFET leakage current
P072	⇔	Sad Sensoface: ISFET operating point
P074	⇔	Sad Sensoface: ORP zero offset
P090	⇔	Error in Buffer Table
P110	User-defined	CIP Counter: Configured number of CIP cycles exceeded: As required, calibrate/adjust sensor or replace.
P111	User-defined	SIP Counter: Configured number of SIP cycles exceeded: As required, calibrate/adjust sensor or replace.
P113	User-defined	Sensor Operating Time: Replace the sensor.
P120	\otimes	Wrong Sensor (Sensor Check)
P121	\otimes	Sensor Error (Factory Settings): Replace the sensor.
P122	⇔	Sensor Memory (Cal Data): The calibration data is defective: Recalibrate/readjust sensor.
P123	⇔	New Sensor, Adjustment Required

No.	Message Type	Message / Notes
P124	\(Sensor Date: The sensor data is implausible. Check and, as necessary, adjust the configuration.
P130	Info	SIP Cycle Counted
P131	Info	CIP Cycle Counted
P200	\otimes	Noise Level at pH Input
P201	\(Cal Temperature: The calibration temperature is impermissible: Check the calibration temperature. Note the information in the Calibration chapter.
P202	Info	Cal: Buffer Unknown
P203	Info	Cal: Identical Buffers
P204	Info	Cal: Buffers Interchanged
P205	Info	Cal: Sensor Unstable: The drift criterion was not adhered to during calibration. Possible causes: improper calibration, sensor cable/connection faulty, sensor worn. Check the sensor and calibration and repeat as required. Otherwise, replace the sensor.
P206	\(\bar{\phi} \)	Cal: Slope: Slope out of permissible range: Repeat calibration/adjustment or replace sensor.
P207	€	Cal: Zero Point: Zero point out of permissible range: Repeat calibration/adjustment or replace sensor.
P208	\otimes	Cal: Sensor Failure (ORP Check): Replace the sensor.
P254	Info	Module Reset

No.	Message Type	pH/pH Calculation Block Messages
A010	\otimes	pH Diff Range: Max/min device limits exceeded: - Check both pH values Check sensor/cable connections.
A011	\otimes	pH Diff LO_LO: Value below configured monitoring limit
A012	A	pH Diff LO: Value below configured monitoring limit
A013	À	pH Diff HI: Value above configured monitoring limit
A014	\otimes	pH Diff HI_HI: Value above configured monitoring limit
A015	\otimes	Temperature Diff Range: Max/min device limits exceeded: - Check both temperature values Check sensor/cable connections.
A016	\otimes	Temperature Diff LO_LO: Value below configured monitoring limit
A017	À	Temperature Diff LO: Value below configured monitoring limit
A018	A	Temperature Diff HI: Value above configured monitoring limit
A019	\otimes	Temperature Diff HI_HI: Value above configured monitoring limit
A020	\otimes	ORP Diff Range: Max/min device limits exceeded: - Check both ORP values Check sensor/cable connections.
A021	\otimes	ORP Diff LO_LO: Value below configured monitoring limit
A022	À	ORP Diff LO: Value below configured monitoring limit
A023	A	ORP Diff HI: Value above configured monitoring limit
A024	\otimes	ORP Diff HI_HI: Value above configured monitoring limit
A045	\otimes	pH Voltage Diff Range: Max/min device limits exceeded: - Check both pH voltage values Check sensor/cable connections.
A046	\otimes	pH Voltage Diff LO_LO: Value below configured monitoring limit
A047	À	pH Voltage Diff LO: Value below configured monitoring limit
A048	A	pH Voltage Diff HI: Value above configured monitoring limit
A049	\otimes	pH Voltage Diff HI_HI: Value above configured monitoring limit
A200	⇔	Calculation Block Configuration

Note: Function check (HOLD) active


Note: The oxygen measurement requires TAN option FW4400-015: Oxygen Measurement in the Saturation and ppb Range




Menu Selection: Parameter Setting ▶ MSU ... Module ▶ Memosens Oxy:

Parameter	Default	Selection / Range
In most Filters		
Input Filter	T	
Pulse Suppression	Weak	Suppression of interference pulses: Off, Weak, Medium, Strong
Sensor Data		
Measure in	Liquids	Liquids, Gases
Relative Humidity	50.0 %	If measuring in gases
Sensoface	Off	Off, On
Sensor Monitoring Details (see page 68)		Slope, Zero Point, Sensocheck, Response Time, Sensor Operating Time, Sensor Wear, SIP Counter
Cal Presettings	<u>.</u>	
Calibration Mode	In Air	Presetting for calibration mode: In Air, In Water, Data Entry, Product, Zero Point, Temperature
Product Calibration	Sat. %Air	If "Product" is selected: Sat. %Air, Conc. (Liquid), Partial Pressure
Calibration Timer		
Monitoring	Off	Off, Auto: 720 h, Individual
Pressure Correction	<u>.</u>	
Ext. Pressure Transmitter		
Pressure Transmitter	Absolute	
Current Input	0 20 mA	0 20 mA / 4 20 mA
Start 0(4) mA	0000 mbar	xxxx mbar
End 20 mA	2000 mbar	xxxx mbar
Pressure During Meas.		
Detection	Manual	Manual (default 1013 mbar), External, AO 1 if PROFINET is active ¹⁾
Pressure During Cal		
Detection	Manual	Manual (default 1013 mbar), External, AO 1 if PROFINET is active 1)
Salinity Correction		
Entry	Salinity	Salinity, Chlorinity, Conductivity (0.00 g/kg or 0.000 µS/cm depending on selection) The calculated salinity is displayed if Chlorinity/ Conductivity are selected.

Parameter	Default	Selection / Range
Messages		
Saturation %Air Messages 1)	Off	Off, Max. Device Limits, Variable Limits
Saturation %O ₂ Messages 1)	Off	Off, Max. Device Limits, Variable Limits
Concentration Messages	Off	Off, Max. Device Limits, Variable Limits
Partial Pressure Messages	Off	Off, Max. Device Limits, Variable Limits
Temperature Messages	Off	Off, Max. Device Limits, Variable Limits





Note: During calibration, the function check (HOLD) operating state is active for the corresponding module channel. The assigned current outputs and relay contacts behave as configured (BASE module).


The calibration data is saved in the Memosens sensor. This means that Memosens sensors can be cleaned, reconditioned, calibrated, and adjusted away from the measurement location, e.g., in a laboratory. Sensors in the system are replaced on-site by adjusted sensors.

Calibration: Detecting deviations without readjustment **Adjustment:** Detecting deviations with readjustment

Adjustment

Adjustment means that the values determined by a calibration are applied to the sensor. The zero point and slope values determined during calibration are entered in the adjustment record:

Diagnostics MSU ... Module Memosens Oxy Cal/Adj Record
These values are not used to calculate the process variables until the calibration
has been terminated with an adjustment.



Oxygen Calibration/Adjustment Explanations

Every oxygen sensor has an individual slope and an individual zero point. Both values are altered, for example, by aging. In order to achieve sufficient accuracy in the oxygen measurement, the meter must be regularly adjusted to reflect the sensor data (adjustment).

The "slope" is the sensor current value with atmospheric oxygen saturation, 25 °C/77 °F and 1013 mbar/14.69 psi: nA/100 %. Only the "nA" measurement symbol appears on the display.

This is technically not a "slope" but rather a calibration point. The value is provided with the intention of enabling the sensor to be compared with the datasheet values.

If the electrolyte, the membrane body, or both are replaced during maintenance of amperometric sensors, this change must be confirmed manually in the Maintenance menu:

Maintenance ▶ MSU ... Module ▶ Memosens Oxy ▶ Membrane Body Replacement

Calibration is required after each membrane body replacement. This entry impacts on the accuracy of the calibration.

Recommendations for Calibration

For best performance, you should always calibrate in air. Compared to water, air is a calibration medium that is easy to handle, stable, and thus safe. In the most cases, however, the sensor must be removed for a calibration in air. In some processes, the sensor cannot be removed for calibration. Here, calibration must be performed directly in the process medium (e.g., with aeration).

Common Combination: Process Variable / Calibration Mode

Measurement	Calibration	
Saturation	Water	
Concentration	Air	

Note: Make sure that the temperature of the sensor and the temperature of the calibration solution are not too far apart. The ideal temperature is 25 °C/77 °F. If there is a temperature difference between the calibration medium and the measured medium, you must keep the sensor in the respective medium for several minutes before and after calibration in order to achieve stable measured values.


The type of calibration pressure detection is preset during parameter setting.


Parameter Setting MSU ... Module Memosens Oxy Pressure Correction Pressure During Cal

Note: Amperometric sensors must be sufficiently polarized prior to calibration/ adjustment. Follow the information on the sensor in the sensor user manual to ensure that the calibration is neither falsified nor unstable.

Function Check (HOLD) During Calibration / Adjustment

Behavior of the signal and switching outputs during calibration / adjustment





Calibration Mode: In Air

The slope is corrected using the saturation value (100 % Air), similar to the air saturation of water. Since this analogy only applies to water-vapor saturated air (100 % relative humidity) and in many cases the calibration air is less humid, the relative humidity of the calibration air must also be specified. If you do not know the exact value of the relative humidity of the calibration air, you can take the following reference values for a sufficiently precise calibration:

- Ambient air: 50 % rel. humidity (average)
- Bottled gas (synthetic air): 0 % rel. humidity

Calibration Procedure

Notes: The sensor membrane must be dry. Be sure to keep temperature and pressure constant during calibration. If there is a temperature difference between calibration and measured media, the sensor requires some equalization time before and after calibration.

When you access calibration, the calibration values selected in the calibration presettings are used. These settings can still be changed in the Calibration menu. If you do not want to perform calibration, use the left softkey to go back one level and abort the calibration.

- 01. Select "In Air" Calibration Mode and press *enter* to confirm.
- 02. Take the sensor out of the medium and clean it.
- 03. Carefully dab the membrane dry with a paper tissue.
- 04. Expose the sensor to air with a known water vapor saturation and press *enter* to confirm.
 - √ Display of selected calibration medium (air)
- 05. Enter relative humidity, e.g.: Ambient air: 50 %, Bottled gas: 0 %
- 06. Enter Cal Pressure: Enter the calibration pressure if "Manual" was configured.
- 07. Start with the right softkey: Next
 - ✓ Drift check is carried out.

The following parameters are displayed: Sensor Current, Calibration Pressure, and Response Time.

08. Exit calibration with the *right softkey: Next*.

✓ The calibration record is displayed. By pressing the *right softkey: Adjust*, the calculated calibration values are used to calculate the measured values in Protos and also saved in the Memosens sensor.

Calibration Mode: In Water

The slope is corrected using the saturation value (100 %) related to saturation with air.

Calibration Procedure

Note: Ensure sufficient sensor incident flow (see the specifications of the oxygen sensor.)

The calibration medium must be in equilibrium with air. Oxygen exchange between water and air is very slow. Therefore, it takes a relatively long time until water is saturated with atmospheric oxygen. If there is a temperature difference between calibration medium and measured medium, you must keep the sensor in the respective medium for several minutes before and after calibration.

When you access calibration, the calibration values selected in the calibration presettings are used. These settings can still be changed in the Calibration menu. If you do not want to perform calibration, use the left softkey to go back one level and abort the calibration.

- 01. Select "In Water" Calibration Mode and press enter to confirm.
- 02. Take the sensor out of the medium and clean it.
- 03. Carefully dab the membrane dry with a paper tissue.
- 04. Expose the sensor to the calibration medium (air-saturated water), ensure sufficient incident flow, and press *enter* to confirm.
 - √ Display of selected calibration medium (air-saturated water)
- 05. Enter Cal Pressure: Enter the calibration pressure if "Manual" was configured.
- 06. Start with the *right softkey: Next*.
 - ✓ Drift check is carried out.

The following parameters are displayed: Sensor Current, Calibration Pressure, and Response Time.

The time for the sensor signal to stabilize can be shortened with the *left soft-key: Exit* (without drift check: reduced accuracy of calibration values). The response time indicates how long the sensor needs until the sensor signal is stable. If the sensor signal or the measured temperature fluctuate considerably

or the sensor is inadequately polarized, the calibration procedure is aborted after around 2 minutes. In this case, calibration needs to be restarted. Return the sensor to the process once this has been successfully completed. Make sure that the temperature of the sensor and the temperature of the calibration solution are not too far apart. The ideal temperature is 25 °C/77 °F.

07. Exit calibration with the *right softkey: Next*.

√ The calibration record is displayed. By pressing the *right softkey: Adjust*, the calculated calibration values are used to calculate the measured values in Protos and also saved in the Memosens sensor.

Calibration Mode: Data Entry

Entry of values for slope and zero point of the sensor, related to 25 °C/77 °F, 1013 mbar/14.69 psi.

Slope = sensor current at 100 % atmospheric oxygen, 25 °C/77 °F, 1013 mbar/14.69 psi

Calibration Procedure

When you access calibration, the calibration values selected in the calibration presettings are used. These settings can still be changed in the Calibration menu. If you do not want to perform calibration, use the left softkey to go back one level and abort the calibration.

- 01. Select "Data Entry" calibration mode and press enter to confirm.
- 02. Remove the sensor and install the premeasured sensor.
- 03. Continue with the *right softkey: Next*.
- 04. Enter the zero point and slope measured values.
- ✓ By pressing the *right softkey: Adjust*, the calculated calibration values are used to calculate the measured values in Protos and also saved in the Memosens sensor.

Calibration Mode: Product

If the sensor cannot be removed – e.g., for sterility reasons – its slope can be determined by "sampling". The current measured value is saved in the device for this purpose. A sample is taken directly afterward at the measuring point. The reference value is entered in the device. Protos calculates the sensor's correction values from the difference between the measured and reference values, and corrects the zero point in the event of small saturation values, the slope in the event of large values.

Calibration Procedure

When you access calibration, the calibration values selected in the calibration presettings are used. These settings can still be changed in the Calibration menu. If you do not want to perform calibration, use the left softkey to go back one level and abort the calibration.

- 01. Select "Product" calibration mode and press enter to confirm.
 - ✓ Process variable Saturation, Concentration, and Partial Pressure configured as in "Cal Presettings".
- 02. Change the process variable, if required.
- 03. Prepare for sampling.


NOTICE! Measure the reference value at temperature and pressure conditions similar to those of the process.

04. Start with the right softkey: Next.

Product calibration is performed in 2 steps.

Step 1:

- 05. Take sample.
 - √ The measured value and temperature at the time of sampling are displayed.
- 06. Save with the *right softkey: Save*.



Step 2: Lab value has been measured.

- 01. Open the Product Calibration menu again.
- 02. Enter reference value ("Lab value").
- 03. Confirm with the *right softkey: Next* or repeat calibration with the *left soft-key: Cancel*.
- ✓ By pressing the *right softkey: Adjust*, the calibration values obtained during calibration for calculation of the process variables are applied to the device and saved in the sensor.

Exception: Sample value can be determined and entered on site:

Left softkey: Entry.

Calibration Mode: Zero Point

For trace measurements below 500 ppb, the zero point should be calibrated. If a zero correction is performed, the sensor should remain in the calibration medium (media containing CO_2 : at least 120 min) for at least 10 to 60 minutes to obtain stable, non-drifting values. During zero correction, a drift check is not performed.

Calibration Procedure

When you access calibration, the calibration values selected in the calibration presettings are used. These settings can still be changed in the Calibration menu. If you do not want to perform calibration, use the left softkey to go back one level and abort the calibration.

- 01. Select "Zero Point" Calibration Mode and press enter to confirm.
- 02. Press the *right softkey: Next*.
 - ✓ Zero point correction is carried out. The measured sensor current is displayed.
- 03. Enter the input current for the zero point.
- 04. Press the right softkey: Next.
- √ The calibration record is displayed. By pressing the *right softkey: Adjust*, the calculated calibration values are used to calculate the measured values in Protos and also saved in the Memosens sensor.

Calibration Mode: Temperature

This function is used to adjust the individual temperature detector tolerances or cable lengths for the purpose of increasing the accuracy of the temperature measurement. The adjustment requires an accurate measurement of the process temperature using a calibrated reference thermometer. The measurement error of the reference thermometer should be less than 0.1 K. Adjustment without an accurate measurement may result in falsification of the displayed measured value.

Calibration Procedure

- 01. Select "Temperature" calibration mode and press enter to confirm.
- 02. Enter the measured process temperature and press *enter* to confirm
 √ The temperature offset is displayed..
- 03. Adjust the temperature detector with the right softkey: Save.

The current adjustment and temperature offset data can be displayed in the Diagnostics menu:


▶ MSU ... Module ▶ Memosens Oxy ▶ Temp. Offset Log

Oxy Maintenance Functions

Note: Function check (HOLD) active

The current outputs and relay contacts behave in accordance with the parameter settings. Since the device is in function check (HOLD) mode, certain media can be used to validate the sensor and check the measured values without affecting the signal outputs.

To end the function check, return to measuring mode.



From the measuring mode: menu key: Menu Selection.

Select Maintenance (maint) using the arrow keys and confirm by pressing


Passcode (factory setting): 2958 Then select the module and sensor.

During maintenance, the sensor monitor allows validation of the sensor by immersing it in a known solution, for example, and checking the values

Membrane Body Replacement

If the electrolyte or the membrane body of a Memosens oxygen sensor are replaced during maintenance of a sensor, this must be confirmed with "Yes" in the Maintenance menu: After confirming the check, the calibration counter for the membrane body "Membrane Cals" is reset to "0". See the Sensor wear monitor figure on page 87.

Oxy Diagnostic Functions



For a detailed description of general diagnostic functions, see the basic unit user manual.

Overview of Oxy Diagnostic Functions

In diagnostics mode, you can access the following submenus without interrupting the measurement:

Diagnostics MSU ... Module:

Module Diagnostics Protos periodically performs a self-test in the back-

ground. The results can be displayed here.

Diagnostics MSU ... Module Memosens Oxy:

Submenus	Description
Sensor Information	The Sensor Information submenu shows data from the currently connected Memosens sensor, e.g., manufacturer, order no., serial no., firmware and hardware version, last calibration, operating time.
Sensor Monitor	The raw measured values are displayed in the sensor monitor for diagnostic purposes.


Oxy Diagnostic Functions


Submenus

Oxy Sensor Diagram

Description

The sensor diagram clearly indicates the status of the parameters in the connected sensor, including the calibration timer. Inactive parameters are shown in gray and set to 100 % (e.g., disabled response time). The parameter values should lie between the outer (100 %) and inner (50 %) polygon. A warning signal flashes if a value drops below the inner polygon (< 50 %).





Critical range - "inner circle": Value outside tolerance


The tolerance limits (radius of the "inner circle") can be individually adjusted. See Parameter Setting • Sensor Data • Sensor Monitoring Details.

The calibration/adjustment record shows the data from the last calibration/adjustment performed on the currently connected sensor. The temp. offset log shows the data from the last temperature equalization performed on the currently connected sensor.

Oxy Cal/Adj Record

Temp. Offset Log

Sensor Wear Monitor



The sensor wear monitor shows the sensor operating time and the maximum temperature during the operating time, as well as wear and the forecast remaining time.

Oxy Messages

⊗ Failure ⚠ Out of Specification ⇔ Maintenance Required

See also "Decommissioning"; p. 159

No.	Message Type	Message / Notes					
D008	\otimes	Factory Settings: Switch device off (approx. 10 s). If the message persists, send in the device.					
D009	\otimes	Firmware Error: Switch device off (approx. 10 s). Reload the firmware. If the message persists, send in the device.					
D010	\otimes	Saturation %Air Range: Max/min limits of range exceeded or sensor not/incorrectly connected, cable faulty.					
D011	\otimes	Saturation %Air LO_LO: Value below configured monitoring limit					
D012	À	Saturation %Air LO: Value below configured monitoring limit					
D013	҈Ѧ	Saturation %Air HI: Value above configured monitoring limit					
D014	\otimes	Saturation %Air HI_HI: Value above configured monitoring limit					
D015	\otimes	Temperature Range: Max/min limits of range exceeded or sensor not/incorrectly connected, cable faulty.					
D016	\otimes	Temperature LO_LO: Value below configured monitoring limit					
D017	Æ	Temperature LO: Value below configured monitoring limit					
D018	À	Temperature HI: Value above configured monitoring limit					
D019	\otimes	Temperature HI_HI: Value above configured monitoring limit					
D020	\otimes	Concentration Range: Max/min limits of range exceeded or sensor not/incorrectly connected, cable faulty.					
D021	\otimes	Concentration LO_LO: Value below configured monitoring limit					
D022	҈Ѧ	Concentration LO: Value below configured monitoring limit					
D023	A	Concentration HI: Value above configured monitoring limit					
D024	\otimes	Concentration HI_HI: Value above configured monitoring limit					
D025	\otimes	Partial Pressure Range: Max/min limits of range exceeded or sensor not/incorrectly connected, cable faulty.					
D026	\otimes	Part. Pressure LO_LO: Value below configured monitoring limit					
D027	À	Part. Pressure LO: Value below configured monitoring limit					


No.	Message Type	Message / Notes				
D028	A	Part. pressure HI: Value above configured monitoring limit				
D029	\otimes	Part. Pressure HI_HI: Value above configured monitoring limit				
D045	\otimes	Saturation $\%O2$ Range: Max/min limits of range exceeded or sensor not/inco rectly connected, cable faulty.				
D046	\otimes	Saturation %O2 LO_LO: Value below configured monitoring limit				
D047	À	Saturation %O2 LO: Value below configured monitoring limit				
D048	À	Saturation %O2 HI: Value above configured monitoring limit				
D049	\otimes	Saturation %O2 HI_HI: Value above configured monitoring limit				
D060	⊗,�	Sad Sensoface: Slope - Readjust sensor Check/refill electrolyte Replace sensor.				
D061	⊗,�	Sad Sensoface: Zero Point - Readjust sensor Check/refill electrolyte Replace sensor.				
D062	User-defined	Sad Sensoface: Sensocheck - Readjust sensor Replace sensor.				
D063	•	Sad Sensoface: Response Time - Check/refill electrolyte Replace sensor.				
D064	•	Sad Sensoface: Calibration Timer				
D070	User-defined	Sad Sensoface: Sensor Wear Sensor is Worn (100 %): - Readjust sensor Check/refill electrolyte Replace sensor.				
D080	€	Sensor Current Range - Check polarization voltage: Parameter Setting • MSU Module • Memosens Oxy • Sensor Data - Refill electrolyte Recalibrate/readjust.				
D113	User-defined	Sensor Operating Time: Replace the sensor.				
D120	\otimes	Wrong Sensor (Sensor Check)				
D121	\otimes	Sensor Error (Factory Settings): Replace the sensor.				
D122	€	Sensor Memory (Cal Data): The calibration data is defective: Recalibrate/readjust sensor.				

No.	Message Type	Message / Notes				
D123	⇔	New Sensor, Adjustment Required				
D124	⇔	Sensor Date: The sensor data is implausible. Check and, as necessary, adjust the configuration.				
D200	•	Temp O2 Conc/Sat: The temperature is outside the valid range for oxygen concentration/saturation.				
D201	€	Cal Temperature: The calibration temperature is impermissible: Check the calibration temperature. Note the information in the Calibration chapter.				
D203	Info	Cal: Identical media				
D204	Info	Cal: Media interchanged				
D205	Info	Cal: Sensor Unstable: The drift criterion was not adhered to during calibration. Possible causes: improper calibration, sensor cable/connection faulty, sensor worn. Check the sensor and calibration and repeat as required. Otherwise, replace the sensor.				
D254	Info	Module Reset				

No.	Message Type	pe Oxy/Oxy Calculation Block Messages				
H010	\otimes	Saturation %AIR Diff Range: Max/min device limits exceeded: - Check both saturation values Check sensor/cable connections.				
H011	\otimes	Saturation %AIR Diff LO_LO: Value below configured monitoring limit				
H012	Â	Saturation %AIR Diff LO: Value below configured monitoring limit				
H013	Â	Saturation %AIR -Diff HI: Value above configured monitoring limit				
H014	\otimes	Saturation %AIR Diff HI_HI: Value above configured monitoring limit				
H015	\otimes	Temperature Diff Range: Max/min device limits exceeded: - Check both temperature values Check sensor/cable connections.				
H016	\otimes	Temperature Diff LO_LO: Value below configured monitoring limit				
H017	À	Temperature Diff LO: Value below configured monitoring limit				
H018	҈Ѧ	Temperature Diff HI: Value above configured monitoring limit				
H019	\otimes	Temperature Diff HI_HI: Value above configured monitoring limit				
H020	\otimes	Conc. (Liquid) Diff Range: Max/min device limits exceeded: - Check both concentration values Check sensor/cable connections.				
H021	\otimes	Conc. (Liquid) Diff LO_LO: Value below configured monitoring limit				
H022	Â	Conc. (Liquid) Diff LO: Value below configured monitoring limit				
H023	҈Ѧ	Conc. (Liquid) Diff. HI: Value above configured monitoring limit				
H024	\otimes	Conc. (Liquid) Diff HI_HI: Value above configured monitoring limit				
H045	\otimes	Saturation %O2 Diff Range: Max/min device limits exceeded: - Check both saturation values Check sensor/cable connections.				
H046	\otimes	Saturation %O2 Diff LO_LO: Value below configured monitoring limit				
H047	A	Saturation %O2 Diff LO: Value below configured monitoring limit				
H048	A	Saturation %O2 Diff HI: Value above configured monitoring limit				
H049	\otimes	Saturation %O2 Diff HI_HI: Value above configured monitoring limit				
H090	\otimes	Conc. (Gas) Diff Range (Measurement in Gases): Max/min device limits exceeded: - Check both concentration values Check sensor/cable connections.				


No.	Message Type	Oxy/Oxy Calculation Block Messages			
H091	\otimes	Conc. (Gas) Diff LO_LO: Value below configured monitoring limit.			
H092	Â	Conc. (Gas) Diff LO: Value below configured monitoring limit			
H093	À	Conc. (Gas) Diff HI: Value above configured monitoring limit			
H094	\otimes	Conc. (Gas) Diff HI_HI: Value above configured monitoring limit			
H200	⇔	Calculation Block Configuration			

Note: Function check (HOLD) active



Menu Selection: Parameter Setting ▶ MSU ... Module ▶ Memosens Cond:

Parameter	Default	Selection / Range		
Input Filter				
Pulse suppression	Off	Enable/disable suppression of interference pulses.		
Sensor Data				
Sensoface	On	On, Off		
Sensor Monitoring Details (see page 95) • Cell Constant • Sensocheck • Sensor Operating Time	Auto Off Off	Auto, Individual Off, On Off, Individual (max. 9999 d)		
• CIP Counter	Off	Off, Individual		
SIP Counter	Off	Off, Individual		
Cal Presettings				
Calibration Mode	Automatic	Presetting for calibration mode: Automatic, Manual, Product, Installation Factor (4-Electrode Sensor), Data Entry, Temperature		
Cal Solution	NaCl Sat	In Automatic mode: Selection of calibration solution NaCl 0.01 m: 1183 µS/cm NaCl 0.1 m: 10.683 mS/cm NaCl Sat 251.3 mS/cm KCl 0.01 m: 1413 µS/cm KCl 0.1 m: 12.88 mS/cm KCl 1 m: 111.80 mS/cm		
Product Calibration	Conductivity	Conductivity, concentration (with TAN option FW4400-009)		
Conductivity	Without TC	Without TC, With TC		
TC Process Medium (see page	97)			
Temperature Compensation	Off	Off, Linear, EN27888, Ultrapure Water (with TAN option FW4400-008)		
Concentration (see page 99)				
Concentration	Off	Off, On		
TDS (see page 97)				
TDS Function	Off	Off, On (preset 1.00)		
USP (see page 98)				
USP Function	Off	Off, On		
Messages (see page 101)				
Messages	Temperature: Max. Device Limits	Conductivity, Resistivity, Concentration, Temperature, Salinity, TDS Can be adjusted for all types of monitoring: Off, Max. Device Limits, Variable Limits)		



CIP Counter/SIP Counter

CIP/SIP cycles are used to clean or sterilize process-wetted parts in the process. Depending on the application, either one chemical (alkaline solution, water) or several chemicals (alkaline solution, water, acidic solution, water) are used.

- CIP temperature > 55 °C / 131 °F
- SIP temperature > 115 °C / 239 °F


The cleaning (cleaning in place) and sterilization (sterilization in place) cycles are counted to measure the load on the sensor, e.g., in biotechnology applications.

Note: If measurements are generally taken at high temperatures (> 55 $^{\circ}$ C / 131 $^{\circ}$ F), the counters should be switched off.

When a CIP/SIP counter is switched on, a maximum number of cycles can be entered. A message can be set to indicate when a counter has reached a specified value.

Note: A CIP or SIP cycle is only entered into the logbook 2 hours after the start to ensure that the cycle is complete.

Presettings for Calibration



Temperature Compensation of Process Medium

The following are available for temperature compensation:

- Off
- Linear (enter temperature coefficient TC)
- EN 27888 (natural waters)
- Ultrapure Water (with different trace impurities)

Trace Impurities in Ultrapure Water (with TAN option FW4400-008)

NaCl	Neutral ultrapure water, for conductivity measurement in water
	processing downstream of gravel bed filter
HCI	Acidic ultrapure water, for conductivity measurement downstream
	of cation filter

NH₃ Ammoniacal ultrapure water

Alkaline ultrapure water

Note: If TC correction for the process medium is enabled, "TC" is shown on the display in measuring mode.

TDS Function

NaOH

TDS (total dissolved solids) = weight of dissolved solids that influence conductivity

The TDS function provides a quick way of determining the evaporation residue of water. A TDS factor must be entered for this purpose.

The factor establishes a simple linear relationship between measured conductivity and evaporation residue. It is dependent on the composition of the medium and must be empirically determined by the user.

USP Function (Cond)

Monitoring Ultrapure Water in the Pharmaceutical Industry

The conductivity of ultrapure water in the pharmaceutical industry can be monitored online in accordance with the guideline "USP" (U.S. Pharmacopeia), Annex 5, Section 645 "Water Conductivity". The conductivity is measured without temperature compensation and is compared with limits. The water is usable without any further testing when the conductivity is below the USP limit.

Configuring the USP Function

The USP value can be configured as a process variable USP% for output (display, current output, limit, measurement recorder)
Settings can be changed in the USP submenu:
Parameter Setting MSU ... Module Memosens Cond USP

Reduced Limit Value: The USP limit can be reduced down to 10 %. **Monitoring:** Select whether and how an exceeded limit is to be displayed.

Off No message, but the parameter is still shown in the

Diagnostics menu.

Failure A failure message is shown in off-limit conditions;

the corresponding NAMUR icon is displayed.

Maintenance A Maint. Required message is shown in off-limit conditions;

the corresponding NAMUR icon is displayed.

USP Function: Specify a Relay Contact

The USP function can also be assigned to a relay contact:

Parameter Setting > BASE Module ... > Contact K... > Usage: USP Output

Display of the USP Function in the Diagnostics Menu

Diagnostics • MSU ... Module • Memosens Cond • USP Function Display of the USP limit, the reduced limit, and conductivity.

Concentration (Cond)

Note: Concentration determination requires the activation of TAN option FW4400-009.

The substance concentration in percent by weight (wt%) is determined for H_2SO_4 , HNO_3 , HCI, NaOH, NaCI, and Oleum. For concentration curves, see page 174.

Conditions for Concentration Determination

The following conditions must be met for a reliable concentration determination:

- For calculation of concentration, the medium to be measured must be a purely binary mixture (e.g., water-hydrochloric acid). Presence of other dissolved substances (e.g., salts) leads to incorrect concentration values.
- In the region of small slopes (e.g., at the range limits) small changes in conductivity can correspond to large changes in concentration. This may lead to an unsteady display of the concentration value.
- As the concentration value is calculated from the measured conductivity and temperature values, accurate temperature measurement is very important. Therefore, you should make sure that conductivity sensor and process medium are in thermal equilibrium.

Settings can be changed in the Concentration submenu:

Parameter Setting MSU ... Module Memosens Cond Concentration

- 01. Concentration: On
- 02. Select the medium:

NaCl (0-28 %), HCl (0-18 %), NaOH (0-24 %), H $_2$ SO $_4$ (0-37 %), HNO $_3$ (0-30 %), H $_2$ SO $_4$ (89-99 %), HCl (22-39 %), HNO $_3$ (35-96 %), H $_2$ SO $_4$ (28-88 %), NaOH (15-50 %), Oleum (12-45 %), table

You can define limits for warning and failure messages for the concentration value:

Parameter Setting ▶ MSU ... Module ▶ Memosens Cond ▶ Messages ▶ Concentration Messages

Concentration (Cond)

Specifying a Concentration Solution for Conductivity Measurement

To specify the customer-specific solution, 5 concentration values A-E are entered in a matrix together with 5 temperature values 1-5. First enter the 5 temperature values, then the corresponding conductivity values for each of the concentrations A-E.

These solutions are then available as "Table" in addition to the permanently set standard solutions.

Settings can be changed in the System Control in the Concentration Table submenu:

Parameter Setting > System Control > Concentration Table:

- 01. Enter temperatures 1 to 5.
- 02. Enter values for concentrations A-E for the respective temperatures.

Notes: The temperature values must be rising (Temp. 1 is the lowest, Temp. 5 the highest temperature).

The concentration values must be rising (Conc. A is the lowest, conc. E the highest concentration).

The table values A1 ... E1, A2 ... E2, etc. must all be rising within the table or all falling.

Points of inflection are not allowed.

Incorrect table entries are indicated by an exclamation point in a red triangle. The table is built up as 5x5 matrix:


	Conc. A	Conc. B	Conc. C	Conc. D	Conc. E
Temp 1	A1	B1	C1	D1	E1
Temp 2	A2	B2	C2	D2	E2
Temp 3	A3	В3	C3	D3	E3
Temp 4	A4	B4	C4	D4	E4
Temp 5	A5	B5	C5	D5	E5

The concentration table is selected in the menu:

Parameter Setting ▶ MSU ... Module ▶ Memosens Cond ▶ Cal Presettings

Calibration Mode: Automatic

Cal Solution: Table



pH Value Calculation (Cond)

Note: Two conductivity sensors required


pH Value Calculation by Means of Dual Conductivity Measurement

When monitoring boiler feed water in power plants, the pH value can be calculated by means of a dual conductivity measurement. For that purpose, the boiler feedwater conductance is measured before and after the cation exchanger. This commonly used method of indirect pH value measurement does not require much maintenance and has the following advantage:

Normal pH measurement in ultrapure water is very critical. Boiler feedwater does not contain many ions. This requires the use of a special electrode, which must be calibrated constantly and the service life of which is generally rather short.

Function

Two conductivity sensors are used to measure the conductivity upstream and downstream of ion exchanger. The concentration of the sodium hydroxide and the pH value are determined from the two calculated conductivity values via a "calculation block" in accordance with the calculation formulas set out below:



Calculating the Concentration of Sodium Hydroxide / the pH Value

$$c(NaOH) = \frac{Cond1- 1/3 Cond2}{243} pH = 11+log[c(NaOH)]$$

For a description of calculation blocks, see the basic unit user manual.

pH Value Calculation (Cond)

Recommended pH Ranges:

 10 ± 0.2 for < 136 bar operating overpressure or 9.5 ± 0.2 for > 136 bar operating overpressure




Figure:

Conditioning of water in natural circulation boilers with sodium hydroxide. Relationship between the pH value and the conductivity measured upstream and downstream of the cation exchanger.

Source: Appendix to VGB guideline for boiler feed water, boiler water, and steam of steam generators above 68 bar permissible operating overpressure (VGB-R 450 L, 1988 edition)

Note: During calibration, the function check (HOLD) operating state is active for the corresponding module channel. The assigned current outputs and relay contacts behave as configured (BASE module).


The calibration data is saved in the Memosens sensor. This means that Memosens sensors can be cleaned, reconditioned, calibrated, and adjusted away from the measurement location, e.g., in a laboratory. Sensors in the system are replaced on-site by adjusted sensors.

Calibration: Detecting deviations without readjustment **Adjustment:** Detecting deviations with readjustment

Adjustment

Adjustment means that the values determined by a calibration are applied to the sensor. The zero point and slope values determined during calibration are entered in the adjustment record:

Diagnostics MSU ... Module Memosens Cond Cal/Adj Record
These values are not used to calculate the process variables until the calibration
has been terminated with an adjustment.



Explanations Regarding Calibration/Adjustment with 2-/4-Electrode Sensors

Each conductivity sensor has an individual cell constant. Depending on the sensor design, the cell constant may vary over a wide range. Because the conductivity value is calculated from the measured conductance and the cell constant, the device must know the cell constant.

During calibration or sensor adjustment, either the known (printed) cell constant of the used conductivity sensor is entered in the device, or it is determined automatically by measuring a calibration solution with known conductivity.

Notes on Calibration

- Use only fresh calibration solutions.
- The used calibration solution must be configured.
- The accuracy of the calibration is crucially dependent on precise acquisition
 of the calibration solution temperature. On the basis of the measured or
 entered temperature, Protos calculates the setpoint of the calibration solution from a stored table.
- Note the response time of the temperature probe.
- To determine the exact cell constant, wait for temperature equalization of the temperature probe and calibration solution before calibration.

Since the cell constant is subject to production-related fluctuations, it is recommended that the removed sensor be calibrated with a calibration solution (e.g., saturated NaCl). The cell constants of the sensors are dependent on the installation geometry – especially in the case of fringe-field sensors:

- If the sensor is installed in a free space (minimum distances exceeded), the cell constant specified in the specifications can be entered directly. "Data Entry" Calibration Mode
- If the installation space is tight (minimum distances are not reached), the sensor must be adjusted in its installed state, as the resulting cell constant has changed. "Product" Calibration Mode

Temperature Compensation during Calibration


The conductivity value of the calibration solution is temperature-dependent. During calibration, the temperature of the calibration solution must therefore be known so that the actual value can be obtained from the conductivity table.


Automatic Temperature Compensation

Protos measures the temperature of the calibration solution using the temperature detector integrated in the Memosens sensor.

Function Check (HOLD) During Calibration / Adjustment

Behavior of the signal and switching outputs during calibration / adjustment





Calibration Mode: Automatic

During automatic calibration, the conductivity sensor is immersed in a standard calibration solution (NaCl or KCl, set during parameter setting in the "Cal Presettings" submenu). On the basis of the measured conductance and temperature, Protos automatically calculates the cell constant. The temperature dependence of the calibration solution is accounted for.

Notes on Calibration

- Use only fresh calibration solutions. The used calibration solution must be configured.
- The accuracy of the calibration is crucially dependent on a precise acquisition of the calibration solution temperature. On the basis of the measured or entered temperature, Protos calculates the setpoint of the calibration solution from a stored table.
- Note the response time of the temperature probe.
- To determine the exact cell constant, wait for temperature equalization of the temperature probe and calibration solution before calibration.
- If the measured conductance or temperature fluctuate greatly, the calibration procedure is aborted after approx. 2 min.
- · Repeat calibration if an error message appears.

Calibration Procedure

When you access calibration, the calibration values selected in the calibration presettings are used. These settings can still be changed in the Calibration menu. If you do not want to perform calibration, use the left softkey to go back one level and abort the calibration.

- 01. Select "Automatic" calibration mode and press *enter* to confirm.

 √ Display of calibration solution as configured in Cal Presettings.
- 02. Change the calibration solution, if required.
- 03. Take the sensor out of the medium and rinse it well in deionized water.
- 04. Immerse the sensor in the calibration solution.
- 05. Start calibration with the right softkey: Next.
 - √ Calibration is performed.

The following parameters are displayed: Calibration Temperature, Solution Table Value (conductivity depending on calibration temperature), and Response Time.

✓ The calibration record is displayed. By pressing the *right softkey: Adjust*, the calculated calibration values are used to calculate the measured values in Protos and also saved in the Memosens sensor.

Calibration Mode: Manual

During calibration with manual entry of the conductivity value of the calibration solution, the sensor is immersed in a calibration solution. Protos determines a conductivity/calibration temperature pair value. The temperature-corrected conductivity value of the calibration solution must then be entered. For this purpose, take the conductivity value that goes with the displayed temperature from the calibration solution TC table. Conductivity intermediate values must be interpolated. Protos automatically calculates the cell constant.

Notes on Calibration

- Use only fresh calibration solutions. The used calibration solution must be configured.
- The accuracy of the calibration is dependent on precise acquisition of the calibration solution temperature. On the basis of the measured or entered temperature, Protos calculates the setpoint of the calibration solution from a stored table.
- Note the response time of the temperature probe.
- To determine the exact cell constant, wait for temperature equalization of the temperature probe and calibration solution before calibration.
- If the measured conductance or temperature fluctuate greatly, the calibration procedure is aborted after approx. 2 min.
- · Repeat calibration if an error message appears.

Calibration Procedure

Calibration MSU ... Module Memosens Cond

When you access calibration, the calibration values selected in the calibration presettings are used. These settings can still be changed in the Calibration menu. If you do not want to perform calibration, use the left softkey to go back one level and abort the calibration.

- 01. Select "Manual" calibration mode and press *enter* to confirm.
- 02. Take the sensor out of the medium, rinse it well in deionized water, and dry it.
- 03. Immerse the sensor in the calibration solution.
- 04. Start calibration with the *right softkey: Next*.
 - √ Calibration is performed.

The following parameters are displayed: Calibration Temperature and Response Time.

- 05. Enter conductivity.
- 06. Continue with the right softkey: Next.

✓ The calibration record is displayed. By pressing the *right softkey: Adjust*, the calibration values obtained during calibration for calculation of the measured values are applied to the device. The calibration values are also stored in the sensor.

Calibration Mode: Product

If the sensor cannot be removed – e.g., for sterility reasons – its cell constant can be determined by "sampling". The currently measured value (conductivity or concentration¹⁾) in the process is stored by Protos for this purpose. Right after this, take a sample from the process. The value of this sample is separately determined under process conditions (same temperature!) wherever possible. The calculated value is entered in the measuring system. Protos calculates the conductivity sensor's cell constant from the deviation between the process measured value and the sample value.

Product Calibration without TC Compensation (With Conductivity)

A sample is taken from the process. The sample's measured value is determined in the laboratory at the temperature at which the sample was taken ("Sample Temperature", see display). It may be necessary to thermostat the sample in the laboratory accordingly. Temperature compensation of the reference measuring devices must be disabled (TC = 0 %/K).

Product Calibration with TC Compensation $T_{ref} = 25 \text{ °C/77 °F (With Conductivity)}$

A sample is taken from the process. During measurement in the laboratory (TC linear), the same values for reference temperature and temperature coefficient must be set in both the reference measuring device and Protos. In addition, the measuring temperature should match the sample temperature (see display) as closely as possible. The sample should be transported in a vacuum flask (Dewar) to ensure this.

NOTICE! Product calibration is only possible if the process medium is stable (no chemical reactions that affect conductivity). At higher temperatures, evaporation may falsify results.

Calibration Procedure

Calibration MSU ... Module Memosens Cond

When you access calibration, the calibration values selected in the calibration presettings are used. These settings can still be changed in the Calibration menu. If you do not want to perform calibration, use the left softkey to go back one level and abort the calibration.

- 01. Select "Product" calibration mode and press enter to confirm.
- 02. Prepare for sampling.
- 03. Start with the right softkey: Next.

Product calibration is performed in 2 steps.

Step 1:

- 04. Take sample.
 - ✓ The measured value and temperature at the time of sampling are displayed.
- 05. Save with the *right softkey: Save*.
 - √ An information window is shown.
- 06. Right softkey: Close
- 07. As required, exit calibration by pressing the *left softkey: Back*.

Note: The icon indicates that product calibration has not yet been completed.

Step 2: Lab value has been measured.

- 08. Open the Product Calibration menu again.
- 09. Right softkey: Next
- 10. Enter the lab value and press enter to confirm.
- 11. Confirm with the *right softkey: Next* or repeat calibration with the *left soft-key: Cancel*.
- ✓ The calibration record is displayed. By pressing the *right softkey: Adjust*, the calibration values obtained during calibration for calculation of the measured values are applied to the device. The calibration values are stored in the sensor.

Exception: Sample value can be determined and entered on site:

- 01. Take sample.
 - √ The measured value and temperature at the time of sampling are displayed.
- 02. Left softkey: Entry
- 03. Enter the lab value and press enter to confirm.
- 04. Confirm with the *right softkey: Next* or repeat calibration with the *left soft-key: Cancel*.
- ✓ The calibration record is displayed. By pressing the *right softkey: Adjust*, the calibration values obtained during calibration for calculation of the measured values are applied to the device. The calibration values are stored in the sensor.

Calibration Mode: Installation Factor

When using a 4-electrode sensor in a tight space, an installation factor can be entered for calibration / adjustment.

Calibration Procedure

Calibration MSU ... Module Memosens Cond

When you access calibration, the calibration values selected in the calibration presettings are used. These settings can still be changed in the Calibration menu. If you do not want to perform calibration, use the left softkey to go back one level and abort the calibration.

Make sure that the sensor is in normal mounting position in the medium.

- 01. Select "Install. Factor" calibration mode and press *enter* to confirm.
- 02. Enter the installation factor.
- 03. Press the right softkey: Next.

✓ The calibration record is displayed. By pressing the *right softkey: Save*, the calibration values obtained during calibration for calculation of the measured values are applied to the device. The calibration values are stored in the sensor.

Calibration Mode: Data Entry

Entry of values for the cell constant of a sensor, related to 25 °C/77 °F.

Calibration Procedure

Calibration MSU ... Module Memosens Cond

When you access calibration, the calibration values selected in the calibration presettings are used. These settings can still be changed in the Calibration menu. If you do not want to perform calibration, use the left softkey to go back one level and abort the calibration.

- 01. Select "Data Entry" calibration mode and press enter to confirm.
- 02. Remove the sensor and install the premeasured sensor.
- 03. Continue with the right softkey: Next.
- 04. Enter the cell constant of the premeasured sensor.
- ✓ The calibration record is displayed. By pressing the *right softkey: Adjust*, the calibration values obtained during calibration for calculation of the measured values are applied to the device. The calibration values are stored in the sensor.

Calibration Mode: Temperature

This function is used to adjust the individual temperature probe tolerances or cable lengths for the purpose of increasing the accuracy of the temperature measurement. The adjustment requires an accurate measurement of the process temperature using a calibrated reference thermometer. The measurement error of the reference thermometer should be less than 0.1 K. Adjustment without an accurate measurement may result in falsification of the displayed measured value.

The adjustment value is stored in the sensor.

Calibration Procedure

Calibration MSU ... Module Memosens Cond

When you access calibration, the calibration values selected in the calibration presettings are used. These settings can still be changed in the Calibration menu. If you do not want to perform calibration, use the left softkey to go back one level and abort the calibration.

- 01. Select "Temperature" calibration mode and press enter to confirm.
- 02. Enter the measured process temperature and press *enter* to confirm.
 √ The temperature offset is displayed.
- 03. Adjust the temperature detector with the *right softkey: Save*.

The current adjustment and temperature offset data can be displayed in the Diagnostics menu:


Diagnostics ▶ MSU ... Module ▶ Memosens Cond ▶ Temp. Offset Log

Cond Maintenance Functions


Note: Function check (HOLD) active

The current outputs and relay contacts behave in accordance with the parameter settings. Since the device is in function check (HOLD) mode, certain media can be used to validate the sensor and check the measured values without affecting the signal outputs.

To end the function check, return to measuring mode.



Cond Diagnostic Functions



For a detailed description of general diagnostic functions, see the basic unit user manual.

Overview of Cond Diagnostic Functions

In diagnostics mode, you can access the following submenus without interrupting the measurement:

Diagnostics ▶ MSU ... Module:

Module Diagnostics Protos periodically performs a self-test in the background.

The results can be displayed here.

Diagnostics MSU ... Module Memosens Cond:

Submenus	Description
Sensor Information	The Sensor Information submenu shows data from the currently connected Memosens sensor, e.g., manufacturer, order no., serial no., firmware and hardware version, last calibration, operating time:
Sensor Monitor	The raw measured values, such as conductance, effective resistance, and temperature are displayed for diagnostic purposes.

Cond Diagnostic Functions

Submenus	Description
Cond Cal/Adj Record	The calibration/adjustment record shows the data
	from the last calibration/adjustment performed on the currently connected sensor.
Temp. Offset Log	The temp. offset log shows the data from the last
	temperature equalization performed on the currently connected sensor.
USP Function	If configured: Display of UPS Limit value, Reduced
	Limit Value, and Conductivity
Sensor Wear Monitor	The sensor wear monitor displays the sensor operating time and the maximum temperature during the operating time.

Cond Messages

Separation Specification Speci



See also "Decommissioning"; p. 159

No.	Message Type	Cond Messages					
C008	\otimes	Factory Settings: Switch device off (approx. 10 s). If the message persists, send in the device.					
C009	\otimes	Firmware Error: Switch device off (approx. 10 s). Reload the firmware. If the message persists, send in the device.					
C010	\otimes	Conductivity Range: Measuring range exceeded. Possible causes: Sensor not/incorrectly connected, cable incorrectly connected/faulty, range incorrectly specified, incorrect cell constant set.					
C011	\otimes	Conductivity LO_LO: Value below configured monitoring limit					
C012	À	Conductivity LO: Value below configured monitoring limit					
C013	A	Conductivity HI: Value above configured monitoring limit					
C014	\otimes	Conductivity HI_HI: Value above configured monitoring limit					
C015	\otimes	Temperature Range: Measuring range exceeded. Possible causes: Sensor not/incorrectly connected, cable incorrectly connected/faulty, range incorrectly specified					
C016	\otimes	Temperature LO_LO: Value below configured monitoring limit					
C017	A	Temperature LO: Value below configured monitoring limit					
C018	A	Temperature HI: Value above configured monitoring limit					
C019	\otimes	Temperature HI_HI: Value above configured monitoring limit					
C020	\otimes	Resistivity Range: Measuring range exceeded. Possible causes: Sensor not/incorrectly connected, cable incorrectly connected/faulty, range incorrectly specified, incorrect cell constant set.					
C021	\otimes	Resistivity LO_LO: Value below configured monitoring limit					
C022	A	Resistivity LO: Value below configured monitoring limit					
C023	A	Resistivity HI: Value above configured monitoring limit					
C024	\otimes	Resistivity HI_HI: Value above configured monitoring limit					
C025	\otimes	Concentration Range: Measuring range exceeded. Possible causes: Sensor not/incorrectly connected, cable incorrectly connected/faulty, range incorrectly specified, incorrect cell constant set.					


No.	Message Type	Cond Messages			
C026	\otimes	Concentration LO_LO: Value below configured monitoring limit			
C027	Â	Concentration LO: Value below configured monitoring limit			
C028	À	Concentration HI: Value above configured monitoring limit			
C029	\otimes	Concentration HI_HI: Value above configured monitoring limit			
C040	\otimes	Salinity Range: Range exceeded. Possible causes: Sensor not/incorrectly connected, cable incorrectly connected/faulty, incorrect cell constant set.			
C041	\otimes	Salinity LO_LO: Value below configured monitoring limit			
C042	A	Salinity LO: Value below configured monitoring limit			
C043	A	Salinity HI: Value above configured monitoring limit			
C044	\otimes	Salinity HI_HI: Value above configured monitoring limit			
C045	\otimes	Conductance Range: Value above range limit. Possible causes: Sensor not/incorrectly connected, incorrect sensor for range, cable faulty (short circuit).			
C060	•	Sad Sensoface: Polarization The sensor is polarized. Sensor is not suitable for the range or the process medium: Connect a suitable sensor.			
C061	€	Sad Sensoface: Cable			
C062	User-defined	Sad Sensoface: Cell constant Incorrect cell constant set, incorrect adjustment: Repeat calibration/adjustment. Replace the sensor as required.			
C070	\otimes	TDS Range: Measuring range exceeded. Possible causes: Sensor not/incorrectly connected, cable incorrectly connected/faulty, incorrect cell constant set.			
C071	\otimes	TDS LO_LO: Value below configured monitoring limit			
C072	A	TDS LO: Value below configured monitoring limit			
C073	A	TDS HI: Value above configured monitoring limit			
C074	\otimes	TDS HI_HI: Value above configured monitoring limit			
C090	User-defined	USP Limit: The configured USP limit was exceeded.			
C091	User-defined	Reduced USP Limit: The configured reduced USP limit was exceeded.			
C110	User-defined	CIP Counter: Configured number of CIP cycles exceeded: As required, calibrate/adjust sensor or replace.			
C111	User-defined	SIP counter: Configured number of SIP cycles exceeded: As required, calibrate/adjust sensor or replace.			
C113	User-defined	Sensor Operating Time: Replace the sensor.			

No.	Message Type	Cond Messages
C120	\otimes	Wrong Sensor (Sensor Check)
C121	\otimes	Sensor Error (Factory Settings): Replace the sensor.
C122	\(\operatorname{\operatornam	Sensor Memory (Cal Data): The calibration data is defective: Recalibrate/readjust sensor.
C123	⇔	New Sensor, Adjustment Required
C124	\(\operatorname{\operatornam	Sensor Date: The sensor data is implausible. Check and, as necessary, adjust the configuration.
C130	Info	SIP Cycle Counted
C131	Info	CIP cycle counted
C200	À	Reference temperature The reference temperature for temperature compensation is invalid.
C201	Â	Temperature compensation
C202	҈Ѧ	TC Range (Maintenance Required): The measured value is at the limit of the permissible compensation range (table).
C203	\otimes	TC Range (Failure): The measured value is outside the permissible compensation range (table).
C204	Info	Cal: Sensor Unstable: The drift criterion was not adhered to during calibration. Possible causes: improper calibration, sensor cable/connection faulty, sensor worn. Check the sensor and calibration and repeat as required. Otherwise, replace the sensor.
C205	Info	Cal: Sensor Failure: Replace the sensor.
C254	Info	Module Reset

No.	Message Type	Cond / Cond Calculation Block Messages
E010	\otimes	Conductivity Diff Range: Max/min device limits exceeded: - Check both conductivity values Check sensor/cable connections.
E011	\otimes	Conductivity Diff LO_LO: Value below configured monitoring limit
E012	Â	Conductivity Diff LO: Value below configured monitoring limit
E013	À	Conductivity Diff HI: Value above configured monitoring limit
E014	\otimes	Conductivity Diff HI_HI: Value above configured monitoring limit
E015	\otimes	Temperature Diff Range: Max/min device limits exceeded: - Check both temperature values Check sensor/cable connections.
E016	\otimes	Temperature Diff LO_LO: Value below configured monitoring limit
E017	À	Temperature Diff LO: Value below configured monitoring limit
E018	À	Temperature Diff HI: Value above configured monitoring limit
E019	\otimes	Temperature Diff HI_HI: Value above configured monitoring limit
E020	\otimes	Resistivity Diff Range: Max/min device limits exceeded: - Check both resistance values Check sensor/cable connections.
E021	\otimes	Resistivity Diff LO_LO: Value below configured monitoring limit
E022	A	Resistivity Diff LO: Value below configured monitoring limit
E023	À	Resistivity Diff HI: Value above configured monitoring limit
E024	\otimes	Resistivity Diff HI_HI: Value above configured monitoring limit
E030	\otimes	RATIO Range: Min/max device limits exceeded: - Check both conductivity values.
E031	\otimes	RATIO LO_LO: Value below configured monitoring limit
E032	A	RATIO LO: Value below configured monitoring limit
E033	A	RATIO HI: Value above configured monitoring limit
E034	\otimes	RATIO HI_HI: Value above configured monitoring limit
E035	\otimes	PASSAGE Range: Min/max device limits exceeded: - Check both conductivity values.
E036	\otimes	PASSAGE LO_LO: Value below configured monitoring limit
E037	Æ	PASSAGE LO: Value below configured monitoring limit
E038	À	PASSAGE HI: Value above configured monitoring limit

No.	Message Type	Cond / Cond Calculation Block Messages			
E039	\otimes	PASSAGE HI_HI: Value above configured monitoring limit			
E045	\otimes	REJECTION Range: Min/max device limits exceeded: Check both conductivity values.			
E046	\otimes	REJECTION LO_LO: Value below configured monitoring limit			
E047	A	REJECTION LO: Value below configured monitoring limit			
E048	À	REJECTION HI: Value above configured monitoring limit			
E049	\otimes	REJECTION HI_HI: Value above configured monitoring limit			
E050	\otimes	DEVIATION Range: Min/max device limits exceeded: - Check both conductivity values.			
E051	\otimes	DEVIATION LO_LO: Value below configured monitoring limit			
E052	À	DEVIATION LO: Value below configured monitoring limit			
E053	Â	DEVIATION HI: Value above configured monitoring limit			
E054	\otimes	DEVIATION HI_HI: Value above configured monitoring limit			
E060	\otimes	pH Range: Range outside the permitted range of the VGB guideline: - Check both conductivity values Check choice of alkalizing agent Check ion exchanger Check both sensors/cables.			
E061	\otimes	pH LO_LO: Value below configured monitoring limit			
E062	À	pH LO: Value below configured monitoring limit			
E063	҈Ѧ	pH HI: Value above configured monitoring limit			
E064	\otimes	pH HI_HI: Value above configured monitoring limit			
E200	€	Calculation Block Configuration			

Note: Function check (HOLD) active




Menu Selection: Parameter Setting ▶ MSU ... Module ▶ ... Condl:

Parameter	Default	Selection / Range						
Input Filter								
Pulse Suppression	Off	Enable/disable suppression of interference pulses.						
Sensor Data (see page 130)	Sensor Data (see page 130)							
Sensoface	On	On, Off						
Sensocheck (Memosens: in "Sensor Monitoring Details")	Off	Off, Failure, Maintenance						
Sensor Monitoring Details (Memosens only) • Cell Factor • Sensocheck • Sensor Operating Time • CIP Counter • SIP Counter	Auto Off Off Off Off	Auto, Individual Off, On Off, Individual (max. 9999 d) Off, Individual Off, Individual						
Temperature Detection (SE670/SE680K only) Measuring Temp Cal Temperature Cal Presettings	Auto Auto	Auto, Manual Auto, Manual						
Calibration Mode	Automatic	Automatic, Manual, Product, Zero Point, Installation						
Calibration wode	Automatic	Factor, Data Entry, Temperature						
Cal Solution	NaCl Sat	NaCl 0.01 m: 1183 μS/cm NaCl 0.1 m: 10.683 mS/cm NaCl Sat 251.3 mS/cm KCl 0.01 m: 1413 μS/cm KCl 0.1 m: 12.88 mS/cm KCl 1m: 111.80 mS/cm						
Product Calibration	Without TC	Without TC, With TC						
TC Process Medium								
Temperature Compensation	Off	Off, Linear, EN27888, Ultrapure Water (with TAN option FW4400-008)						
Concentration (see page 135)								
Concentration	Off	Off, On						
TDS (see page 133)								
TDS Function	Off	Off, On (preset 1.00)						
USP (see page 134)	•							
USP Function	Off	Off, On						

Menu Selection: Parameter Setting ▶ MSU ... Module ▶ ... Condl:

Parameter	Default	Selection / Range
Messages		
Messages	Temperature: Max. Device Limits	Conductivity, Resistivity, Concentration, Temperature, Salinity, TDS Can be adjusted for all types of monitoring: Off, Max. Device Limits, Variable Limits)



Sensoface


Parameter Setting MSU ... Module ... Condl Sensor Data:

The Sensoface icons provide the user with diagnostic information on wear and required maintenance of the sensor. In measuring mode, an icon (happy, neutral, or sad smiley) is shown on the display to reflect the continuous monitoring

Sensoface monitors the toroidal conductivity sensor on the basis of the

Cell Factor, Zero Point, and, if Sensocheck is enabled: Primary/Secondary Coil

Additionally with Memosens sensors: Number of CIP and SIP cycles compared to the specified "Sensor Monitoring Details".



Sensocheck

Monitoring of primary and secondary coils. Additionally with Memosens sensors: Monitoring of the cell factor compared to the specified "Sensor Monitoring Details".

Enable or disable Sensocheck to generate Failure or Maint. Required messages.

Memosens: Parameter Setting ▶ MSU ... Module ▶ Memosens Condl ▶ Sensor Data Sensor Monitoring Details

SE670/SE680K: Parameter Setting MSU ... Module Condl Sensor Sensor Data

Menu

Condl Parameter Setting

CIP Counter/SIP Counter

CIP/SIP cycles are used to clean or sterilize process-wetted parts in the process. Depending on the application, either one chemical (alkaline solution, water) or several chemicals (alkaline solution, water, acidic solution, water) are used.

- CIP temperature > 55 °C / 131 °F
- SIP temperature > 115 °C / 239 °F

The cleaning (cleaning in place) and sterilization (sterilization in place) cycles are counted to measure the load on the sensor, e.g., in biotechnology applications.

Note: If measurements are generally taken at high temperatures (> 55 $^{\circ}$ C / 131 $^{\circ}$ F), the counters should be switched off.

When a CIP/SIP counter is switched on, a maximum number of cycles can be entered. A message can be set to indicate when a counter has reached a specified value.

Note: A CIP or SIP cycle is only entered into the logbook 2 hours after the start to ensure that the cycle is complete.

Action

Presettings for Calibration

Display

The calibration presettings can be 7 defined in the parameter settings or IB Cal Presettings (Administrator) adjusted directly in the Calibration Calibration Mode Manual menu prior to the calibration. Product Calibration * Product Conductivity Zero Point Parameter Setting ▶ ... Condl... ▶ Installation Fact Cal Presettings: Data Entry Back Calibration Mode: Presetting of calibration mode, e.g., Automatic, Manual, Product, Zero Point, Installation Factor, Data Entry, Temperature Further options are available depending on the calibration mode. Automatic: Selection of calibration solution Product Calibration: Conductivity/concentration¹⁾ Conductivity: Selection with/without temperature compensation Concentration: Selection of medium

Temperature Compensation of Process Medium

The following are available for temperature compensation:

- Off
- Linear (enter temperature coefficient TC)
- EN 27888 (natural waters)
- Ultrapure water (with different trace impurities)

Trace Impurities in Ultrapure Water (with TAN option FW4400-008)

NaCl	Neutral ultrapure water, for conductivity measurement in water
	processing downstream of gravel bed filter
HCl	Acidic ultrapure water, for conductivity measurement downstream
	of cation filter

NH₃ Ammoniacal ultrapure water

NaOH Alkaline ultrapure water

Note: If TC correction for the process medium is enabled, "TC" is shown on the display in measuring mode.

TDS Function

TDS (total dissolved solids) = weight of dissolved solids that influence conductivity

The TDS function provides a quick way of determining the evaporation residue of water. A TDS factor must be entered for this purpose.

The factor establishes a simple linear relationship between measured conductivity and evaporation residue. It is dependent on the composition of the medium and must be empirically determined by the user.

USP Function (Condl)

Monitoring Ultrapure Water in the Pharmaceutical Industry

The conductivity of ultrapure water in the pharmaceutical industry can be monitored online in accordance with the guideline "USP" (U.S. Pharmacopeia), Annex 5, Section 645 "Water Conductivity". The conductivity is measured without temperature compensation and is compared with limits. The water is usable without any further testing when the conductivity is below the USP limit.

Configuring the USP Function

The USP value can be configured as a process variable USP% for output (display, current output, limit, measurement recorder)
Settings can be changed in the USP submenu:
Parameter Setting MSU ... Module ... Condl USP

Reduced Limit Value: The USP limit can be reduced down to 10 %. **Monitoring:** Select whether and how an exceeded limit is to be displayed.

Off No message, but the parameter is still shown in the

Diagnostics menu.

Failure A failure message is shown in off-limit conditions;

the corresponding NAMUR icon is displayed.

Maintenance A Maint. Required message is shown in off-limit conditions;

the corresponding NAMUR icon is displayed.

USP Function: Specify a Relay Contact

The USP function can also be assigned to a relay contact:

Parameter Setting > BASE Module ... > Contact K... > Usage: USP Output

Display of the USP Function in the Diagnostics Menu

Diagnostics MSU ... Module ... Condl USP Function Display of the USP limit, the reduced limit, and conductivity.

Concentration (Condl)

Note: Concentration determination requires the activation of TAN option FW4400-009.

The substance concentration in percent by weight (wt%) is determined for H₂SO₄, HNO₃, HCl, NaOH, NaCl, and Oleum (see page 174ff).

Conditions for Concentration Determination

The following conditions must be met for a reliable concentration determination:

- For calculation of concentration, the medium to be measured must be a purely binary mixture (e.g., water-hydrochloric acid). Presence of other dissolved substances (e.g., salts) leads to incorrect concentration values.
- In the region of small slopes (e.g., at the range limits) small changes in conductivity can correspond to large changes in concentration. This may lead to an unsteady display of the concentration value.
- As the concentration value is calculated from the measured conductivity and temperature values, accurate temperature measurement is very important. Therefore, you should make sure that conductivity sensor and process medium are in thermal equilibrium.

Settings can be changed in the Concentration submenu:

Parameter Setting MSU ... Module ... Condl Concentration

- 01. Concentration: On
- 02. Select the medium:

```
NaCl (0-28 %), HCl (0-18 %), NaOH (0-24 %), \rm H_2SO_4 (0-37 %), HNO_3 (0-30 %), \rm H_2SO_4 (89-99 %), HCl (22-39 %), HNO_3 (35-96 %), \rm H_2SO_4 (28-88 %), NaOH (15-50 %), Oleum (12-45 %), table
```

You can define limits for warning and failure messages for the concentration value:

Parameter Setting MSU ... Module ... Condl Messages Concentration Messages

Concentration (Condl)

Specifying a Concentration Solution for Conductivity Measurement

To specify the customer-specific solution, 5 concentration values A-E are entered in a matrix together with 5 temperature values 1-5. First enter the 5 temperature values, then the corresponding conductivity values for each of the concentrations A-E.

These solutions are then available as "Table" in addition to the permanently set standard solutions.

Settings can be changed in the System Control in the Concentration Table submenu:

Parameter Setting > System Control > Concentration Table:

- 01. Enter temperatures 1 to 5.
- 02. Enter values for concentrations A-E for the respective temperatures.

Notes: The temperature values must be rising (Temp. 1 is the lowest, Temp. 5 the highest temperature).

The concentration values must be rising (Conc. A is the lowest, conc. E the highest concentration).

The table values A1 ... E1, A2 ... E2, etc. must all be rising within the table or all falling.

Points of inflection are not allowed.

Incorrect table entries are indicated by an exclamation point in a red triangle. The table is built up as 5x5 matrix:

	Conc. A	Conc. B	Conc. C	Conc. D	Conc. E
Temp 1	A1	B1	C1	D1	E1
Temp 2	A2	B2	C2	D2	E2
Temp 3	A3	В3	C3	D3	E3
Temp 4	A4	B4	C4	D4	E4
Temp 5	A5	B5	C5	D5	E5


The concentration table is selected in the menu:

Parameter Setting MSU ... Module ... Condl Cal Presettings

Calibration Mode: Automatic

Cal Solution: Table

Back



Note: During calibration, the function check (HOLD) operating state is active for the corresponding module channel. The assigned current outputs and relay contacts behave as configured (BASE module).

The calibration data is saved in the Memosens sensor. This means that Memosens sensors can be cleaned, reconditioned, calibrated, and adjusted away from the measurement location, e.g., in a laboratory. Sensors in the system are replaced on-site by adjusted sensors.


Calibration: Detecting deviations without readjustment **Adjustment:** Detecting deviations with readjustment

Adjustment

Adjustment means that the values determined by a calibration are applied to the sensor. The zero point and slope values determined during calibration are entered in the adjustment record:

Diagnostics MSU ... Module ... Condl Cal/Adj Record

These values are not used to calculate the process variables until the calibration has been terminated with an adjustment.



Explanations Regarding Calibration/Adjustment with Toroidal Sensors

Each inductive (toroidal) conductivity sensor has an individual cell factor. The cell factor may vary depending on the sensor design. Because the conductivity value is calculated from the measured conductance and the cell factor, the measuring system must know the cell factor. During calibration or sensor adjustment, either the known (printed) cell factor of the used toroidal conductivity sensor is entered in the measuring system, or it is determined automatically by measuring a calibration solution with known conductivity.

Notes on Calibration

- · Use only fresh calibration solutions.
- The used calibration solution must be configured.
- The accuracy of the calibration is crucially dependent on a precise acquisition of the calibration solution temperature. On the basis of the measured or entered temperature, Protos calculates the setpoint of the calibration solution from a stored table.
- Note the response time of the temperature probe.
- To determine the exact cell factor, wait for temperature equalization of the temperature probe and calibration solution before calibration.
- If the measured conductance or temperature fluctuate greatly, the calibration procedure is aborted after approx. 2 min. Repeat calibration if an error message appears.

Since the cell factor is subject to production-related fluctuations, it is recommended that the removed sensor be calibrated with a calibration solution (e.g., saturated NaCl).

 If the installation space is tight (minimum distances are not reached), the sensor must be adjusted in its installed state, as the resulting cell factor has changed.

Calibration Mode: "Product Calibration".

Temperature Compensation during Calibration

The conductivity value of the calibration solution is temperature-dependent. During calibration, the temperature of the calibration solution must therefore be known so that the actual value can be obtained from the conductivity table.

Automatic Temperature Compensation


During automatic logging of the calibration temperature, Protos measures the temperature of the calibration solution using the temperature probe integrated in the Memosens sensor.


If the sensor does not have an integrated temperature detector:

- Connect an external temperature probe and select it in the Parameter Setting menu.
- Set the manual temperature for calibration.

Function Check (HOLD) During Calibration / Adjustment

Behavior of the signal and switching outputs during calibration / adjustment





Calibration Mode: Automatic

During automatic calibration:, the conductivity sensor is immersed in a standard calibration solution (NaCl or KCl, set during parameter setting). On the basis of the measured conductance and temperature, Protos automatically calculates the cell factor.

The temperature dependence of the calibration solution is accounted for.

Notes on Calibration

- Use only fresh calibration solutions. The used calibration solution must be configured.
- The accuracy of the calibration is crucially dependent on a precise acquisition of the calibration solution temperature. On the basis of the measured or entered temperature, Protos calculates the setpoint of the calibration solution from a stored table.
- Note the response time of the temperature probe.
- To determine the exact cell factor, wait for temperature equalization of the temperature probe and calibration solution before calibration.
- If the measured conductance or temperature fluctuate greatly, the calibration procedure is aborted after approx. 2 min.
- Repeat calibration if an error message appears.

Calibration Procedure

Calibration MSU ... Module ... Condl

When you access calibration, the calibration values selected in the calibration presettings are used. These settings can still be changed in the Calibration menu. If you do not want to perform calibration, use the left softkey to go back one level and abort the calibration.

- 01. Select "Automatic" calibration mode and press *enter* to confirm.
 - √ Display of calibration solution as configured in Cal Presettings.
- 02. Change the calibration solution, if required.
- 03. Take the sensor out of the medium, rinse it well in deionized water, and dry it.
- 04. Immerse the sensor in the calibration solution.
- 05. Start calibration with the right softkey: Next.
 - √ Calibration is performed.

The following parameters are displayed: Calibration Temperature, Solution Table Value (conductivity depending on calibration temperature), and Response Time.

✓ By pressing the *right softkey: Adjust*, the calibration values obtained during calibration for calculation of the process variables are applied to the device.

Calibration Mode: Manual

During calibration with manual entry of the conductivity value of the calibration solution, the sensor is immersed in a calibration solution. Protos determines a conductivity / calibration temperature pair value. The temperature-corrected conductivity value of the calibration solution must then be entered. For this purpose, take the conductivity value that goes with the displayed temperature from the calibration solution TC table. Conductivity intermediate values must be interpolated.

Protos automatically calculates the cell factor.

Notes on Calibration

- Use only fresh calibration solutions. The used calibration solution must be configured.
- The accuracy of the calibration is crucially dependent on a precise acquisition of the calibration solution temperature. On the basis of the measured or entered temperature, Protos calculates the setpoint of the calibration solution from a stored table.
- Note the response time of the temperature probe.
- To determine the exact cell factor, wait for temperature equalization of the temperature probe and calibration solution before calibration.
- If the measured conductance or temperature fluctuate greatly, the calibration procedure is aborted after approx. 2 min.
- Repeat calibration if an error message appears.

Calibration Procedure

Calibration MSU ... Module ... Condl

When you access calibration, the calibration values selected in the calibration presettings are used. These settings can still be changed in the Calibration menu. If you do not want to perform calibration, use the left softkey to go back one level and abort the calibration.

- 01. Select "Manual" calibration mode and press *enter* to confirm.
- 02. Take the sensor out of the medium and rinse it well in deionized water.
- 03. Immerse the sensor in the calibration solution.
- 04. Start calibration with the right softkey: Next.
 - √ Calibration is performed.
 - The following parameters are displayed: Calibration Temperature and Response Time.
- 05. Enter conductivity.
- 06. Continue with the right softkey: Next.
- ✓ By pressing the *right softkey: Adjust*, the calibration values obtained during calibration for calculation of the process variables are applied to the device.

Calibration Mode: Product

If the sensor cannot be removed – e.g., for sterility reasons – its cell factor can be determined by "sampling". The currently measured value (conductivity or concentration¹⁾) in the process is stored by Protos for this purpose. Right after this, take a sample from the process. The value of this sample is separately determined under process conditions (same temperature!) wherever possible. The calculated value is entered in the measuring system. Protos calculates the conductivity sensor's cell factor from the deviation between the process measured value and the sample value.

Product Calibration without TC Compensation (With Conductivity)

A sample is taken from the process. The sample's measured value is determined in the laboratory at the temperature at which the sample was taken ("Sample Temperature", see display). It may be necessary to thermostat the sample in the laboratory accordingly. Temperature compensation of the reference measuring devices must be disabled (TC = 0 %/K).

Product Calibration with TC Compensation $T_{ref} = 25 \, ^{\circ}\text{C}/77 \, ^{\circ}\text{F}$ (With Conductivity)

A sample is taken from the process. During measurement in the laboratory (TC linear), the same values for reference temperature and temperature coefficient must be set in both the reference measuring device and Protos. In addition, the measuring temperature should match the sample temperature (see display) as closely as possible. The sample should be transported in a vacuum flask (Dewar) to ensure this.

NOTICE! Product calibration is only possible if the process medium is stable (no chemical reactions that affect conductivity). At higher temperatures, evaporation may falsify results.

Calibration Procedure

Calibration MSU ... Module Memosens Condl

When you access calibration, the calibration values selected in the calibration presettings are used. These settings can still be changed in the Calibration menu. If you do not want to perform calibration, use the left softkey to go back one level and abort the calibration.

- 01. Select "Product" calibration mode and press enter to confirm.
- 02. Prepare for sampling.
- 03. Start with the right softkey: Next.

Product calibration is performed in 2 steps.

Step 1:

- 04. Take sample.
 - ✓ The measured value and temperature at the time of sampling are displayed.
- 05. Save with the right softkey: Save.
 - √ An information window is shown.
- 06. Right softkey: Close
- 07. As required, exit calibration by pressing the left softkey: Back.

Note: The icon indicates that product calibration has not yet been completed.

Step 2: Lab value has been measured.

- 08. Open the Product Calibration menu again.
- 09. Right softkey: Next
- 10. Enter the lab value and press enter to confirm.
- 11. Confirm with the *right softkey: Next* or repeat calibration with the *left soft-key: Cancel*.
- ✓ The calibration record is displayed. By pressing the *right softkey: Adjust*, the calibration values obtained during calibration for calculation of the measured values are applied to the device. The calibration values are stored in the sensor.

Exception: Sample value can be determined and entered on site:

- 01. Take sample.
 - √ The measured value and temperature at the time of sampling are displayed.
- 02. Left softkey: Entry
- 03. Enter the lab value and press enter to confirm.
- 04. Confirm with the *right softkey: Next* or repeat calibration with the *left soft-key: Cancel*.
- ✓ The calibration record is displayed. By pressing the *right softkey: Adjust*, the calibration values obtained during calibration for calculation of the measured values are applied to the device. The calibration values are stored in the sensor.

Calibration Mode: Zero Point

Zero point correction

Calibration Procedure

Calibration MSU ... Module ... Condl

When you access calibration, the calibration values selected in the calibration presettings are used. These settings can still be changed in the Calibration menu. If you do not want to perform calibration, use the left softkey to go back one level and abort the calibration.

- 01. Select "Zero Point" calibration mode.
- 02. Take the sensor out of the medium, rinse it in deionized water, and dry it. The sensor should be dry, since zero calibration is performed in air.
- O3. Press the *right softkey: Next*.
 ✓ Zero point correction is carried out. The permissible zero deviation is model-dependent; on the SE 670 sensor, for example, it is ± 0.050 mS/cm.
- 04. Press the right softkey: Next.
- ✓ The calibration record is displayed. By pressing the *right softkey: Adjust*, the calibration values obtained during calibration for calculation of the measured values are applied to the device. The calibration values are stored in the sensor.

Calibration Mode: Installation Factor

When using a Memosens sensor in a tight space, an installation factor can be entered for calibration / adjustment.

Calibration Procedure

Calibration MSU ... Module Memosens Condl

When you access calibration, the calibration values selected in the calibration presettings are used. These settings can still be changed in the Calibration menu. If you do not want to perform calibration, use the left softkey to go back one level and abort the calibration.

Make sure that the sensor is in normal mounting position in the medium.

- 01. Select "Install. Factor" calibration mode and press enter to confirm.
- 02. Enter the installation factor.
- 03. Press the right softkey: Next.

✓ The calibration record is displayed. By pressing the *right softkey: Save*, the calibration values obtained during calibration for calculation of the measured values are applied to the device. The calibration values are stored in the sensor.

Calibration Mode: Data Entry

Entry of values for the cell factor and zero point of a sensor, related to $25 \, ^{\circ}\text{C}/77 \, ^{\circ}\text{F}$.

If concentration measurement is activated (TAN option FW-E009), the concentration is also shown in this menu and directly adjusted with the cell factor. This makes direct calibration to the concentration value possible.

Calibration Procedure

Calibration MSU ... Module ... Condl

When you access calibration, the calibration values selected in the calibration presettings are used. These settings can still be changed in the Calibration menu. If you do not want to perform calibration, use the left softkey to go back one level and abort the calibration.

- 01. Select "Data Entry" calibration mode and press enter to confirm.
- 02. Remove the sensor and install the premeasured sensor.
- 03. Continue with the *right softkey: Next*.
- 04. Enter the cell factor of the premeasured sensor.
- ✓ By pressing the *right softkey: Adjust*, the calibration values obtained during calibration for calculation of the process variables are applied to the device. The calibration values are stored in the sensor.

Calibration Mode: Temperature

This function is used to adjust the individual temperature detector tolerances or cable lengths for the purpose of increasing the accuracy of the temperature measurement. The adjustment requires an accurate measurement of the process temperature using a calibrated reference thermometer. The measurement error of the reference thermometer should be less than 0.1 K. Adjustment without an accurate measurement may result in falsification of the displayed measured value.

Calibration Procedure

Calibration MSU ... Module ... Condl

When you access calibration, the calibration values selected in the calibration presettings are used. These settings can still be changed in the Calibration menu. If you do not want to perform calibration, use the left softkey to go back one level and abort the calibration.

- 01. Select "Temperature" calibration mode and press enter to confirm.
- 02. Enter the measured process temperature and press *enter* to confirm.
 √ The temperature offset is displayed.
- 03. Adjust the temperature detector with the right softkey: Save.

When using the SE670 or SE680K sensors, the current adjustment and temperature offset data can be displayed in the Diagnostics menu:


Diagnostics MSU ... Module Sensor Condl Temp. Offset Log

Condl Maintenance Functions


Note: Function check (HOLD) active

The current outputs and relay contacts behave in accordance with the parameter settings. Since the device is in function check (HOLD) mode, certain media can be used to validate the sensor and check the measured values without affecting the signal outputs.

To end the function check, return to measuring mode.



Condl Diagnostic Functions



For a detailed description of general diagnostic functions, see the basic unit user manual.

Overview of Condl Diagnostic Functions

In diagnostics mode, you can access the following submenus without interrupting the measurement:

Diagnostics MSU ... Module:

Module Diagnostics Protos periodically performs a self-test in the background. The results can be displayed here.

Diagnostics ▶ MSU ... Module ▶ ...Condl:

Submenus	Description
Sensor Information	The Sensor Information submenu shows data from the currently connected Memosens sensor, e.g., manufacturer, order no., serial no., firmware and hardware version, last calibration, operating time.
Sensor Monitor	The raw measured values are displayed in the sensor monitor for diagnostic purposes.
Condl Cal/Adj Record	The calibration/adjustment record shows the data from the last calibration/adjustment performed on the currently connected sensor.
Temp. Offset Log (SE670/SE680K only)	The temp. offset log shows the data from the last temperature equalization performed on the currently connected sensor.

Condl Messages

Condl Messages



See also "Decommissioning"; p. 159

No.	Message Type	Condl Messages
T008	\otimes	Factory Settings: Switch device off (approx. 10 s). If the message persists, send in the device.
T009	\otimes	Firmware Error: Switch device off (approx. 10 s). Reload the firmware. If the message persists, send in the device.
T010	\otimes	Conductivity Range: Measuring range exceeded. Possible causes: Sensor not/incorrectly connected, cable incorrectly connected/faulty, range incorrectly specified, incorrect cell factor set.
T011	\otimes	Conductivity LO_LO: Value below configured monitoring limit
T012	À	Conductivity LO: Value below configured monitoring limit
T013	À	Conductivity HI: Value above configured monitoring limit
T014	\otimes	Conductivity HI_HI: Value above configured monitoring limit
T015	\otimes	Temperature Range: Measuring range exceeded. Possible causes: Sensor not/incorrectly connected, cable incorrectly connected/faulty, range incorrectly specified
T016	\otimes	Temperature LO_LO: Value below configured monitoring limit
T017	À	Temperature LO: Value below configured monitoring limit
T018	A	Temperature HI: Value above configured monitoring limit
T019	\otimes	Temperature HI_HI: Value above configured monitoring limit
T020	\otimes	Resistivity Range: Measuring range exceeded. Possible causes: Sensor not/incorrectly connected, cable incorrectly connected/faulty, range incorrectly specified, incorrect cell factor set.
T021	\otimes	Resistivity LO_LO: Value below configured monitoring limit
T022	A	Resistivity LO: Value below configured monitoring limit
T023	A	Resistivity HI: Value above configured monitoring limit
T024	\otimes	Resistivity HI_HI: Value above configured monitoring limit
T025	\otimes	Concentration Range: Measuring range exceeded. Possible causes: Sensor not/incorrectly connected, cable incorrectly connected/faulty, range incorrectly specified, incorrect cell factor set.

Condl Messages

No.	Message Type	Condl Messages
T026	\otimes	Concentration LO_LO: Value below configured monitoring limit
T027	À	Concentration LO: Value below configured monitoring limit
T028	҈Ѧ	Concentration HI: Value above configured monitoring limit
T029	\otimes	Concentration HI_HI: Value above configured monitoring limit
T040	\otimes	Salinity Range: Range exceeded. Possible causes: Sensor not/incorrectly connected, cable incorrectly connected/faulty, incorrect cell factor set.
T041	\otimes	Salinity LO_LO: Value below configured monitoring limit
T042	À	Salinity LO: Value below configured monitoring limit
T043	A	Salinity HI: Value above configured monitoring limit
T044	\otimes	Salinity HI_HI: Value above configured monitoring limit
T045	\otimes	Conductance Range: Value above range limit. Possible causes: Sensor not/incorrectly connected, incorrect sensor for range, cable faulty (short circuit).
T060	\(\operatorname{\operatornam	Sad Sensoface: Primary Coil Sensor faulty: Replace the sensor.
T061	⊕	Sad Sensoface: Secondary Coil Sensor faulty: Replace the sensor.
T063	⊗	Sad Sensoface: Zero Point Adjust the sensor zero point.
T064	User-defined	Sad Sensoface: Cell Factor Incorrect cell factor set, incorrect adjustment: Repeat calibration/adjustment. Replace the sensor as required.
T070	\otimes	TDS Range: Measuring range exceeded. Possible causes: Sensor not/incorrectly connected, cable incorrectly connected/faulty, range incorrectly specified, incorrect cell factor set.
T071	\otimes	TDS LO_LO: Value below configured monitoring limit
T072	A	TDS LO: Value below configured monitoring limit
T073	A	TDS HI: Value above configured monitoring limit
T074	\otimes	TDS HI_HI: Value above configured monitoring limit
T090	User-defined	USP Limit: The configured USP limit was exceeded.
T091	User-defined	Reduced USP Limit: The configured reduced USP limit was exceeded.
T110	User-defined	CIP Counter: Configured number of CIP cycles exceeded: As required, calibrate/adjust sensor or replace.
T111	User-defined	SIP Counter: Configured number of SIP cycles exceeded: As required, calibrate/adjust sensor or replace.

Condl Messages

No.	Message Type	Condl Messages
T113	User-defined	Sensor Operating Time: Replace the sensor.
T120	\otimes	Wrong Sensor (Sensor Check)
T121	\otimes	Sensor Error (Factory Settings): Replace the sensor.
T122	©	Sensor Memory (Cal Data): The calibration data is defective: Recalibrate/readjust sensor.
T123		New Sensor, Adjustment Required
T124	⇔	Sensor Date: The sensor data is implausible. Check and, as necessary, adjust the configuration.
T130	Info	SIP Cycle Counted
T131	Info	CIP Cycle Counted
T200	À	Reference temperature The reference temperature for temperature compensation is invalid.
T201	A	Temperature compensation
T202	À	TC Range (Maintenance Required): The measured value is at the limit of the permissible compensation range (table).
T203	\otimes	TC Range (Failure): The measured value is outside the permissible compensation range (table).
T204		Sensor Coding
T205	Info	Cal: Sensor Unstable: The drift criterion was not adhered to during calibration. Possible causes: improper calibration, sensor cable/connection faulty, sensor worn. Check the sensor and calibration and repeat as required. Otherwise, replace the sensor.
T254	Info	Module Reset

Decommissioning

Returns

If required, send the product in a clean condition and securely packed to your local contact, see www.knick.de.

Disposal

The local codes and regulations must be observed when disposing of the product.

Specifications

Sensor input	Interface for Memosens I, II, III (channels A, B, C)				
·	Channel B: Add-on function FW4400-014				
	Channels B+C: Add-on function FW4400-018				
Power supply	$U = 2.99 \dots 3.22 \text{ V}, I_{\text{max}} = 6 \text{ mA}$				
Explosion protection (MSU 4400X-180)	For entity parameters, see attachment to certificates or control drawings.				
Interface	RS-485				
Transfer rate	9,600 Bd				
Max. cable length	100 m				
Current input	0/4 20 mA / 100 Ω				
	e.g., for external pressure signal with OXY				
Start/end of scale	Can be configured within range				
Characteristic curve	Linear				
Measurement error	< 1% current value + 0.1 mA				
	(± 1 count, plus sensor error)				
General data					
RoHS conformity	According to EU directive 2011/65/EU				
EMC	EN 61326-1, EN 61326-2-3, NAMUR NE 21				
Emitted interference	Industrial applications ¹⁾ (EN 55011 Group 1 Class A)				
Immunity to interference	Industrial applications				
Lightning protection	to EN 61000-4-5, Installation class 2				
Rated operating conditions (module installed)					
Ambient temperature	Safe area: –20 55 °C/-4 131 °F Ex: –20 50 °C/-4 122 °F				
Relative humidity	5 95 %				
Climatic class	3K5 according to EN 60721-3-3				
Location class	C1 according to EN 60654-1				

¹⁾ This equipment is not designed for domestic use, and is unable to guarantee adequate protection of the radio reception in such environments.

Specifications

Tuesday Adams and Adams an	20 70 % / 4 150 %		
Transport / storage temperature	-20 /U °C / -4 158 °F		
Screw clamp connectors	Tightening torque 0.5 0.6 Nm		
	Single or stranded wires 0.2 2.5 mm ²		
Wiring	Stripping length max. 7 mm		
	Temperature resistance > 75 °C / 167 °F		
Power supply (KBUS)	6.8 8.0 V / ≤ 75 mA		

Mettler-Toledo Buffer Table

°C	рН				
0	2,03	4,01	7,12	9,52	
5	2,02	4,01	7,09	9,45	
10	2,01	4,00	7,06	9,38	
15	2,00	4,00	7,04	9,32	
20	2,00	4,00	7,02	9,26	
25	2,00	4,01	7,00	9,21	
30	1,99	4,01	6,99	9,16	
35	1,99	4,02	6,98	9,11	
40	1,98	4,03	6,97	9,06	
45	1,98	4,04	6,97	9,03	
50	1,98	4,06	6,97	8,99	
55	1,98	4,08	6,98	8,96	
60	1,98	4,10	6,98	8,93	
65	1,99	4,13	6,99	8,90	
70	1,99	4,16	7,00	8,88	
75	2,00	4,19	7,02	8,85	
80	2,00	4,22	7,04	8,83	
85	2,00	4,26	7,06	8,81	
90	2,00	4,30	7,09	8,79	
95	2,00	4,35	7,12	8,77	

Buffer Table, Knick CaliMat

°C	рН				
Order No.	CS-P0200A/	CS-P0400A/	CS-P0700A/	CS-P0900A/	CS-P1200A/
0	2.01	4.05	7.09	9.24	12.58
5	2.01	4.04	7.07	9.16	12.39
10	2.01	4.02	7.04	9.11	12.26
15	2.00	4.01	7.02	9.05	12.13
20	2.00	4.00	7.00	9.00	12.00
25	2.00	4.01	6.99	8.95	11.87
30	2.00	4.01	6.98	8.91	11.75
35	2.00	4.01	6.96	8.88	11.64
40	2.00	4.01	6.96	8.85	11.53
50	2.00	4.01	6.96	8.79	11.31
60	2.00	4.00	6.96	8.73	11.09
70	2.00	4.00	6.96	8.70	10.88
80	2.00	4.00	6.98	8.66	10.68
90	2.00	4.00	7.00	8.64	10.48

Buffer Table, DIN 19267

°C	рН				
0	1,08	4,67	6,89	9,48	13,95*
5	1,08	4,67	6,87	9,43	13,63*
10	1,09	4,66	6,84	9,37	13,37
15	1,09	4,66	6,82	9,32	13,16
20	1,09	4,65	6,80	3,27	12,96
25	1,09	4,65	6,79	9,23	12,75
30	1,10	4,65	6,78	9,18	12,61
35	1,10	4,65	6,77	9,13	12,45
40	1,10	4,66	6,76	9,09	12,29
45	1,10	4,67	6,76	9,04	12,09
50	1,11	4,68	6,76	9,00	11,98
55	1,11	4,69	6,76	8,96	11,79
60	1,11	4,70	6,76	8,92	11,69
65	1,11	4,71	6,76	8,90	11,56
70	1,11	4,72	6,76	8,88	11,43
75	1,11	4,73	6,77	8,86	11,31
80	1,12	4,75	6,78	8,85	11,19
85	1,12	4,77	6,79	8,83	11,09
90	1,13	4,79	6,80	8,82	10,99
95	1,13*	4,82*	6,81*	8,81*	10,89*

^{*} extrapoliert / extrapolated / extrapolée

NIST Standard Buffer Table (DIN 19266: 2000-01)

°C	рН			
0				
5	1.668	4.004	6.950	9.392
10	1.670	4.001	6.922	9.331
15	1.672	4.001	6.900	9.277
20	1.676	4.003	6.880	9.228
25	1.680	4.008	6.865	9.184
30	1,685	4.015	6.853	9.144
37	1,694	4.028	6.841	9.095
40	1.697	4.036	6.837	9.076
45	1.704	4.049	6.834	9.046
50	1.712	4.064	6.833	9.018
55	1.715	4.075	6.834	9.985
60	1.723	4.091	6.836	8.962
70	1.743	4.126	6.845	8.921
80	1.766	4.164	6.859	8.885
90	1.792	4.205	6.877	8.850
95	1.806	4.227	6.886	8.833

Note:

The actual pH(S) values of the individual batches of the reference materials are documented in a certificate of an accredited laboratory. This certificate is supplied with the respective buffers. Only these pH(S) values shall be used as standard values for the secondary reference buffer materials. Correspondingly, this standard does not include a table with standard pH values for practical use. The table above only provides examples of pH(PS) values for orientation.

Buffer Table, NIST Technical Buffers

°C	рН			
0	4.00	7.14	10.30	
5	4.00	7.10	10.23	
10	4.00	7.04	10.11	
15	4.00	7.04	10.11	
20	4.00	7.02	10.05	
25	4.01	7.00	10.00	
30	4.01	6.99	9.96	
35	4.02	6.98	9.92	
40	4.03	6.98	9.88	
45	4.05	6.98	9.85	
50	4.06	6.98	9.82	
55	4.07	6.98	9.79	
60	4.09	6.99	9.76	
65	4.09 *	6.99 *	9.76 *	
70	4.09 *	6.99 *	9.76 *	
75	4.09 *	6.99 *	9.76 *	
80	4.09 *	6.99 *	9.76 *	
85	4.09 *	6.99 *	9.76 *	
90	4.09 *	6.99 *	9.76 *	
95	4.09 *	6.99 *	9.76 *	

^{*} Values complemented

Buffer Table, Hamilton

°C	рН				
0	1,99	4,01	7,12	10,19	12,46
5	1,99	4,01	7,09	10,19	12,46
10	2,00	4,00	7,06	10,15	12,34
15	2,00	4,00	7,04	10,11	12,23
20	2,00	4,00	7,02	10,06	12,11
25	2,00	4,01	7,00	10,01	12,00
30	1,99	4,01	6,99	9,97	11,90
35	1,98	4,02	6,98	9,92	11,80
40	1,98	4,03	6,97	9,86	11,70
45	1,97	4,04	6,97	9,83	11,60
50	1,97	4,06	6,97	9,79	11,51
55	1,97	4,08	6,98	9,77	11,51
60	1,97	4,10	6,98	9,75	11,51
65	1,97	4,13	6,99	9,74	11,51
70	1,97	4,16	7,00	9,73	11,51
75	1,97	4,19	7,02	9,73	11,51
80	1,97	4,22	7,04	9,73	11,51
85	1,97	4,26	7,06	9,74	11,51
90	1,97	4,30	7,09	9,75	11,51
95	1,97	4,35	7,09	9,75	11,51

Buffer Table, Kraft

℃	рН				
0 5	2.01 2.01	4.05 4.04	7.13 7.07	9.24 9.16	11.47* 11.47
10	2.01	4.02	7.05	9.11	11.31
15	2.00	4.01	7.02	9.05	11.15
20	2.00	4.00	7.00	9.00	11.00
25	2.00	4.01	6.98	8.95	10.85
30	2.00	4.01	6.98	8.91	10.71
35	2.00	4.01	6.96	8.88	10.57
40	2.00	4.01	6.95	8.85	10.44
45	2.00	4.01	6.95	8.82	10.31
50	2.00	4.00	6.95	8.79	10.18
55	2.00	4.00	6.95	8.76	10.18*
60	2.00	4.00	6.96	8.73	10.18*
65	2.00	4.00	6.96	8.72	10.18*
70	2.01	4.00	6.96	8.70	10.18*
75	2.01	4.00	6.96	8.68	10.18*
80	2.01	4.00	6.97	8.66	10.18*
85	2.01	4.00	6.98	8.65	10.18*
90	2.01	4.00	7.00	8.64	10.18*
95	2.01	4.00	7.02	8.64	10.18*

^{*} Values complemented

Buffer Table, Hamilton A

℃	рН				
0	1.99	4.01	7.12	9.31	11.42
5	1.99	4.01	7.09	9.24	11.33
10	2.00	4.00	7.06	9.17	11.25
15	2.00	4.00	7.04	9.11	11.16
20	2.00	4.00	7.02	9.05	11.07
25	2.00	4.01	7.00	9.00	11.00
30	1.99	4.01	6.99	8.95	10.93
35	1.98	4.02	6.98	8.90	10.86
40	1.98	4.03	6.97	8.85	10.80
45	1.97	4.04	6.97	8.82	10.73
50	1.97	4.05	6.97	8.78	10.67
55	1.98	4.06	6.98	8.75	10.61
60	1.98	4.08	6.98	8.72	10.55
65	1.98	4.10	6.99	8.70	10.49
70	1.99	4.12	7.00	8.67	10.43
75	1.99	4.14	7.02	8.64	10.38
80	2.00	4.16	7.04	8.62	10.33
85	2.00	4.18	7.06	8.60	10.28
90	2.00	4.21	7.09	8.58	10.23
95	2.00	4.24	7.12	8.56	10.18

Buffer Table, Hamilton B

_°C	рН				
0	1.99	4.01	6.03	9.31	11.42
5	1.99	4.01	6.02	9.24	11.33
10	2.00	4.00	6.01	9.17	11.25
15	2.00	4.00	6.00	9.11	11.16
20	2.00	4.00	6.00	9.05	11.07
25	2.00	4.01	6.00	9.00	11.00
30	1.99	4.01	6.00	8.95	10.93
35	1.98	4.02	6.00	8.90	10.86
40	1.98	4.03	6.01	8.85	10.80
45	1.97	4.04	6.02	8.82	10.73
50	1.97	4.05	6.04	8.78	10.67
55	1.98	4.06	6.06	8.75	10.61
60	1.98	4.08	6.09	8.72	10.55
65	1.98	4.10	6.11	8.70	10.49
70	1.99	4.12	6.13	8.67	10.43
75	1.99	4.14	6.15	8.64	10.38
80	2.00	4.16	6.18	8.62	10.33
85	2.00	4.18	6.21	8.60	10.28
90	2.00	4.21	6.24	8.58	10.23
95	2.00	4.24	6.27	8.56	10.18

Buffer Table, HACH

°C	рН		
0	4,00	7,118	10,30
5	4,00	7,087	10,23
10	4,00	7,059	10,17
15	4,00	7,036	10,11
20	4,00	7,016	10,05
25	4,01	7,000	10,00
30	4,01	6,987	9,96
35	4,02	6,977	9,92
40	4,03	6,970	9,88
45	4,05	6,965	9,85
50	4,06	6,964	9,82
55	4,07	6,965	9,79
60	4,09	6,968	9,76
65	4,10	6,980	9,71
70	4,12	7,000	9,66
75	4,14	7,020	9,63
80	4,16	7,040	9,59
85	4,18	7,060	9,56
90	4,21	7,090	9,52
95	4,24	7,120	9,48

Buffer Table, Ciba

°C	pН			
0	2,04	4,00	7,10	10,30
5	2,09	4,02	7,08	10,21
10	2,07	4,00	7,05	10,14
15	2,08	4,00	7,02	10,06
20	2,09	4,01	6,98	9,99
25	2,08	4,02	6,98	9,95
30	2,06	4,00	6,96	9,89
35	2,06	4,01	6,95	9,85
40	2,07	4,02	6,94	9,81
45	2,06	4,03	6,93	9,77
50	2,06	4,04	6,93	9,73
55	2,05	4,05	6,91	9,68
60	2,08	4,10	6,93	9,66
65	2,07*	4,10*	6,92*	9,61*
70	2,07	4,11	6,92	9,57
75	2,04*	4,13*	6,92*	9,54*
80	2,02	4,15	6,93	9,52
85	2,03*	4,17*	6,95*	9,47*
90	2,04	4,20	6,97	9,43
95	2,05*	4,22*	6,99*	9,38*


^{*} extrapolated

Buffer Table, Reagecon


°C	рН				
0°C	*2,01	*4,01	*7,07	*9,18	*12,54
5°C	*2,01	*4,01	*7,07	*9,18	*12,54
10°C	2,01	4,00	7,07	9,18	12,54
15°C	2,01	4,00	7,04	9,12	12,36
20°C	2,01	4,00	7,02	9,06	12,17
25°C	2,00	4,00	7,00	9,00	12,00
30°C	1,99	4,01	6,99	8,95	11,81
35°C	2,00	4,02	6,98	8,90	11,63
40°C	2,01	4,03	6,97	8,86	11,47
45°C	2,01	4,04	6,97	8,83	11,39
50°C	2,00	4,05	6,96	8,79	11,30
55°C	2,00	4,07	6,96	8,77	11,13
60°C	2,00	4,08	6,96	8,74	10,95
65°C	*2,00	*4,10	*6,99	*8,70	*10,95
70°C	*2,00	*4,12	*7,00	*8,67	*10,95
75°C	*2,00	*4,14	*7,02	*8,64	*10,95
80°C	*2,00	*4,16	*7,04	*8,62	*10,95
85°C	*2,00	*4,18	*7,06	*8,60	*10,95
90°C	*2,00	*4,21	*7,09	*8,58	*10,95
95°C	*2,00	*4,24	*7,12	*8,56	*10,95

^{*} Values complemented


Conductivity Concentration Curves




Conductivity Concentration Curves



Conductivity Concentration Curves





MSU 4400(X)-180 Module

Α

Adjustment (COND) 105 Adjustment (CONDI) 139 Adjustment (ORP) 50 Adjustment (OXY) 71 Adjustment (pH) 28

Application in hazardous locations 8

Autoclaving counter (pH/ORP) 21

Automatic buffer recognition (Calimatic) 34

Automatic calibration with standard calibration solution (COND) 110 Automatic calibration with standard calibration solution (CONDI) 144

В

Buffer Tables 162

C

Calibration by data entry of premeasured sensors (COND) 117
Calibration by data entry of premeasured sensors (CONDI) 152
Calibration by data entry of premeasured sensors (OXY) 80
Calibration by data entry of premeasured sensors (pH) 40

Calibration by sampling (COND) 114 Calibration by sampling (CONDI) 148

Calibration by sampling (OXY) 81

Calibration by sampling (pH) 38

Calibration (COND) 104

Calibration (CONDI) 138

Calibration in air (OXY) 76

Calibration in water (OXY) 78

Calibration (ORP) 49

Calibration (OXY) 70

Calibration (pH) 27

Calibration, preset values (COND) 94

Calibration, preset values (CONDI) 128

Calibration, preset values (ORP) 47

Calibration, preset values (OXY) 66

Calibration, preset values (pH) 18

Calibration solution (CONDI parameter setting) 128

Calibration solution (COND parameter setting) 94

Calibration with installation factor (COND) 116

Calimatic 34

MSU 4400(X)-180 Module

Cell constant, sensor monitoring 94 Cell factor, sensor monitoring 128 CIP counter (COND) 96 CIP counter (CONDI) 132 CIP counter (pH/ORP) 21 Combo sensor, pH/ORP, calibration 33 Combo sensor, pH/ORP, parameter setting 17 Concentration, CONDI parameter setting 128 Concentration, COND parameter setting 94 Concentration curves 174 Concentration determination 99 Concentration table (COND) 99 Concentration table (CONDI) 135 Cond calibration/adjustment record 121 Condl calibration/adjustment record 155 Condl messages 156 Condl parameters, range and specification 128 Condl sensor parameter setting 127 Cond messages 122 Cond parameter setting 93 Cond parameters, range and specification 94 Cond sensor wear monitor 121 Connection data 160 Contacting conductivity measurement 93 Corrective maintenance 8

D

Data entry of premeasured sensors 80
Data entry of premeasured sensors (COND) 117
Data entry of premeasured sensors (CONDI) 152
Data entry of premeasured sensors (pH) 40
Decommissioning 159
Device limits, messages (CONDI) 101
Device limits, messages (CONDI) 137
Device limits (messages) (OXY) 69
Device limits, messages (pH) 26
Diagnostic functions (CONDI) 120
Diagnostic functions (CONDI) 155
Diagnostic functions (ORP) 59

MSU 4400(X)-180 Module

Diagnostic functions (OXY) 86 Diagnostic functions (pH) 44 Disposal 159

Ε

Electrostatic discharge (ESD) 13

EMC 160

EN27888, temperature compensation (COND) 94

EN27888, temperature compensation (CONDI) 128

Error messages (COND) 122

Error messages (CONDI) 156

Error messages (OXY) 88

Error messages (pH, ORP) 61

Explosion protection, safety instructions 8

F

Firmware version 9

G

Generating messages (measuring module) (COND) 101 Generating messages (measuring module) (CONDI) 137 Generating messages (measuring module) (OXY) 69 Generating messages (measuring module) (pH) 26

Н

Hardware/firmware version 9 Hydrochloric acid, concentration curve 175

ı

Inductive conductivity measurement 127
Inserting the module 13
Installation factor, calibration (CONDI) 151
Installation, insert module 13
Intended use 6
ISFET sensor, calibration 33
ISFET sensor, operating point 41
ISFET sensor, parameter setting 17

K

KCI, cal presettings (COND) 94 KCI, cal presettings (CONDI) 128

MSU 4400(X)-180 Module

M

Maintenance menu (COND) 119
Maintenance menu (CONDI) 154
Maintenance menu (ORP) 58
Maintenance menu (OXY) 85
Maintenance menu (pH) 43
Manual calibration (COND) 112
Manual calibration (CONDI) 146
Manual entry of buffer values (calibration) 36
Membrane body replacement 85
Memosens OXY 65
Messages, ORP sensor monitoring 48
Messages, pH sensor monitoring 20
Module compatibility 9
Module configuration 15
Module firmware 9

N

NaCI, cal presettings (COND) 94 NaCI, cal presettings (CONDI) 128 Nitric acid, concentration curve 174

0

Oleum, concentration curve 176
One-point calibration 30
ORP adjustment 54
ORP calibration/adjustment record 59
ORP check 56
ORP data entry 53
ORP parameter, range and specification 47
ORP parameter setting 46
Oxy calibration/adjustment record 87
Oxygen measurement 65
Oxy messages 88
Oxy parameter setting 65
Oxy sensor diagram 87
Oxy sensor wear monitor 87

Ρ

MSU 4400(X)-180 Module

Package contents 7 Parameter setting 14 pH calibration/adjustment record 45 pH, ORP messages 61 pH parameter, range and specification 18 pH parameter setting 17 pH sensor diagram 45 pH sensor wear monitor 45 pH value calculation 102 Pressure correction (OXY) 66 Product calibration (COND) 114 Product calibration (CONDI) 148 Product calibration (OXY) 81 Product calibration (pH) 38 R Rated operating conditions 160 Removal from operation 159 Returns 159 S Safety instructions 8 Saline solution, concentration curve 176 Salinity correction (OXY) 66 Sensocheck (COND) 95 Sensocheck (CONDI) 131 Sensocheck (ORP) 48 Sensocheck (OXY) 68 Sensocheck (pH) 20 Sensoface (COND) 95 Sensoface (CONDI) 131 Sensoface (ORP) 48 Sensoface (OXY) 68 Sensoface (pH) 20 Sensor data (CONDI) 130 Sensor data (ORP) 48 Sensor data (OXY) 68 Sensor data (pH) 20 Sensor monitor, diagnostics (COND) 120 Sensor monitor, diagnostics (CONDI) 155

MSU 4400(X)-180 Module

Sensor monitor, diagnostics (ORP) 59 Sensor monitor, diagnostics (pH) 44 Sensor monitoring details (COND) 95 Sensor monitoring details (CONDI) 131 Sensor monitoring details (ORP) 48 Sensor monitoring details (OXY) 68 Sensor monitoring details (pH) 20 Sensor monitor, maintenance (COND) 119 Sensor monitor, maintenance (CONDI) 154 Sensor monitor, maintenance (ORP) 58 Sensor monitor, maintenance (OXY) 85 Sensor monitor, maintenance (pH) 43 Serial number, display 9 SIP counter (COND) 96 SIP counter (CONDI) 132 SIP counter (pH) 21 Sodium hydroxide solution, concentration curve 175 Specifications 160 Sulfuric acid, concentration curve 174 Supplemental directives 2 т Table of contents 3 Technical data 160 Temperature compensation (COND) 94 Temperature compensation (CONDI) 128 Temperature compensation during calibration (COND) 107 Temperature compensation during calibration (CONDI) 141 Temperature compensation during calibration (pH) 31 Temperature compensation of process medium (pH) 24 Temperature probe adjustment (COND) 118 Temperature probe adjustment (CONDI) 153 Temperature probe adjustment (ORP) 57 Temperature probe adjustment (OXY) 84 Temperature probe adjustment (pH) 42 Temp. offset log (COND) 121 Temp. offset log (CONDI) 155 Temp. offset log (ORP) 60 Temp. offset log (OXY) 87

MSU 4400(X)-180 Module

Temp. offset log (pH) 45
Terminal assignments 11
Terminal plate 10
Terminal plate adhesive labels 10
Three-point calibration 31
Two-point calibration 30

U

Ultrapure water, monitoring (COND) 98 Ultrapure water, monitoring (CONDI) 134 Unical, wiring 11 USP function (COND) 98 USP function (CONDI) 134

W

Warnings 2 Wiring 11

Z

Zero correction (CONDI) 150 Zero correction (OXY) 83



Knick Elektronische Messgeräte GmbH & Co. KG

Headquarters

Beuckestraße 22 • 14163 Berlin

Germany

Phone: +49 30 80191-0

Fax: +49 30 80191-200 info@knick.de

www.knick.de

Local Contacts

www.knick-international.com

Translation of the original instructions Copyright 2022 • Subject to change

Version: 2

This document was published on October 11, 2022 The latest documents are available for download on our website under the corresponding product description.



L00772 TA-201.180-KNEN02