The Art of Measuring.

Knick >

SE 605 H

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
Aderbelegung Memosens-Kabel

Safety Instructions	10
Short Description	11
Applications	11
Dimension Drawings	11
Product Code / Calibration	12
Product Code	12
Calibration and Adjustment	12
Specifications (Sensor)	13
Pressure / Temperature Diagram	13
Disposal	13
Specifications (Cable)	14
Memosens Connecting Cable (Accessory)	14
Product Code	14
Wire Assignment of Memosens Cable	14

▲ WARNING! Failure to observe this warning may result in serious injury.

The safety alert symbol on the nameplate means: **Read these instructions for use, observe the Specifications, and follow the Safety Instructions.**

▲ WARNING! Risk of injury from process medium

- Only correct installation (see dimension drawing) ensures a tight assembly.
- Do not dismantle the sensor as long as the process medium is under pressure!

All Applications

Hazards due to pressure, temperature, aggressive media or explosive atmosphere are possible, depending on the location of use. Therefore, the installation, operation, and servicing of the sensor shall only be carried out by suitably trained personnel authorized by the operating company.

Hazardous Areas

Comply with all applicable local codes and standards for the installation of electrical equipment in hazardous locations. For orientation, please refer to IEC 60079-14, EU directives 2014/34/EU and 1999/92/EC (ATEX), NFPA 70 (NEC), ANSI/ISA-RP12.06.01. The electrical and thermal parameters of the sensors must be adhered to.

Memosens Ex sensors are marked by an orange-red ring. Combined with a model CA/MS-***X** or model CA/MS-***X**-L measuring cable or a certified measuring cable which is identical in hardware and function, the sensor may be connected to a suitable measuring device, as described in the Certificates BVS 15 ATEX E141 X and IECEx BVS 15.0114X.

Hazardous Areas: Electrical and Thermal Parameters

Certificate Number:

BVS 16 ATEX E 037 X
IECEx BVS 16.0030X

Marking:

 II 1G
Ex ia IIC T3/T4/T6 Ga

Thermal Parameters:

Temperature class	Ambient temperature range Ta	Maximum permissible process temperature
T6	-20 °C < Ta < +70 °C	70 °C
T4	-20 °C < Ta < +120 °C	120 °C
T3	-20 °C < Ta < +135 °C	135 °C

Special Conditions




- The cable and the sensor shall only be used within the ambient temperature range specified for the temperature class.
- The measuring cable including its connecting head must be protected from electrostatic charging if it passes through areas of Zone 0 (category 1G).
- The Memosens sensors shall not be operated in electrostatically critical processing conditions. Intense vapor or dust flows directly impacting on the connection system shall be avoided.
- Metallic process connection parts must be mounted at the installation site so that they are electrostatically conductive (< 1 MΩ).

Applications

The SE 605H Memosens® is a coaxial 2-electrode conductivity sensor made of stainless steel. It has a large measuring range from ultrapure water to 1000 $\mu\text{S}/\text{cm}$. The outer and inner electrodes can be detached for easy cleaning and simple O-ring replacement. A temperature detector is integrated in the sensor.

The sensor is particularly suitable for monitoring ultrapure water in the pharmaceutical industry, food industry, and biotechnology.

Dimension Drawings

⚠ WARNING! Risk of injury

Never dismantle a sensor with clamp union as long as it is installed.

Replacing the Sensor Gaskets

1. Remove the sensor.
2. Loosen the set screws.
3. Screw off the outer electrode with the Tri-Clamp.
4. Screw off the inner electrode.
5. Replace the gasket.
6. Screw on the inner electrode.
7. Screw on the outer electrode with the Tri-Clamp.
8. Tighten the set screws.
9. Install the sensor.

Product Code


SE 605H -		MS				/	0	0	0
Explosion protection	Without	N							
	ATEX	X							
Process connection	Ingold socket 25 mm (G 1¼")		H	0					
	Ingold socket 25 mm (G 1¼"), 45 mm		H	Z					
	Clamp 1 ½"		J	1					
Sealing material	Elastomeric ring Set F, FKM-FDA				F				
	Elastomeric ring Set E, EPDM-FDA				E				
	Elastomeric ring Set H, FFKM-FDA				H				
	Elastomeric ring Set V, FKM-FDA-USP VI				V				
	Elastomeric ring Set U, EPDM-FDA-USP VI				U				
	Elastomeric ring Set W, FFKM-FDA-USP VI				W				
Certificates	Without								0
	Inspection Certificate 3.1 according to EN 10204								3
	FDA - USP VI								U
	Surface Ra < 0.4 µm								4
	Surface Ra < 0.8 µm								8

Calibration and Adjustment

Remove sensor, clean, rinse, and dry it, then immerse it in a calibration solution.
 Suitable calibration solution CS-C15K/500 (15 µS/cm), CS-C147K/500 (147 µS/cm).
 Perform a product calibration if the sensor can or shall not be removed.
 Observe the user manual of the measuring instrument.

Cell constant:	0.021 /cm
Range:	0 ... 600 $\mu\text{S}/\text{cm}$ ($\pm 2\%$) 600 ... 1000 $\mu\text{S}/\text{cm}$ ($\pm 4\%$)
Material:	<ul style="list-style-type: none"> • Sensor body and electrodes: Stainless steel, 1.4435 • Insulator: PEEK • Gaskets: See product code
Temperature detector:	NTC 30 k Ω
Temperature:	Medium: -30 ... 135 °C, Environment: -25 ... 80 °C
Pressure with Ingold socket:	Max. 25 bar (-30 ... 70 °C); max. 10 bar (135 °C)
Pressure with clamp union:	Max. 16 bar (-20 ... 110 °C); max. 10 bar (135 °C)
Process connection:	See product code
Electrical connection:	Memosens® connector
Dimensions:	See dimension drawing
Weight:	Approx. 0.5 kg

Pressure / Temperature Diagram

Disposal

Observe the applicable local or national regulations concerning the disposal of "waste electrical and electronic equipment".

Memosens Connecting Cable (Accessory)

Socket material	PEEK
Sheath material	TPE, max.130 °C, permanent
Cable diameter	6.3 mm
Wire assignment	See figure below
Length	up to 100 m
Ingress protection	IP 68 (10 m water column, 25 °C, 45 days, 1 mol/l KCl)

Product Line


Cable type	Cable length	Order number
Memosens cable	3 m	CA/MS-003NAA
	5 m	CA/MS-005NAA
	10 m	CA/MS-010NAA
	20 m	CA/MS-020NAA
Memosens cable, Ex*	3 m	CA/MS-003XAA
	5 m	CA/MS-005XAA
	10 m	CA/MS-010XAA
	20 m	CA/MS-020XAA

Other lengths available on request.

*) The cable is Ex-certified, ATEX II IG Ex ia IIC T3/T4/T6.

The Type-Examination Certificate is enclosed with each hazardous-area sensor.

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DE Copyright 2019 • Änderungen vorbehalten
GmbH & Co. KG