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!

NOTICE

Observe the following before initial use
To ensure proper functioning, read these instructions carefully before putting the sensor into operation. Operation and servicing of the sensor should be carried out by specially trained personnel.

Precautions
Despite its exterior robustness, the sensor is a high-precision optical instrument. Therefore, you should take the following precautions when handling the sensor:

- Do not touch the optical parts unnecessarily.
- Avoid any mechanical stress of the optical parts (pressure, scratches).
- Do not the expose the interior of the sensor cap to direct sunlight.

Removing the sensor
Do not turn the sensor in the fitting because this would detach the sensor cap. This would impair the sensor sealing against penetration of process medium and the sensor would be damaged.
To remove the sensor:

- Loosen the PG13.5 nut without turning the sensor.
- Pull the sensor out of the fitting without turning it.

2 Intended Use

The Knick Series SE 740 sensors with integrated wear monitoring are intended for optical oxygen measurement by means of quenching. A temperature detector is integrated. SE 740 sensors are designed for the following applications, among others: Food, pharmaceuticals, fermentations and other biotechnological processes. The sensors are sterilizable and they have a hygienic form factor.

3 Package Contents

- The package content includes:
- Sensor
 - Instructions for use
 - Test report 2.2
 - Inspection certificate 3.1
 - EU Declaration of Conformity

4 Product Marking

The rating plate is printed on the sensor and includes the following information:

- Manufacturer
- Sensor model/length (e.g. SE 740/1)
- Pressure range
- Temperature range
- CE marking

In addition, each sensor is provided with a serial number to allow identification / retraceability. The sensor cap has an individual serial number.

5 Installation and Commissioning

- Inspect the sensor for mechanical damage of the body, the PG 13.5 thread or the sensor cap.
- Connect the sensor to the analyzer. The sensor is immediately ready to use.

6 Measuring

- Take off the transport protection.
- Immerse the sensor sufficiently deep into the process medium (the sensor cap should always be completely wetted). The sensor signal stabilizes within a few minutes. The oxygen sensor enables precise measurements independent of the flow rate.

7 Calibration

The SE 740 oxygen sensor is factory-calibrated. This allows initial start-up. We recommend calibrating regularly and prior to batch processes. Please refer to the user manual of the analyzer for information on calibration.

8 Troubleshooting

Error	Cause	Remedy
No measured values	Sensor not connected to meter	Establish the connection
	Cable defective	Change cable
Measured value too high / too low or error message	Coating on optical surface	Clean the outside of the sensor
	Service life of sensor cap has ended	Check and replace the sensor cap if required
	Dirt inside the sensor cap.	Clean the inside of glass window and sensor cap. Tighten the sensor cap as far as possible,
	Sensor cap untight or defective. O-ring defective.	replace sensor cap and O-ring if necessary.

9 Maintenance and Cleaning

To replace the sensor cap, proceed as follows:

- Clean the outside of the sensor.
- Take off the sensor cap as shown.

Please note: If your fingers do not get a good grip or if the sensor cap is very tight, pull the included silicone tube over the cap.

- If required, clean the cap including the thread.
- Examine the O-ring and replace if required.
- Replace the sensor cap and screw it tight by hand.

The SE 740 sensors withstand the cleaning procedures normally used in biotech industries. Very frequent cleaning with steam or hot alkaline solutions may shorten the service life of the sensor cap.

10 Specifications

Sensor	
Measuring range	4 ppb ... 25 ppm 0.05 ... 300% saturation
Response time (25 °C, air -> N2)	98 % full scale < 30 s
Pressure range P _{rel}	-1 ... 12 bar
Temperature range	-10 ... 130 °C no oxygen measurement at 85 °C and above
O-ring material	EPDM FDA
Process-wetted sensor parts	stainless steel 1.4435, silicone FDA
Temp detector	Yes
Process adaptation	Pg13.5
Sensor connector	M12, 8 pins
Length (1/2/3/4)	120/215/325/425 mm

11 Spare Parts, Cables

Designation	Order code
Sensor cap for SE 740	ZU 0914
Cable, 5 m, M12, 8 pins	CA/M12-005N485
Cable, 10 m, M12, 8 pins	CA/M12-010N485

12 Warranty

12-month warranty after date of shipment in the event of manufacturing faults.

13 Disposal

Please observe the applicable local or national regulations concerning the disposal of "waste electrical and electronic equipment".

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English

- M12 connection, 8 pins
- Serial number / length
- Pg 13.5 thread
- Compression ring
- O-ring
- Sensor body
- Sensor cap with serial number
- Transport protection

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Elektronische Messgeräte GmbH & Co. KG