



**pH Portables and pH Sensors – Product Range**




		Order No.
<p>Portavo 907 pH</p> 	<p>Portavo 907 for measurement with analog pH sensors or digital sensors for pH, conductivity or oxygen, incl. Paraly SW 112 configuration software with USB connector cable.</p>	<p>907 Multi pH</p>
<p>Portavo 904 pH</p> 	<p>Portavo 904 for measurement with analog or digital Memosens pH sensors, incl. Paraly SW 112 configuration software with USB connector cable.</p>	<p>904 pH</p>
<p>Portavo 902 pH</p> 	<p>Portavo 902 for measurement with analog or digital Memosens pH sensors</p>	<p>902 pH</p>
<p>pH set</p>		
	<p>pH set consisting of Portavo 902 pH and one SE 102 NMS digital Memosens pH sensor with cable</p>	<p>902 pH</p>
<p>pH/Pt 1000 sensor</p> 	<p>Digital Memosens pH sensor plastic body, ceramic junction, length 120 mm</p>	<p>SE 101 NMS</p>
<p>pH/Pt 1000 sensor</p> 	<p>Digital Memosens pH sensor glass body, ceramic junction, length 110 mm</p>	<p>SE 102 NMS</p>
<p>Memosens cable</p> 	<p>Measuring cable for digital sensors with Memosens connector, 1.5 m</p>	<p>CA/MS-001XFA</p>

**Accessories – Product Range**








pH/Pt 1000 sensor		Order No.
	Plastic body, fiber junction, length 120 mm	SE 101 N
pH/Pt 1000 sensor		
	Glass body, ceramic junction, length 110 mm	SE 102 N
pH puncture sensor		
	Plastic body, open junction, length 65 / 25 mm	SE 104 N
pH/Pt 1000 sensor		
	for measurements in hazardous areas Zone 0, incl. equipotential bonding cable, glass body, ceramic junction, length 105 mm	ZU 6979
Electrode quivers		
	5 quivers, as replacement, for leak-proof storage of sensors	ZU 0929
Pt 1000 temperature detector		
	for fast response temperature measurements: Monel 2.4360, -10 ... +100 °C, accuracy class A to IEC 751	ZU 6959
Pt 1000 temperature detector		
	for fast response temperature measurements, with tilted tip for puncture measurements in semi-solid substances: Monel 2.4360, -10 ... +100 °C, accuracy class A to IEC 751	ZU 0156

**Accessories and Buffer Solutions – Product Range**

Robust field case		Order No.
	For meter and sensor	ZU 0934
Li-ion battery		Order No.
	Li-ion battery (USB chargeable, for Portavo 904 and 907 only)	ZU 0925

CaliMat pH Buffer Solutions		Quantity	Order No.
	pH 2.00 (20 °C)	250 ml	CS-P0200/250
	pH 4.00 (20 °C)	250 ml	CS-P0400/250
		1000 ml	CS-P0400/1000
	pH 7.00 (20 °C)	250 ml	CS-P0700/250
		1000 ml	CS-P0700/1000

**Buffer Solutions Product Range – Continued**

	pH 9.00 (20 °C)	250 ml 1000 ml	CS-P0900/250 CS-P0900/1000
	pH 12.00 (20 °C)	250 ml	CS-P1200/250
	Set pH 4.00 (20 °C)	3 x 250 ml	CS-PSET4
	Set pH 7.00 (20 °C)	3 x 250 ml	CS-PSET7
	Set pH 9.00 (20 °C)	3 x 250 ml	CS-PSET9
	Set pH 4.00 / 7.00 / 9.00 (20 °C)	3 x 250 ml	CS-PSET479
	KCl solution	250 ml	ZU 0062