

MemoSuite Advanced

ChangeLog

Dokument-Version 1.0

Knick Elektronische Messgeräte GmbH & Co. KG

www.knick.de

Content

Content.....	2
1 Versions-Historie	2
1.1 Version 2.0.0.....	2
1.2 Version 2.1.0.....	2
1.3 Version 2.2.0.....	3
1.4 Version 2.3.0.....	3
1.5 Version 2.4.0.....	3
1.6 Version 2.5.0.....	3
1.7 Version 2.5.1.....	3
1.8 Version 2.5.2.....	3

1 Versions-Historie

1.1 Version 2.0.0

New: Editing Memosens measuring point TAG and number

New: Importing MemoLogs from Portavo

New: Hiding (deactivating) sensors

New: Database backup and restore

New: Multi-sensor operation

New: Support for Hamilton Memosens sensors

New: Calibration comments in calibration reports

New: French user language

New: Calibration reports for external calibrations

New: Chinese user language

1.2 Version 2.1.0

New: Windows 10 support

1.3 Version 2.2.0**Adjustable limits for sensor evaluation during calibration/adjustment**

The limits of zero point and slope within which a calibration/adjustment is considered successful can now be set by the user (for the pH calibration types Calimat, Manual, and Reference).

1.4 Version 2.3.0

New: Measurement reports (PDF)

New: Configurable limits for zero point/slope evaluation (GMP calibration)

New: Sensor TAG in the calibration protocol

1.5 Version 2.4.0

New: Login and access permissions via ActiveDirectory

1.6 Version 2.5.0

New: Conductivity measurement value, temperature-compensated

1.7 Version 2.5.1

Update of the MySQL DB used to 5.7.44 (security update)

New: Windows 11 support

1.8 Version 2.5.2

Bug fix: Generation of calibration reports under Windows 11